[Kurzanleitung](#) | [Short instructions](#) | [Instruction brève](#) | [Guía rápida](#)

EASYCAL™ 5

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Graph der Messwerte

Graph showing measurement values (µl) over time. The y-axis ranges from 85 to 115 µl. The x-axis shows measurement numbers 0, 2, and 3. A dashed horizontal line is at 100 µl. Three data points are plotted at approximately 100, 100, and 100 µl.

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

The screenshot shows the EASYCAL™ 5 software interface. At the top, it displays the date and time: 2021-11-11 11:03:25 - DI / ISO-69 - 704780 BRAND GMBH + CO KG Transferpette® S Variabel 95 - 1000µl. The interface is divided into several sections:

- Umgebungsdaten (Environmental Data):** Shows temperature (21 °C), humidity (40 %rH), and pressure (900 mibar).
- Soll-Wert (Target Value):** 1000 µl.
- Anzahl (Count):** 2 / 10.
- Letzte Messungen (Last Measurements):** A list of previous measurements with icons for volume, mass, and time.
- Neue Messung (New Measurement):** A form to enter a new measurement.
- Prüfmittel (Test Material):** A table with columns for Nummer (Number), Zeitpunkt (Time Point), Messwert (Measurement Value), and Volumen (µl). The table shows two entries: 1 (11.11.2021 14:48:43, 1.000 g, 1003 µl) and 2 (11.11.2021 14:49:26, 1.000 g, 1103 µl).
- Graph der Messwerte (Measurement Value Graph):** A scatter plot showing the measurement values for the two entries.
- Statistik Gerät (Device Statistics):** Shows the mean value (1053,16 ± 0) and a range from 992 to 1008.
- Statistik Auftrag (Order Statistics):** A table of statistical data:

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The screenshot shows a dialog box titled 'Messungen verwerfen' (Reject Measurements). The dialog contains the following text and elements:

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info@brand.de

www.brand.de

More operating manuals and translations?

Please refer to <http://www.brand.de/om> or use the following Quick Response Code:

The original operating manual is in German. Other languages are translations of the original operating manual.

Subject to technical changes, errors, and misprints.

Contents

1	Preliminary remarks	35
2	Installation.....	36
3	First steps in EASYCAL™ 5.....	45
3.1	Tutorial.....	45
3.1.1	Set up balance as test equipment.....	45
3.1.2	Creating order	45
3.1.3	Determine test equipment	50
3.1.4	Performing the test	50
3.1.5	Complete test	58

1 Preliminary remarks


This tutorial is intended to make it easier for you to use EASYCAL™ 5 as part of the 30-day test period. Using an exemplary process, you will learn how to fill different areas of EASYCAL™ 5, how to carry out the first device test and how to print out a certificate.

Please note that in the course of this introduction, sample data will be generated in the database, which cannot be deleted later. We recommend coordinating this with the regulatory requirements of your quality management system.

2 Installation


Prerequisite:

- ✓ You have already received the EASYCAL™ 5 installation file, EasyCal-Setup.exe, from the BRAND home-page and saved it locally.
 - ✓ You have also received a test license (ticket) by E-mail.
 - ✓ You have local admin rights and an existing internet connection.
- a. Start the installation by opening the EasyCal-Setup.exe and follow the instructions in the installation wizard.


→ EASYCAL™ 5 and all software components required for operation are installed on your PC.

- b. After completing the installation, start EASYCAL™ 5 by ticking the box and clicking on Complete.


c. You will be asked to activate EASYCAL™ 5.


d. Click on Activate now.


→ Your standard web browser opens with the page <https://license.brand.de/index.php>.

e. Enter your ticket (= license key) here.

f. Click on Continue.

g. Click on Activate License.

Activate Licenses

To activate your licenses:

1. Select the licenses you want to activate.
2. Select a license file to load the licenses into.
3. Click "Activate Selected Licenses Now".

<input checked="" type="checkbox"/>	Name	Activated On	Serial number	Status
<input checked="" type="checkbox"/>	Easycal 5 Demo Lizenz (30 Tage)	-		Available

Select license file

select / create automatically

[Activate Selected Licenses Now](#) [Manual license transfer](#)

h. If you do not change any of the settings, click on Activate Selected Licenses Now.

↪ The online license transfer is carried out.

Online License Transfer

Starting license transfer.
 Downloading license template.
 Registering license template.
 Creating license request.
 Downloading license update.
 Importing license update into target.
 Creating receipt.
 Uploading receipt.

License transfer completed successfully!

- i. Click on OK.

BRAND. For lab. For life.

My Licenses


Your license code contains licenses for the following products. If you want to activate those licenses, click "Activate Licenses". If you want to relocate licenses from one computer to another, click "Deactivate / Relocate Licenses".

Name	Activated On	Serial number	Status
Easycal 5 Demo Lizenz (30 Tage)	2021-11-15 10:28:54	130-1018824202	Activated

Deactivate / Relocate Licenses

- j. The license is now activated. You can close the browser window.

- k. Go to the product activation window.

Product activation required

Thank you for installing EASYCAL (V. 5.0.14.211).


EASYCAL has not been activated yet. In order to use the product, it is necessary that you activate it with a license. This requires an Internet connection on this or another computer.

Activate now

Check again

Cancel and close

- l. Click "Check Again".


Period of use not yet started


The period of use of the license for EASYCAL has not yet been started

OK Cancel


- m. Click on OK to start the usage period. An information window appears.


- n. Click on OK.
- o. At the same time, EASYCAL™ 5 is started and a notification appears.


- p. Acknowledge the notification by clicking Yes so that the database is initially filled.
- q. A notification will appear after a short time.


r. Close your entries by clicking Continue.

- s. The new data is imported and the progress is displayed.


- ↳ After a while, a notification will appear that the update is complete.

- t. Acknowledge these by clicking on OK.

- ↳ EASYCAL™ 5 will start. Various messages will appear to show the progress.

- ↳ A warning appears that you are using a demo version. Read the notes.

- u.** After you acknowledged the warning you can now test EASYCAL™ 5.
- ↪ You will be automatically logged in as "EASYCAL Admin". Note that EASYCAL™ 5 will always start in single user mode after installation. You can change this later in the "Global Settings", but shouldn't do so at this point in time.

3 First steps in EASYCAL™ 5.

3.1 Tutorial

This example shows the process of a test – from the creation of an order to the test to the printing of the certificate. The section refers to master data available in the demo license and serves as a first introduction to the program.

Login

If the program EASYCAL™ 5 is started, a log-in dialog appears. Here, the user logs in with user name and password.

3.1.1 Set up balance as test equipment

In order set up a balance as a piece of test equipment, proceed as follows. The set-up is always done in the order: Set up test equipment in the master data, set up test equipment at the workstation, and select test equipment for the orders.

Prerequisite:

- ✓ You have connected a balance to the workstation.
 - ✓ You know the connection parameters of the balance.
 - ✓ You know to which COM port the balance is connected at the workstation.
- a. Name balance (see Administration > Master data maintenance > Test equipment > Test equipment).
 - b. Set up location for test equipment (if not already done).
 - c. Configure the balance as a test equipment (see Test equipment device type). In doing so, you specify both the function and the connection of the balance to the workstation (see Measuring function, Communication).

Exit this menu. Change to the Settings menu.

- d. Register the balance as test equipment at the workstation (see Settings > Workstation settings > Test equipment).
- e. Select balance and the connection for the workstation. You can also test the connection to the balance.

Exit this menu. Change to the Orders and calibration slips menu.

- f. You can select the connected balance in the order when selecting the test equipment. In the tutorial (see Determine test equipment, p. 50), in the user manual under Select test equipment for the test.

3.1.2 Creating order

- a. Navigate from Main menu to Orders and calibration slips.

Orders and calibration slips

Test type filter: Not specified | Search: Type to search | from 1/11/2022 to 1/26/2022

Order overview | Calibration slips

Order I	Test type	Order number	Article	Status
2	SPEC	2022-01-25-15-2	EC-12484 - Burette BLAUBRAND®, straight stopcock, Schilling, class AS; 0.02 - 10 ml	Created
1	STD	123456	EC-704106 - Transferpette® Fixed volume; 5 µl	Finished

Buttons: Back, Individual device test, Create order, Open

b. Click on Create order.

Select order type

Batch test

Use when testing random samples (of any size). Batch certificates are issued.

Individual test

Use for 100% test or test of a single device. Each device receives its own certificate.

Button: Back

c. Click on Individual test.

d. For Test type, select the value STD – Standard values.

Test types organize your test procedures. You cannot create additional test types; instead, you must choose from the test types already stored. The test plans are linked to the test type as well as the certificate templates. In the user administration, the user must be given the respective access right for the corresponding test type. Because you are currently logged in as Global Admin, you have access to all data.

Individual test - Create order

Test type: **STD - standard values** | Number of devices: **1** | Customer: [dropdown]

Order no.: [input] | Order no.: [input]

Article code: [input]  | OC-Number: [input]

Article: [input] | Prepare calibration slips

Pre- and post-test required

Notes on the order: [text area] | End customer: [dropdown]

Volume ranges: [input] 

Test plan: [dropdown] 

Back | Start test

e. Click on the button  at Item code in order to display already created articles.

Individual test - Create order

Test type: **STD - standard values** | Number of devices: **1** | Customer: [dropdown]


Order no.: [input] | Order no.: [input]

Article code: [input]  | OC-Number: [input]

Article: [input] | Prepare calibration slips

Pre- and post-test required

Notes on the order: [text area] | End customer: [dropdown]

Volume ranges: [input] 

Test plan: [dropdown] 

Back | Start test

Select article

Search: [input]


Article code	Vei	Manufacturer	Product	Product variation	Volume range	Status
10_copy	1	BRAND GMBH + CO KG	Dispensette® S Organic	Fixed volume, Special fixed volume	0.05 - 5 ml	Released
10002	0	BRAND GMBH + CO KG	Burette	SILBERBRAND, ohne Hahn, class B, Economy	0.1 - 25 ml	Released
10012	0	BRAND GMBH + CO KG	Compact Burette	SILBERBRAND, class B, Economy	0.1 - 25 ml	Released
100210	6	Eppendorf AG	Research® plus	Fixed volume, medium grey	10 µl	Released
100211	7	Eppendorf AG	Research® plus	Fixed volume, light grey	20 µl	Released
100211_copy	0	Eppendorf AG	Research® plus	Fixed volume, light grey	20 µl	Released
100212	6	Eppendorf AG	Research® plus	Fixed volume, yellow	10 µl	Released
100212_copy	0	Eppendorf AG	Research® plus	Fixed volume, yellow	10 µl	Released
100213	6	Eppendorf AG	Research® plus	Fixed volume, yellow	20 µl	Released
100214	7	Eppendorf AG	Comforpette®	Fixed volume	20 µl	Released
100215	6	Eppendorf AG	Research® plus	Fixed volume, yellow	50 µl	Released
100216	6	Eppendorf AG	Research® plus	Fixed volume, yellow	100 µl	Released
100217	6	Eppendorf AG	Research® plus	Fixed volume, yellow	200 µl	Released
100218	6	Eppendorf AG	Research® plus	Fixed volume, blue	200 µl	Released
100219	6	Eppendorf AG	Research® plus	Fixed volume, blue	250 µl	Released
100220	6	Eppendorf AG	Research® plus	Fixed volume, blue	500 µl	Released
100221	6	Eppendorf AG	Research® plus	Fixed volume, blue	1000 µl	Released
100310	6	Eppendorf AG	Research® plus	Adjustable volume, dark grey	0.1 - 2.5 µl	Released
100311	8	Eppendorf AG	Research® plus	Adjustable volume, dark grey	0.5 - 10 µl	Released

Daniel Hing | Last status: The application is prepared


f. Enter Transferpette® S 1000 in the search field. As you type, you will see the article selection change. This behavior is the same in all search fields of EASYCAL™ 5. We recommend always using the search function first because EASYCAL™ 5 comes with an extensive master database.

The prefix EC- indicates that it is a prefabricated and unchangeable data set. If you would like to change

these records, copy the record, and edit and save it under a different name. For more information, see “Help” in EASYCAL™ 5.

- g. Select the article code EC-704780. Confirm your selection with a double click.
- h. For Test plan, select the appropriate test plan. As a rule, you will find 1–2 entries here.
- i. The general conditions for the test as well as the defect limits of the assigned article are stored in the test plan.


- j. Click on Start test.

A notification that no order number has been assigned appears.


If there are specifications in your company/quality management process for organizing tests, select no, and assign an order number.

If there are no specifications for organizing tests, select Yes and EASYCAL™ 5 will assign a unique order number

We recommend deciding on a uniform approach at an early stage.

k. Continue the tutorial by clicking on Yes. The job number is created based on the current time stamp.

l. You can enter a serial number for the identification of the device. This is necessary for real tests. However, it is not necessary for the tutorial.

m. Click on Continue. The test order has been created.

The window for starting the actual test will open.

3.1.3 Determine test equipment

Air temperature	Thermometer		351000-71
Temperature of the medium	Thermometer		351000-71
Relative humidity	Hygrometer		600029-10
Air pressure	Barometer		600029-10
Weight	Single-channel balance	Calibrate now	20210708-RZ-01
Weight	Single-channel balance	Calibrate now	20210708-RZ-01
Weight	Single-channel balance	Calibrate now	20210708-RZ-01


Prerequisite:

- ✓ The order is created.
- a. The necessary parameters as well as the test equipment type and the test equipment (test equipment number) are preselected for the test equipment. Because only one piece of test equipment is created, we leave it at this selection.

3.1.4 Performing the test


- a. Continue by clicking on Next characteristic.

You see the central input screen of the measurement functionality of EASYCAL™ 5.

- b. Enter fictitious weight values for Neue Messung. Start with 0, 1, and click on “Apply measured value”. A window in which you must enter or complete the measured values appears.

The selection of the separator (comma or dot) is done via the language selection in the software.

- c. Enter fictitious values for the environment parameters as follows:

Enter measured values ✕

Temperature difference	0.00 °C		0.5 °C		
Weight	<input type="text" value="0.1"/>	g		<input type="checkbox"/>	<input checked="" type="checkbox"/> Manual
Air temperature	<input type="text" value="21"/>	°C	(20°C - 25°C)	<input type="checkbox"/>	<input checked="" type="checkbox"/> Manual
Temperature of the medium	<input type="text" value="21"/>	°C	(20°C - 25°C)	<input type="checkbox"/>	<input checked="" type="checkbox"/> Manual
Relative humidity	<input type="text" value="40"/>	%rh	(40%rh - 65%rh)	<input type="checkbox"/>	<input checked="" type="checkbox"/> Manual
Air pressure	<input type="text" value="961"/>	mbar	(960mbar - 1025mbar)	<input type="checkbox"/>	<input checked="" type="checkbox"/> Manual

Cancel
Apply

d. Click on Apply. The first measurement is hereby completed. The measured value and the volume calculated from it appear in the table.

Test equipment > 2. TV > 1. TV > NV

Number	Time	Measured value	Volume [µl]	Sta
1	11/12/2021 2:06:42 PM	0.10000 g	100.3	✓

e. Now enter the value 0.1 g as the measured value eight times in the same way. The table is filled up.

BRAND - EASYCAL™ 5

2021-11-12-14-10-58-DI / ISO-74 - - 704780 BRAND GMBH + CO KG Transferpette® S Adjustable volume 95 - 1000µl

Ambient data

21 °C 21 °C 40 %RH 961 mbar

Request ambient data

Target value: 100 µl

Number: 9 / 10

Letzte Messungen: 0.1 g, 1003.01211312 µl/g, 100.3 µl, 2:13:39 PM

Neue Messung: g, µl/g, µl,

Test equipment disconnected

Fast measurement
 Automatic printing
 Automatic tare
 Allow all measured values
 Check measured values before entering

Test equipment > 2. TV > 1. TV > NV

Number	Time	Measured value	Volume [µl]	Sta
2	11/12/2021 2:13:27 PM	0.10000 g	100.3	✓
3	11/12/2021 2:13:29 PM	0.10000 g	100.3	✓
4	11/12/2021 2:13:30 PM	0.10000 g	100.3	✓
5	11/12/2021 2:13:32 PM	0.10000 g	100.3	✓
6	11/12/2021 2:13:34 PM	0.10000 g	100.3	✓
7	11/12/2021 2:13:35 PM	0.10000 g	100.3	✓
8	11/12/2021 2:13:37 PM	0.10000 g	100.3	✓
9	11/12/2021 2:13:39 PM	0.10000 g	100.3	✓

Graph of the measured values

Statistics of device

100.3 ± 0

Accuracy: 0% / 8% ✓

Coefficient of variation: 0% / 3% ✓

Mean: 100.3 µl

Standard deviation: 0 µl

Expanded measurement uncertainty: 0 µl

Statistics of order

Accuracy: 0%

Mean: 100.3 µl

Standard deviation of the mean: 0 µl

Coefficient of variation of accuracy: 0%

Mean value CV: 0%

CPK: 0

f. The light blue color in the section above the table shows that the order is being processed. Click the icon in order show and hide the extended certificate information.

g. Enter 0.1 again as the measured value, and pay attention to the progression display above the measured value table.

BRAND - EASYCAL™ 5

2021-11-12-14-10-58-DI / ISO-74 - - 704780 BRAND GMBH + CO KG Transferpette® S Adjustable volume 95 - 1000µl

Ambient data

21 °C 21 °C 40 %rH 961 mbar

Request ambient data

Target value
100 µl

Number
10 / 10

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0.1		g
1003.01211312	1003.01211312	µl / g
100.3		µl
2:17:42 PM		

11 / 11

Fast measurement

Automatic printing

Automatic tare

Allow all measured values

Check measured values before entering

Discard entry Undo discard

Test equipment: **2. TV** ✓ 1. TV 1. TV NV

Number	Time	V
3	11/12/2021 2:13:29 PM	0.
4	11/12/2021 2:13:30 PM	0.
5	11/12/2021 2:13:32 PM	0.
6	11/12/2021 2:13:34 PM	0.
7	11/12/2021 2:13:35 PM	0.
8	11/12/2021 2:13:37 PM	0.
9	11/12/2021 2:13:39 PM	0.
10	11/12/2021 2:17:42 PM	0.

Graph of the measured values

Statistics of device


100.3 ± 0

92 100 108

Back

Previous device Previous characteristic Next characteristic

- h. With input of the 10th measured value, all required measurements for characteristic 2. TV have been performed. The characteristic gets a green check. You can enter further measured values; however, this is not required.


i. Click on 1. TV in the progress bar in order to jump to the next characteristic.

j. For Neue Messung, enter the following values in sequence:

- 0.498
- 0.499
- 0.500

Repeat twice.

➔ You can see the position of the individual measured values in relation to defined error limits (solid lines) and the calculated mean value in the display Graph of the measured values. Under Statistics of device, you can see the calculated actual values for accuracy and coefficient of variation compared with the defined limits.

- k. Enter 0.498 again as the measured value, confirm the environmental parameters, and click NV in the progress display. Then enter 1 as the measured value.


The screenshot shows the EASYCAL™ 5 software interface. At the top, it displays the date and time: 2021-11-12-14-10-58-DI / ISO-74 - - 704780 BRAND GMBH + CO KG Transferpette® 5 Adjustable volume 95 - 1000µl. The ambient data section shows 21 °C, 21 °C, 40 %H, and 961 mibar. The target value is 1000 µl, and the number of measurements is 1 / 10. The last measurement is 1003.01 µl. The measured value is 1003.01 µl. The graph shows the measured value of 1003.01 µl. The statistics of device show Accuracy: 0.3 % / 0.8 %, Coefficient of variation: 0 % / 0.3 %, Mean: 1003.01 µl, Standard deviation: 0 µl, and Expanded measurement uncertainty: 0 µl. The statistics of order show Accuracy: 0.3 %, Mean: 1003.01 µl, Standard deviation of the mean: 0 µl, Coefficient of variation of accuracy: 0 %, Mean value CV: 0 %, and CPK: 0.

- l. We assume that you have made a typing error. To do this, enter 1.1 as the measured value, and confirm the environmental parameters. The following warning appears:

The warning dialog box has a title bar 'Warning' and a close button. It contains a yellow warning triangle icon and the text: 'The following errors occurred during the measurement Channel 1: The statistic limits were exceeded.' There is an 'OK' button at the bottom.


- m. Acknowledge this warning by clicking OK.

You can see that the different areas of the input screen react according to the new measured value. You cannot delete this measured value; however, you can discard it. The value remains visible but is no longer used for statistics.


- n. Click on the number of the measurement in the table, and then click on “Discard entry”.

A message appears.


- o. Specify a reason (e.g., Identified as outliers.). You can also select a reason from the drop-down menu. These texts can be maintained in the master data. Confirm your selection by clicking on Continue.

p. The measured value 1 . 1 is not considered further in the evaluation.

q. Continue entering measured values until 10/10 is displayed for “Number”. When the last required measured value is entered, the color in the section above the table changes from light blue to yellow. This signals to you that the order is complete.

Furthermore, all characteristics in the progress bar are marked with a green check. The discarded measured value is highlighted in pink.

3.1.5 Complete test

Test equipment >> 2. TV >> 1. TV >> NV

Number	Time	Measured value	Volume [µl]	Sta
1	11/12/2021 2:39:32 PM	1.00000 g	1003.01	✓
2	11/12/2021 2:40:16 PM	1.10000 g	1103.31	✓
3	11/12/2021 3:29:39 PM	1.00000 g	1003.01	✓
4	11/12/2021 3:29:41 PM	1.00000 g	1003.01	✓
5	11/12/2021 3:29:43 PM	1.00000 g	1003.01	✓
6	11/12/2021 3:29:44 PM	1.00000 g	1003.01	✓
7	11/12/2021 3:29:46 PM	1.00000 g	1003.01	✓
8	11/12/2021 3:29:47 PM	1.00000 g	1003.01	✓
9	11/12/2021 3:29:50 PM	1.00000 g	1003.01	✓
10	11/12/2021 3:29:50 PM	1.00000 g	1003.01	✓

Statistics of device

Accuracy: 0.3 % / 0.8 % ✓
 Coefficient of variation: 0 % / 0.3 % ✓
 Mean: 1003.01 µl
 Standard deviation: 0 µl
 Expanded measurement uncertainty: 0 µl

Statistics of order

Accuracy: 0.3 %
 Mean: 1003.01 µl
 Standard deviation of the mean: 0 µl
 Coefficient of variation of accuracy: 0 %
 Mean value CV: 0 %
 CPK: 0

Buttons: Back, Previous device, Previous characteristic, Next characteristic, Additional device, Device evaluation, **Complete test**, Print preview

a. Click on Complete test.

Order 2021-11-12-14-10-58-DI

Test type: STD - standard values | Article: 704780 - Transferpette® S Adjustable volume: 95 - 1000 µl
 Order no.: 2021-11-12-14-10-58-DI | Test plan: 4437.0 - EC-STD-ISO-PIpip-var-LHDiv
 Order type: Individual test | Tester: D. Ihrig
 Lot size: 1
 Random sample size: 1 / 1
 Status: Finished | Order no.: | Customer: | OC-Number: | End customer: |

Test characteristic

Test characteristic	Description	Target va	Mean	Tol target	Tol OK	Standard devi	A% target	A%	A% OK	CV% targ	CV%	CV% OK	QS	CPK	Min CPK	CPK OK
1	2. TV	100µl	100.3 µl	OK	0 µl	8 %	0 %	OK	3 %	0 %	0 %	OK	0	0	OK	
2	1. TV	500µl	500.4 µl	OK	0.88 µl	1.6 %	0.1 %	OK	0.6 %	0.18 %	0.18 %	OK	0	0	OK	
3	NV	1000µl	1003.01 µl	OK	0 µl	0.8 %	0.3 %	OK	0.3 %	0 %	0 %	OK	0	0	OK	


Device/calibration slip overview

Calibration slip type	Internal no.	Calibration slip no.	Serial number	Article code	Conducted tests	Status	Device status	Tester	Linked calibration slip
Individual test	13014	ISO-74		704780	30	Finished		D. Ihrig	

Buttons: Evaluation/release, **Print internal log**, Print certificates, New calibration slip


Buttons: Back, Back to the order overview, Save

b. You will now see the order with a summary of the measurements for the respective test characteristics. Click on Print internal log.

c. You will see the certificate editor in view mode. Here you can print or export the stored report.

d. Exit the certificate editor by clicking Back.

e. Click on Evaluation/release.

Evaluation/release
✕

Current status Finished

New status Released

Remarks

Apply to all

Cancel

Save

f. In the window that opens, select Released, and confirm the entry with Saving...

BRAND - EASYCAL™ 5
✕

Order 2021-11-12-14-10-58-DI

Test type	STD – standard values	Article	704780 - Transferpette® S Adjustable volume; 95 - 1000 µl	Notes on the order		Communication with the customer	
Order no.	2021-11-12-14-10-58-DI	Test plan	4437.0 - EC-STD-ISO-Pip-Var-LHDIV				
Order type	Individual test	Tester	D. Ihrig				
Lot size	1	Order no.	<input type="text"/>	Customer	<input type="text"/>		
Random sample size	1 / 1	OC-Number	<input type="text"/>	End customer	<input type="text"/>		
Status	Released						

Test characteristic

Test characteristic	Description	Target	va	Mean	Tol	target	Tol	OK	Standard dev	A%	target	A%	A% OK	CV% targ	CV%	CV% OK	QS	CPK	Min CPK	CPK OK	CPK OK
1	2. TV	100µl		100.3 µl				OK	0 µl	8 %	0 %	OK	3 %	0 %	0 %	OK	0	0		OK	
2	1. TV	500µl		500.4 µl				OK	0.88 µl	1.6 %	0.1 %	OK	0.6 %	0.18 %	0.18 %	OK	0	0		OK	
3	NV	1000µl		1003.01 µl				OK	0 µl	0.8 %	0.3 %	OK	0.3 %	0 %	0 %	OK	0	0		OK	

Device/calibration slip overview

Calibration slip type	Internal no.	Calibration slip no.	Serial number	Article code	Conducted tests	Status	Device status	Tester	Linked calibration slip
Individual test	13014	ISO-74		704780	30	Released		D. Ihrig	

Back

Evaluation/release


Print internal log

Print certificates


New calibration slip

Save


g. Click on Print certificates.


h. Click the check-box, and then click Print preview.


→ The certificate editor opens in Print preview mode. You can still make settings in the right section and enter export annotations depending on the selected content structure.

- i. Here you can also export the certificate as a PDF or print a hard copy. By clicking on Back, you leave the Print preview.

- j. Click on Back to the order overview
- k. From here you can create additional jobs and start tests. By clicking on Back, you can return to the main menu.

➤ This is where this tutorial ends.

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BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

Do you need more operating manuals and translations?
Please refer to <http://www.brand.de/om> or use the following Quick
Response Code:

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The screenshot shows the 'Imprimer les certificats' dialog box. It has a title bar with a close button. Below the title, there is a checkbox labeled 'Tout sélectionner / désélectionner' which is checked. A list box below it contains one item: '19 STD-19', which is also checked. At the bottom of the dialog, there are four buttons: 'Imprimer' (highlighted with a red box), 'Enregistrer', 'Enregistrer et imprimer', and 'Interrompre'. A fifth button, 'Aperçu avant impression', is located at the bottom right of the dialog.

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The screenshot displays the EASYCAL™ 5 software interface in a print preview mode. The main window shows a report titled 'EASYCAL™ 5 Resultados de la auditoría'. The report content is as follows:

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The interface includes a left sidebar with a table of contents, a right sidebar with settings (Ajustes), and a bottom toolbar with options like 'Imprimir', 'Guardar como PDF', and 'Guardar e imprimir'.

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