[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |  
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |  
[Руководство по эксплуатации](#) | [사용지침](#) | [Bruksanvisning](#) | [使用説明書](#) |  
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |  
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

## seripettor® seripettor® pro

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# Table of contents

<b>1 Introduction .....</b>	<b>48</b>
1.1 Scope of delivery .....	48
1.2 Terms of use .....	48
<b>2 Safety Instructions .....</b>	<b>50</b>
2.1 General safety instructions.....	50
2.2 Function.....	50
2.3 Limitations of Use .....	51
2.4 Operating Limitations.....	51
2.5 Operating Exclusions .....	52
2.6 Storage Conditions .....	52
2.7 Recommended Application Range .	52
2.8 Dispenser Selection Chart .....	54
<b>3 Functional and operating elements .....</b>	<b>56</b>
<b>4 Assembly .....</b>	<b>57</b>
4.1 First Steps.....	57
<b>5 Operation .....</b>	<b>61</b>
5.1 Dispensing .....	61
5.2 Dispensing with the flexible dis- charge tube .....	62
5.3 Dispensing with a microfilter .....	64
5.4 Bottle stand .....	65
<b>6 Error limits .....</b>	<b>66</b>
<b>7 Checking the Volume (Calibration) .....</b>	<b>67</b>
<b>8 Cleaning .....</b>	<b>68</b>
8.1 Basic cleaning .....	68
8.2 Cleaning and replacing valves.....	69
8.3 Dislodging a stick valve ball .....	72
8.4 Replacing the dispensing cartridge	72
<b>9 Autoclaving .....</b>	<b>74</b>
9.1 Preparing for sterile work .....	75
<b>10 Troubleshooting .....</b>	<b>76</b>
<b>11 Product markings .....</b>	<b>77</b>
<b>12 Ordering Information .....</b>	<b>78</b>
<b>13 Accessories/spare parts .....</b>	<b>79</b>
<b>14 Repairs.....</b>	<b>86</b>
14.1 Sending for repair .....	86
<b>15 Calibration service .....</b>	<b>88</b>
<b>16 Warranty.....</b>	<b>89</b>
<b>17 Disposal .....</b>	<b>90</b>

# 1 Introduction

## 1.1 Scope of delivery

### seripettor®

Actuating unit seripettor®, valve block with discharge tubes (pre-assembled) for GL 45 threaded bottles, aspiration tube, 2x dosing unit, 3x bottle adapter (PP) and this operating manual.

### seripettor® pro

Actuating unit seripettor® pro, valve block with discharge tubes (pre-assembled) for GL 45 threaded bottles, telescoped aspiration tube, 2x dosing unit, 3x bottle adapter (PP) and this operating manual.

	Adapter for bottle thread	Filling tube length
seripettor®	GL 32-33, GL 38, S 40	250 mm
seripettor® pro 2 ml, 10 ml	GL 32-33, GL 38, S 40	125 - 240 mm
seripettor® pro 25 ml	GL 32-33, GL 38, S 40	170 - 330 mm

## 1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: [www.brand.de](http://www.brand.de).

### 1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

### 1.2.2 Symbols

Symbol	Meaning
	Hazardous area

## 1.2.3 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

## 2 Safety Instructions

### 2.1 General safety instructions

#### Please read carefully!

The instrument seripettor® can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations (e.g., wear protective clothing, eye protection, and protective gloves).
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Comply with the operating exclusions(see Operating Exclusions, p. 52)! If in doubt, contact the manufacturer or supplier.
6. Always perform work in a manner that does not endanger yourself or other people. When dispensing, the discharge tube must always point away from you or any other person. Avoid splattering. Only use suitable vessels.
7. Never press down the piston when the discharge tube is sealed with the closure cap/screw cap.
8. During operation, the piston locking mechanism must be closed.
9. Never remove the discharge tube while the dispensing cylinder is filled.
10. Reagents can accumulate in the closure cap/screw cap of the discharge tube. Therefore, clean regularly.
11. For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
12. Never carry the mounted instrument by the actuating unit or the valve block, see First Steps, p. 57. Breakage or loosening of instrument parts can lead to personal injury from chemicals.
13. Never use force. Always press the piston down gently when dispensing.
14. Use only original accessories and original replacement parts. Do not make any technical modifications. Do not dismantle the instrument any further than is described in the operating manual!
15. Always check that the instrument is in proper working condition before use. The user can come into contact with media if the instrument has been insufficiently cleaned or inspected. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult Troubleshooting. Contact the manufacturer, if necessary.

### 2.2 Function

The seripettor® bottle-top dispenser is used for dispensing liquids directly from the storage bottle and is available in two versions:

seripettor®

seripettor® pro

## 2.2.1 Operation

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

### seripettor®

FEP, PP, PE, EPDM

### seripettor® pro

PP, PE, borosilicate glass, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, PFA, FEP, PTFE, Pt-Ir

## 2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- Use between +15°C and +40°C of instrument and reagent (seripettor®: agar culture medium up to max. 60 °C)
- Vapor pressure up to 500 mbar
- Density up to 2.2 g/cm<sup>3</sup>
- Kinematic viscosity\*:
  - 2 ml instrument: 300 mm<sup>2</sup>/s
  - 10 ml instrument: 150 mm<sup>2</sup>/s
  - 25 ml instrument: 75 mm<sup>2</sup>/s

\* dynamic viscosity [mPas] = kinematic viscosity [mm<sup>2</sup>/s] x density [g/cm<sup>3</sup>]

## 2.4 Operating Limitations

- Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions). If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately (Cleaning, p. 68).
- When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
- The instrument is designed for general laboratory applications. Compatibility of the instrument for a specific application (e.g., trace material analysis, food sector etc.) must be checked by the

user. Approvals for specific applications, e.g. for production and administration of food, pharmaceuticals or cosmetics are not available.

## 2.5 Operating Exclusions

### 2.5.1 seripettor®

Never use seripettor® for:

- Liquids attacking FEP, PP, PE or EPDM
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Non-polar solvents such as hydrocarbons and halogen hydrocarbons
- Concentrated or oxidizing acids
- Explosive liquids
- Carbon disulfide

### 2.5.2 seripettor® pro

Never use seripettor® pro for:

- Liquids attacking PP, PE, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide\*)
- Liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Liquids which are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>)
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Non-polar solvents such as hydrocarbons and halogen hydrocarbons
- Concentrated or oxidizing acids (except for HCl)
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated carbon solutions) because solid particles may clog or damage the instrument

\* Dissolved sodium azide permitted up to a concentration of max. 0.1%.

## 2.6 Storage Conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

## 2.7 Recommended Application Range

For guidelines on selecting the right dispenser, please observe the corresponding Operating Exclusions and the 'Dispenser selection chart' on the next page.

### 2.7.1 seripettor® application range

#### Aqueous solutions

Biological buffer solutions and detergents, anti-foaming agents, culture media, vitamin solutions, etc. as well as hydrogen peroxide can be dispensed in everyday routines.

Agar culture media can be dispensed up to max. 60 °C.

## Acids

Low-concentrated or diluted, non-oxidizing acids can be filled.

## Alkalis

Can be used for dispensing alkaline media such as NaOH, KOH and ammonia.

## Polar solvents

E.g., ethanol, methanol, acetylacetone, etc.

## 2.7.2 seripettor® pro application range

The seripettor® pro bottle-top dispenser expands the range of applications. It enables dispensing of acids (e.g., concentrated HCl), polar solvents (e.g., acetone) and UV-sensitive media

## 2.8 Dispenser Selection Chart

### NOTICE

**Not suitable for hydrofluoric acid (HF)**  
seripettor® and seripettor® pro are not suitable for dosing hydrofluoric acid (HF).

	seripettor®	seripettor® pro
1-butanol	—	✓
Acetaldehyde	—	✓
Acetic acid, 100% (glacial acetic acid)	—	✓
Acetic acid, 5%	✓	✓
Acetic acid, 96%	—	✓
Acetone	—	✓
Acetonitrile	—	✓
Acetophenone	✓	—
Acetylacetone	✓	✓
Adipic acid	✓	✓
Agar (60°C)	✓	—
Allyl alcohol	✓	✓
Aluminum chloride	✓	✓
Amino acids	✓	✓
Ammonia, 30%	✓	✓
Ammonium chloride	✓	✓
Ammonium fluoride	✓	✓
Ammonium sulfate	✓	✓
Amyl alcohol (pentanol)	✓	✓
Aniline	—	✓
Barium chloride	✓	✓
Benzaldehyde	—	✓
Benzoic acid methyl ester	—	✓
Benzyl alcohol	—	✓
Benzyl chloride	—	✓
Benzylamine	—	✓
Boric acid, 10%	✓	✓
BSA serum	✓	✓
Butanediol	✓	✓
Butylamine	—	✓
Calcium carbonate	✓	✓
Calcium chloride	✓	✓
Calcium hydroxide	✓	✓
Calcium hypochlorite	—	✓
Chloroacetaldehyde, 45	—	✓
Chloroacetic acid	—	✓
Chromic acid, 50%	—	✓
Copper sulfate	✓	✓
Cumene (isopropylbenzene)	—	✓
Diethylene glycol	✓	✓
Dimethyl sulfoxide (DMSO)	—	✓
Dimethylaniline	—	✓
Ethanol	✓	✓
Ethidium bromide solution	✓	✓
Ethyl methyl ketone	—	✓

	seripettor®	seripettor® pro
Formaldehyde, 40%	✓	✓
Formamide	✓	✓
Formic acid, 100%	—	✓
Glycerin	✓	✓
Glycol (ethylene glycol)	✓	✓
Glycolic acid, 50%	✓	✓
Guanidine hydrochloride	✓	✓
HEPES buffer	✓	✓
Hexanoic acid	✓	✓
Hexanol	—	✓
Hydrobromic acid	—	✓
Hydrochloric acid, 37	—	✓
Hydrogen peroxide, 35	✓	—
Hydroiodic acid	✓	✓
Isoamyl alcohol	—	✓
Isobutanol (isobutyl alcohol)	✓	✓
Isopropanol (2-propanol)	✓	✓
Lactic acid	✓	✓
LB media	✓	✓
McCoy's 5A	✓	✓
MEM	✓	✓
Methanol	✓	✓
Methyl propyl ketone	—	✓
Mineral oil (engine oil)	—	✓
Monochloroacetic acid	—	✓
n-amyl acetate	—	✓
n-Butyl acetate	—	✓
Nitric acid, 10%	—	✓
Octoxinol 9 (TRITON™ x-100)	✓	✓
Oxalic acid	✓	✓
PBS buffer	✓	✓
Perchloric acid	—	✓
Phenol	—	✓
Phosphoric acid, 85	—	✓
Piperidine	—	✓
Polysorbate (TWEEN®)	✓	✓
Potassium chloride	✓	✓
Potassium dichromate	✓	✓
Potassium hydroxide	✓	✓
Potassium hydroxide in ethanol	✓	✓
Potassium permanganate	✓	✓
Propionic acid	✓	✓
Propylene glycol (propane-diol)	✓	✓
Pyridine	—	✓
Pyruvic acid	✓	✓
Ringer's solution	✓	✓
RPMI 1640	✓	✓
Salicylaldehyde	—	✓
Salicylic acid	✓	✓
SDS (sodium dodecyl sulfate)	✓	✓
Silver acetate	✓	✓
Silver nitrate	✓	✓
Sodium acetate	✓	✓
Sodium chloride	✓	✓
Sodium dichromate	✓	✓
Sodium fluoride	✓	✓

	<b>seripettor®</b>	<b>seripettor® pro</b>
Sodium hydroxide, 30%	✓	✓
Sodium hypochlorite, 20% (active chlorine approx. 10%)	—	✓
Sulfuric acid, 10%	✓	✓
Tartaric acid	—	✓
TE buffer	✓	✓
TRIS buffer	✓	✓
Tris-buffered saline with Tween20	✓	✓
Urea	✓	✓
Zinc chloride, 10%	✓	✓
Zinc sulfate, 10%	✓	✓

## Legend:

✓ = The device is suitable for the medium


— = The device is not suitable for the medium

This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. In addition to the chemicals listed above, a number of organic or inorganic saline solutions (e.g. biological buffers), biological detergents as well as media for cell cultures can be dispensed. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

Last updated: 1023/11

# 3 Functional and operating elements

seripettor® and seripettor® pro

- 1 Actuating unit
- 2 Piston locking mechanism
- 3 Volume adjustment
- 4 Dispensing cartridge (piston/cylinder)
- 5 Valve block
- 6 Valve block adapter (GL 45 bottle thread)
- 7 Filling tube
- 8 Discharge tube with integrated discharge valve
- 9 Adapter for discharge tube
- 10 Screw cap
- 11 Closure caps
- 12 Discharge tube

## 4 Assembly

### 4.1 First Steps


#### ⚠ WARNING

#### Follow during every use, especially with hazardous media

- Wear protective clothing, eye protection and appropriate hand protection!
- Never press down the piston when the screw cap is screwed on!
- Avoid splashing the reagent!
- To avoid splashes dispense slowly.
- Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use and Operating exclusions, p. 52.

#### 1. Screwing on the dispensing cartridge

- a. Screw dispensing cartridge into valve block.


#### 2. Installing the actuating unit

- a. Open piston lock on the actuating unit.

- b. Put actuating unit over the dispensing cartridge and screw it down to the valve block.



- c. Depress actuating unit up to the lower stop and keep holding it using one hand and with the other hand, close piston locking on the actuating unit.
- d. Release actuating unit and check to see whether the piston in the dispensing cartridge has moved upward.

**NOTICE**


Check that the discharge tube is securely seated. Re-tighten cap nut(s) after two days of use.

### 3. Installing the filling tube

- a. seripettor®: Cut the filling tube off according to the bottle height and attach it.
- a'. seripettor® pro: Adjust the length of the telescopic filling tube to the bottle height and attach it.

## 4. Installing the instrument on a bottle

- a. Screw the instrument (GL 45 threads) onto the reagent bottle.
- b. To avoid tipping over, use a bottle stand for small bottles.

### NOTICE


For bottles with other thread sizes, select a suitable adapter, see Accessories/spare parts, p. 79

## 5. Transporting the instrument

- a. When installed on a reagent bottle, always carry the instrument as shown in the figure and always store it upright!

## 6. Priming the instrument

- a.** Hold the discharge tube and remove the closure cap/screw cap.

- b.** To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel. For priming, let the piston move up approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until there are no more air bubbles under the piston. Bubbles less than 1 mm are permitted.

### NOTICE

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splattering.

# 5 Operation

## 5.1 Dispensing

### ⚠ WARNING


#### Follow during every use, especially with hazardous media

- Wear protective clothing, eye protection and appropriate hand protection!
- Never press down the piston when the screw cap is screwed on!
- Avoid splashing the reagent!
- To avoid splashes dispense slowly.
- Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use and Operating exclusions, p. 52.

### 1. Selecting the volume

#### NOTICE

Liquid is dispensed when setting the volume!

- a. Hold the discharge tube and remove the closure cap/screw cap.
- b. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.
- c. With one hand, press the actuating unit downward and hold it.
- d. Place the receiving vessel under the discharge tube orifice.
- e. Loosen the volume selector thumb screw by a half turn, slide the pointer vertically to the desired volume and re-tighten the volume thumb screw.

## 2. Dispensing

- a. Hold the discharge tube and remove the closure cap/screw cap.
- b. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- c. Push the piston down to the lower stop evenly and with minimal effort and then let it move up slowly.
- d. Wipe off the discharge tube against the inner wall of the receiving vessel.
- e. Seal the dispensing tube with the closure cap/screw cap.

## 5.2 Dispensing with the flexible discharge tube

The flexible discharge tube can optionally be used for serial dispensing (see Accessories/spare parts, p. 79). The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. The tubing should lie neatly in loops and must not be twisted.


The parts which come into contact with media consist of borosilicate glass, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, PTFE, platinum-iridium and PP.

Never use the flexible discharge tube for:


- Liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Peroxides, as these are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>)

The operating exclusions for the respective instrument in use also apply.


## Installing the flexible discharge tube

- a.** Screw the tube holder to the valve block and install the collection tube.  
 seripettor®: Before installing the flexible discharge tube, remove the built-in elastomer sealing ring and replace it with the provided PTFE sealing ring.

- b.** Check that the flat seal is inserted in the adapter.

- c.** Slide the discharge tube with flexible tubing onto the adapter and tighten the cap nut by hand. Then slide the adapter onto the valve block and again tighten the cap nut by hand. Use a bottle holder.

## 5.3 Dispensing with a microfilter

### Discharge tube with Luer-Lock attachment for microfilters

The discharge tube with Luer-Lock attachment for microfilters can optionally be used for sterile filtration of media, see Accessories/spare parts, p. 79.

The discharge tube with Luer-Lock attachment allows a microfilter to be connected for sterile filtration. The parts which come into contact with media consist of: borosilicate glass, I<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium and PP (Luer-Lock attachment).

Therefore, never use for:

- Liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Peroxides, as these are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>)

The operating exclusions for the instrument and microfilter in use also apply. Only use microfilters with Luer-Lock. Never exceed the maximum operating force of 100 N.

### Assembly

- Dismantle the standard discharge tube from the seripetto® pro and slide the discharge tube with Luer-Lock approx. 2 mm onto the adapter for discharge tubes.
- Then firmly tighten the cap nut by hand. Commercially available sterile filters can be installed onto the Luer-Lock attachment. After connecting the filter, check that it is firmly seated.

#### NOTICE


Observe the regulations for handling sterile media. The increased flow resistance can lead to leaking liquid in the dispensing cartridge. To keep any liquid leaks to a minimum, we recommend dispensing with little force and using a filter with a large filter surface. Please follow the recommendations of the filter manufacturer in terms of media compatibility. Watch out for any increase in pressure during use, and change the filter in a timely manner before it becomes blocked.

## 5.4 Bottle stand

For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over (Accessories/spare parts).

### Assembly

- a. Position the mounting plate at the appropriate height.
- b. As shown, insert the instrument firmly into the holder until it audibly clicks into place.
- c. Then lock the holder into place with the screw.

## 6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. Testing takes place according DIN EN ISO 8655-6 with a completely filled instrument and with uniform and smooth dispensing.

### Error limits

Nominal volume ml	A* ± %	A* ± µl	CV* %	CV* µl
2	1.2	24	0.2	4
10	1.2	120	0.2	20
25	1.2	300	0.2	50

\*A = Accuracy, CV = Coefficient of Variation

### Partial volume

The percentage values for A and CV are relative to the nominal volume ( $V_N$ ) and must be converted for partial volumes ( $V_T$ ).

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

e.g.	Volume	A* ± %	A* ± µl	CV* %	CV* µl
$V_N$	10	1.2	120	0.2	20
$V_T = 50\% N$	5	2.4	120	0.4	20
$V_T = 10\% N$	1	12	120	2.0	20

\*A = Accuracy, CV = Coefficient of Variation

### NOTICE

The maximum error for a single measurement can be calculated from the sum of error limits  $EL = A + 2 \times CV$  (e.g. for 10 ml:  $120 \mu\text{l} + 2 \times 20 \mu\text{l} = 160 \mu\text{l}$ ).

## 7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3–12 months. This time frame should be adjusted to correspond with individual requirements.

Gravimetric volume testing according to DIN EN ISO 8655-6 (for measurement conditions, see Error limits, p. 66) is performed as follows:

### 1. Prepare the instrument

Clean the instrument (Cleaning, p. 68), fill with distilled water, and prime carefully.

### 2. Check the volume

- a. 10 dispensing operations with distilled H<sub>2</sub>O in 3 volume ranges (100%, 50%, 10%) are recommended
- b. For discharge, depress the piston slowly and steadily without force until the lower stop
- c. Wipe off the tip of the dispensing tube.
- d. Weigh the dispensed quantity on an analytical balance. (Please refer to the user manual of the scale manufacturer.)
- e. Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

### Calculation (for nominal volume)

$x_i$  = weighing results

$n$  = number of weighings

$V_0$  = nominal volume

$Z$  = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy\*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation\*:

$$CV\% = \frac{100 s}{\bar{V}}$$

Standard deviation\*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

\*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

#### NOTICE

Test Instructions (SOPs) are available for download from [www.brand.de](http://www.brand.de).

# 8 Cleaning

## ⚠ WARNING

### Components filled with reagent

The cylinder, valves, telescopic filling tube and discharge tube contain reagent!

- > Never remove the discharge tube while the dispensing cylinder is filled.
- > Point the valves and tube openings away from your body.
- > Wear protective clothing, eye protection and protective gloves!

The instrument must be cleaned in the following situations to assure correct operation:

- immediately when the piston is difficult to move
- prior to a reagent change
- prior to long term storage
- prior to dismantling the instrument
- prior to changing the valve
- prior to autoclaving
- regularly when using liquids which form deposits (e.g., crystallizing liquids)
- regularly when liquids accumulate in the screw cap.

## 8.1 Basic cleaning

### 1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

### 2. Rinsing the instrument


- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.


## 8.2 Cleaning and replacing valves

### 8.2.1 seripettor®

#### Filling valve


- a. Remove filling tube.
- b. Unscrew the filling valve with a coin and extract the filling valve body.

- c. If the sealing ring is contaminated or damaged, carefully remove it using a pair of angled tweezers.
- d. Clean any contaminated parts (e.g., in an ultrasonic bath).
- e. Insert cleaned or new sealing ring.


**f.**

- f.** Screw in the cleaned or new filling valve first by hand and then tighten it with a coin.

## Discharge valve


**e.**

- a.** Loosen the screw cap of the dispensing tube.
- b.** Extract the discharge valve.
- c.** If the sealing ring is contaminated or damaged, carefully remove it using a pair of angled tweezers.
- d.** Clean any contaminated parts (e.g., in an ultrasonic bath).
- e.** Install the dispensing tube with the cleaned or new discharge valve.

## 8.2.2 seripettor® pro


### Filling valve

- a. Remove filling tube.
- b. Use the assembly tool to unscrew the filling valve.


- c. If the sealing ring is contaminated or damaged, carefully remove it using a pair of angled tweezers.
- d. Clean any contaminated parts (e.g., in an ultrasonic bath).
- e. Insert cleaned or new sealing ring.

- f. Screw in the cleaned or new filling valve first by hand and then tighten it with the assembly tool.

## Discharge valve

The discharge valve is built into dispensing tube.

- a. Dismantle the dispensing tube and clean in an ultrasonic bath if necessary.
- b. Install cleaned or new dispensing tube.

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200  $\mu$ l plastic pipette tip (see the figure at the side).

## 8.3 Dislodging a stick valve ball


If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.


In this case, dislodge the valve ball using light pressure; for example, with a 200  $\mu$ l plastic pipette tip (see the figure at the side).

## 8.4 Replacing the dispensing cartridge


The dispensing cartridge is a wear part. In case of sluggishness or leakage between cylinder and piston, the dispensing cartridge must be replaced. To prevent injuries due to chemical exposure, cleaning must be performed before replacing the dispensing cartridge (see Cleaning, p. 68).


- a. Screw the instrument onto an empty bottle and empty it completely by dispensing multiple times.
- b. With one hand, press the actuating unit downward all the way and hold it. Loosen the volume selector thumb screw by a half turn, slide the pointer vertically to the upper stop and re-tighten the volume thumb screw.

**c.** Release the piston lock.

**d.** Unscrew and remove the actuating unit.

- e.** Unscrew the dispensing cartridge and replace it with a new one, see Accessories/spare parts, p. 79.
- f.** Screw the actuating unit on and engage the piston lock.

## 9 Autoclaving

### NOTICE

#### The following are not autoclavable:

Dosing unit and actuator unit

The valve block, aspiration tube, dosing cannula and valves are autoclavable at 121 °C at 2 bar and a holding time of at least 15 minutes according to DIN EN 285.


- Unscrew actuating unit and dispensing cartridge, see Replacing the dispensing cartridge, p. 72.
- Loosely seal the valve block with the closure cap (see Accessories/spare parts, p. 79).
- Insert the autoclavable filling tube (see Accessories/spare parts, p. 79 with the sealing ring up into filling valve.
- Loosely screw the sealed valve block with filling tube and dispensing tube onto the filled bottle.

e.

- Seal the ventilation opening with an autoclavable membrane filter (0.2 µm).

f.

- Open dispensing tube, loosen cap nut.
- Autoclave.


### NOTICE

#### Allow parts to cool slowly after autoclaving

Allow to cool in the closed autoclave to avoid deformation due to too fast cooling (cooling time approx. 2 hours, exception: seripettor® with agar culture medium up to max. 60 °C). After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them.

It is the user's responsibility to ensure effective autoclaving.

## 9.1 Preparing for sterile work

- a. After autoclaving, seal the dispensing tube and tighten the screws.
- b. Under sterile conditions, remove the closure cap.
- c. Screw in the sterile dispensing cartridge (see Accessories/spare parts, p. 79).
- d. Fasten the actuating unit.

# 10 Troubleshooting

Problem	Possible cause	Corrective action
Piston sluggish	Crystal deposits	Replace dispensing cartridge, see Replacing the dispensing cartridge, p. 72).
	Forbidden liquid dispensed	Clarify whether the medium is permitted to be dispensed, see Dispenser Selection Chart.
Liquid is above the dosing piston	Dispensing cartridge leaky	Replace dispensing cartridge, see Replacing the dispensing cartridge, p. 72.
Filling not possible	Valve incorrectly installed	Correctly install valve, see Cleaning and replacing valves, p. 69.
	Volume adjusted to minimum setting	Set the desired volume, see Dispensing, p. 61.
	Filling valve stuck	Unscrew the filling valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Dislodging a stick valve ball, p. 72). If necessary, replace the filling valve.
	Dispensing cartridge leaky	Replace dispensing cartridge, see Replacing the dispensing cartridge, p. 72).
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block, clean it and replace if necessary. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it.
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Valve screw connections loose	Tighten the valves firmly with the assembly tool.
	The instrument has not been primed	Prime the instrument, see First Steps, p. 57
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Perform cleaning, see Cleaning, p. 68. Tighten the valves using the assembly tool.
Dispensed volume too low	Discharge tube is loose or damaged	Correctly install the discharge tube. Replace deformed or damaged discharge tube.
	Filling tube is loose or damaged	Perform cleaning, see Cleaning, p. 68. Push the filling tube on firmly. If that is insufficient, cut off approx. 1 cm of tube at the upper end or replace the filling tube.
	Filling valve is loose or damaged	Perform cleaning, see Cleaning, p. 68. Tighten the filling valve. If necessary, replace the valve with sealing ring, see Cleaning and replacing valves, p. 69.
	Frequent dispensing > 40 °C	Replace dispensing cartridge, see Replacing the dispensing cartridge, p. 72).

# 11 Product markings

Symbol or number	Meaning
	General warning symbol
	Refer to the operating manual
	Use eye protection
	Use hand protection
	Wear protective clothing
XXZXXXXX	Serial number

## 12 Ordering Information

### seripettor®

Volume ml	Graduation ml	A* $\leq \pm \%$	$\mu\text{l}$	CV* $\leq \%$	$\mu\text{l}$	Order no.
0.2-2	0.04	1.2	24	0.2	4	<a href="#">4720120</a>
1 - 10	0.2	1.2	120	0.2	20	<a href="#">4720140</a>
2.5-25	0.5	1.2	300	0.2	50	<a href="#">4720150</a>

### seripettor® pro

Volume ml	Graduation ml	A* $\leq \pm \%$	$\mu\text{l}$	CV* $\leq \%$	$\mu\text{l}$	Order no.
0.2-2	0.04	1.2	24	0.2	4	<a href="#">4720420</a>
1 - 10	0.2	1.2	120	0.2	20	<a href="#">4720440</a>
2.5-25	0.5	1.2	300	0.2	50	<a href="#">4720450</a>

\*Final test values in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20°C) of instrument, ambient environment and distilled water, and with uniform, jolt-free handling. A = Accuracy, VC = Variation Coefficient.

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\* buttress thread

\*\* GL, S: the numbers correspond to the outer diameter of the bottle thread in mm.

#### NOTICE

##### Template for adapter

An adapter template is available for download at [www.brand.de](http://www.brand.de) in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop [store.brand.de](http://store.brand.de) in the “Selection Charts” section.

### Dispensing cartridges

#### NOTICE

The dispensing cartridges are not autoclavable.

Dispensing cartridges for seripettor® and seripettor® pro. Packaged individually.

Model	Material	Pkg unit	Order no.
2 ml	Piston (PE), cylinder (PP).	3	<a href="#">704500</a>
10 ml	Piston (PE), cylinder (PP).	3	<a href="#">704502</a>
25 ml	Piston (PE), cylinder (PP).	3	<a href="#">704504</a>
2 ml, sterile***	Piston (PE), cylinder (PP).	7	<a href="#">704507</a>
10 ml, sterile***	Piston (PE), cylinder (PP).	7	<a href="#">704506</a>
25 ml, sterile***	Piston (PE), cylinder (PP).	5	<a href="#">704508</a>

\*\*\*packaged individually

### Discharge tube for seripettor®

Discharge tube for seripettor® with closure cap, EPDM discharge valve and sealing ring.

Model	Pkg unit	Order no.
2 ml, fine tip	1	<a href="#">704518</a>
10 + 25 ml, standard	1	<a href="#">704520</a>

### Discharge tube for seripettor® pro


Discharge tube for seripettor® pro with integrated discharge valve and flat seal.

Model	Pkg unit	Order no.
2 ml	1	<a href="#">707915</a>
10 ml	1	<a href="#">707916</a>
25 ml	1	<a href="#">707918</a>

Discharge tube adapter sold separately.

## Discharge tube with Luer-Lock attachment for microfilter

Model*	Material	Pkg unit	Order no.
2 ml, 10 ml, 25 ml	FEP/PP	1	<a href="#">707928*</a>

\* Not suitable for HF and peroxide

Discharge tube adapter sold separately.

## Valve set for seripettor®

Description	Order no.
1 discharge valve with sealing ring, 1 filling valve with sealing ring and filling valve body with seal.	<a href="#">6790</a>

## Discharge valve for seripettor® pro

Description	Model	Pkg unit	Order no.
Filling valve with sealing ring	2 + 10 ml	1	<a href="#">6697</a>
	25 ml	1	<a href="#">6698</a>

Filling valve adapter sold separately.

## Sealing rings for seripettor®

Description	Material	Pkg unit	Order no.
Sealing ring	EPDM	5 pc per set	<a href="#">6788</a>


## Flexible discharge tubing for seripettor® and seripettor® pro

Description	Model	Pkg unit	Order no.
PTFE, spiraled tube, approx. 800 mm long, with safety handle.	2 + 10 ml	1	<a href="#">704522</a>
	25 ml	1	<a href="#">704523</a>


\* Not suitable for peroxides

## Membrane filter


Description	Material	Pkg unit	Order no.
PTFE membrane filter (0.2 µm), autoclavable (121°C).	PTFE	10 pc in PE bag	<a href="#">26535</a>

## Adapter for seripettor® pro


Description	Material	Pkg unit	Order no.
Discharge tube adapter (1) for seripettor® pro, ETFE	ETFE	1	<a href="#">6208</a>
Filling valve adapter (2) for seripettor® pro, ETFE	ETFE	1	<a href="#">6707</a>

## Filling valve for seripettor®

Description	Material	Length (mm)	Pkg unit	Order no.
Autoclavable model with additional sealing ring	PP	250	2	<a href="#">704532</a>
	PP	500	2	<a href="#">704534</a>
	PP	250, with sealing ring	1	<a href="#">704536</a>
	PP	500, with sealing ring	1	<a href="#">704538</a>

## Telescoping filling tube for seripettor® pro

Material	Model	Length (mm)	Pkg unit	Order no.
FEP	2 + 10 ml*	70 - 140	1	<a href="#">708210</a>
		125 - 240	1	<a href="#">708212</a>
		195 - 350	1	<a href="#">708214</a>
		250 - 480	1	<a href="#">708216</a>
	25 ml**	170 - 330	1	<a href="#">708218</a>
		250 - 480	1	<a href="#">708220</a>

\* Outer-Ø: 6 mm \*\* Outer-Ø: 7.6 mm

## Actuating unit for seripettor®

Description	Material	Model	Pkg unit	Order no.
Stainless steel lifting spring	PC	2 ml	1	<a href="#">704541</a>
	PC	10 ml	1	<a href="#">704542</a>
	PC	25 ml	1	<a href="#">704544</a>

## Actuating unit for seripettor® pro

Description	Material	Model	Pkg unit	Order no.
Stainless steel lifting spring made of Hastelloy® (rust-free).	PPO. PEI (UV protection)	2 ml	1	<a href="#">704551</a>
	PPO. PEI (UV protection)	10 ml	1	<a href="#">704548</a>
	PPO. PEI (UV protection)	25 ml	1	<a href="#">704549</a>

## Valve block

Description	Model	Pkg unit	Order no.
Valve block for discharge tubing for seripettor® and seripettor® pro	2 + 10 ml	1	<a href="#">6792</a>
	25 ml	1	<a href="#">6794</a>

## Bottle stand

Description	Pkg unit	Order no.
Bottle stand PP. Support rod, 325 mm, base plate 220 x 160 mm.	Pack of 1	<a href="#">704275</a>

## Sealing cap for valve block

Description	Model	Pkg unit	Order no.
Sealing cap for valve block for seripettor® and seripettor® pro	2 + 10 ml	1	<a href="#">704552</a>
	25 ml	1	<a href="#">704554</a>

# 14 Repairs

## 14.1 Sending for repair

### NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

### Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

### Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

### Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

### Contact addresses

#### Germany:

BRAND GMBH + CO KG  
Otto-Schott-Straße 25  
97877 Wertheim (Germany)  
T +49 9342 808 0  
F +49 9342 808 98000  
[info@brand.de](mailto:info@brand.de)  
[www.brand.de](http://www.brand.de)

#### USA and Canada:

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)  
T +1-860-767 2562  
F +1 - 860 - 767 2563  
[info@brandtech.com](mailto:info@brandtech.com)  
[www.brandtech.com](http://www.brandtech.com)

#### India:

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791

#### China:

BRAND (Shanghai) Trading Co., Ltd.  
Rm 201-202, North Tower,  
No. 199 Kaibin Rd, Xuhui District, Shanghai  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
[info@brand.com.cn](mailto:info@brand.com.cn)

info@brand.co.in  
www.brand.co.in

www.brand.cn.com

## 15 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing instructions are available for download on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or our accredited calibration laboratory. Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days. The instruments are supplied with a detailed calibration certificate or a calibration slip in accordance with DIN EN ISO/IEC 17025. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on [www.brand.de](http://www.brand.de) (see Service and Support).

### For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

## 16 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

### **USA and Canada:**

Find more warranty information on [www.brandtech.com](http://www.brandtech.com).

## 17 Disposal

Before disposal, observe the relevant national disposal regulations, and ensure that the product is disposed of properly.

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