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PLT unit

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F I R S T C L A S S · B R A N D

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35 Operating Manual

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	Page
Safety Instructions	36
CE Marking	37
Functions and Limitations of use	37
Operating exclusions	37
Operating elements	38
First Steps	39
Initializing	40
Perform self-test	40
Quick test	40
Full test	42
Test preparation	44
Dynamic or static test?	44
Test with or without pipette tip?	44
Single-channel or multi-channel adapter selection	45
Changing the adapter and filter	46
Testing single-channel pipettes	47
Testing multi-channel pipettes	50
Switch the instrument off	53
Adjustments	54
Selecting the language	54
Adjusting the display contrast	55
Setting the inspection date	56
Setting the Auto-Power-Off Time	57
Selecting the pressure units (hPa, mbar, torr)	58
Restoring the factory settings	59
Firmware update via the USB interface	60
Leak rate/volume correlation table	61
Ordering Information · Accessories · Spare Parts	62
Troubleshooting	63
Repair and Warranty Information	65
Disposal	66

Safety Instructions

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

Please read the following carefully!

1. Every user must read and understand this operating manual before operation.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Never use the instrument in an atmosphere with a danger of explosion.
4. The instrument should only be used for the leak testing of pipettes that operate on the air displacement principle. Observe operating exclusions (s. page 37)! If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered.
6. When testing air displacement pipettes with mounted tips, employ only unused tips.
7. Never use force on the instrument!
8. Use only original manufacturer's accessories and spare parts. Do not attempt to make any technical alterations.
9. Only authorized service personnel may repair or service the instrument.
10. Always check the instrument for visible damage before use. If there is a sign of a potential malfunction, immediately stop the test. Consult the 'Troubleshooting' section of this manual (see page 63), and contact the manufacturer if needed.
11. The AC adapter has to be protected against moisture and must be used only for this instrument.

Warning!

Improper use of the instrument or the battery (short-circuit, mechanical damage, overheating etc.) may cause the explosion of the battery.

The BRAND leak testing instrument is a differential pressure measuring instrument for determining the leak rate of an air-displacement pipette. Static and dynamic testing permit the localization of leaks. Dynamic testing is preferable.

Testing is possible both with tips (for testing the instrument/tip interface and the entire system) and without tips (for testing the instrument). The leak rate is directly related to the volume loss of a pipette

The leak rate QL is a measure of the quantity of material (mass) that flows through a leak per unit time.

Note:

The leak tester cannot substitute for gravimetric testing for the monitoring of measuring instruments. The instrument should be used to check pipettes on a daily basis as a safeguard for the periods between calibrations.

Limitations of Use

The instrument offers a rapid and simple test of leakage in air-displacement pipettes within the following limitations:

- Operating temperature 15 °C to 35 °C
- Rel. Humidity max 80%
- For use at elevations up to 2000 m above sea level
- Storage temperature 0 °C to 50 °C

Important!

When the instrument is relocated to a new environment or even during transport, significant changes in the temperature and atmospheric humidity can affect the measurements.

To avoid error messages, the PLT unit must be allowed to acclimate in the testing laboratory for a period of 2 hours before use!

The leak rate is the ratio of the pV value of a gas and the period of time during which the gas flows through the leak. The pV value is the product of the pressure and the volume of a certain quantity of a gas at the prevailing temperature. In ideal gases, the pV value at a given temperature is a measure of the quantity of material or the mass of the gas. The leak rate depends on the type of gas, absolute pressure, pressure difference, and temperature.

For the pipette test, hPa ml/s is a suitable unit for the leak rate QL. A leak rate of 1 hPa ml/s at an air pressure of 1000 hPa means a volume loss of 1 µl/s.

CE Marking

This sign certifies that the product meets the requirements of the EEC directive and has been tested according to the specified test methods.

Operating Exclusions

Users themselves must check the suitability of the instrument for the intended application.

Do not use the instrument to aspirate liquids! Avoid aggressive vapors (corrosion hazard)! Further information available upon request.

Hollow vessels that are not specified for vacuum use should never be evacuated, since the vessel can burst!

Battery Specifications

PLT unit:

Input: DC 6.5 V = 800 mA

AC adapter:

Input: AC 100 V - 240 V, 50/60 Hz, 0,5 A

Output: DC 6,5 V = 800 mA

Protection class: IP20

Operating Elements

According to ISO 8655, air-displacement pipettes should be checked regularly. The BRAND PLT unit offers a safeguard for the periods between calibration intervals.

Back of the instrument

AC adapter socket

USB port

"Back" button

The function of the "Back" button is varied by pressing for a short or long time. Pressing the button for a short time leads back to a sub-menu in the next higher level of the menu hierarchy. Pressing the button for a long time leads back directly to "Ready" on the Start Menu.

Select knob

In many cases, pressing the left "Start/Execute" button as well as the central knob can be used to confirm a selection. For consistency, the left "Start/Execute" button is used for confirmation in this Operating Manual.

Powering on the instrument

Is everything in the package?

The package contains the PLT unit, one adapter each for testing air-displacement pipettes with tips (mounted) or without tips, 2 blind plugs, 3 replacement PE filters for the pipette adapter, an AC adapter, a quality certificate, and this Operating Manual.

What to do?	How to do it	Keys to press	Display readout
-------------	--------------	---------------	-----------------

1. Connecting the AC adapter

Use only the original AC adapter!


The instrument's AC adapter and AC adapter plug must be freely accessible, and separatable from the power grid at any time.

Connect the AC adapter to the socket in the back of the instrument, without pressing any buttons on the instrument.

After the AC adapter is connected, the display shows the software version number, the creation date, and the BRAND logo.

During initial start-up, the display brings up the "Language" indicator after 3 sec.

2. Selecting the language


The language is selected by turning the select knob.

3. Confirming the language

Press the left "OK" button. The "Settings" menu appears in the language selected.

The "Back" button leads to the Main Menu.

Note:

After start-up, one can select either to perform the self-test (page 40) or directly begin a pipette test (page 44).

Perform self-test

The self-test determines if there are any leaks in the BRAND PLT unit. Two different tests are available. The Quick Test lasts only about 2 minutes, and should be performed on a daily basis. The Full Test lasts about 12 minutes. It includes the additional differentiated testing of the internal tanks, and should be performed on a weekly basis.


Quick test

What to do?	How to do it	Keys to press	Display readout
1. Calling up the Main Menu	The Main Menu is called up by pressing or turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Main Menu</p> <p>Pipette Type</p> <p>Self Test</p> <p>Switch OFF</p> <p>Settings...</p> <hr/> <p>select back</p> </div>
2. Setting up the self-test	Set to "Self Test" by turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Main Menu</p> <p>Pipette Type</p> <p>Self Test</p> <p>Switch OFF</p> <p>Settings...</p> <hr/> <p>select back</p> </div>
3. Confirming the self test	Press the "Execute" button. The "Self Test" submenu is called up. The last test to be performed appears enframed.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Self Test</p> <p>Quick Test</p> <p>Full Test</p> <hr/> <p>select back</p> </div>
4. Selecting the Quick Test	Set to "Quick Test" by turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Self Test</p> <p>Quick Test</p> <p>Full Test</p> <hr/> <p>select back</p> </div>
5. Confirming the Quick Test	Press the "Start/Execute" button. The "Start Menu" submenu is called up.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Self Test</p> <p>Please insert blind plug</p> <hr/> <p>start back</p> </div>
6. Seal the system	Press the provided blind plug firmly into the adapter (for the multi-channel adapter, seal all 4 channels).		

What to do?	How to do it	Keys to press	Display readout
-------------	--------------	---------------	-----------------

7. Starting the Quick Test


Press the "Start/Execute" button. The five components - sensor, pump, valve, filter, and measurement chamber - are tested one after the next. The progress bar in the display indicates the test status. The LED indicator alternately blinks red and green during the test.


(pressure value is only an example)

8. Aerating the system

After completion of the component tests, the display shows a prompt to remove the plug.



Remove the plug and press the "Execute" button. The system is aerated and the indicator will show the current pressure in the measurement chamber.



(pressure value is only an example)

9. Ending the self-test

After aeration of the system, the display changes automatically and indicates whether the unit passed the test. In the case of an error, the system stops the procedure and indicates an error code (see page 64), and the system is aerated.






The test can be repeated by pressing the "Start/Execute" button. The "Back" button leads directly to the Main Menu.






Perform self test

The self test determines if there are any leaks in the BRAND PLT unit. Two different tests are available. The Quick Test lasts only about 2 minutes, and should be performed on a daily basis. The Full Test lasts about 12 minutes. It includes the additional differentiated testing of the internal tanks, and should be performed on a weekly basis.

Full test

What to do?	How to do it	Keys to press	Display readout
1. Calling up the Main Menu	The 'Main Menu' is called up by pressing or turning the select knob.		
2. Setting up the self-test	Set to 'Self Test' by turning the select knob.		
3. Confirming the self-test	Press the 'Execute' button. The 'Self Test' submenu is called up. The last test to be performed appears enframed.		
4. Selecting the Full Test	Set to 'Selecting the Full Test' by turning the select knob.		
5. Confirming the Full Test	Press the 'Start/Execute' button. The 'Start Menu' submenu is called up.		
6. Seal the system	Press the provided plug firmly into the adapter (for the multi-channel adapter, seal all 4 channels).		

What to do?	How to do it	Keys to press	Display readout
<p>7. Starting the Full Test</p>	<p>Press the 'Start/Execute' button. The six components - sensor, pump, valve, filter, measurement chamber, and the tanks - are tested one after the next. The progress bar in the display indicates the test status. The LED indicator alternately blinks red and green during the test.</p>		 <p>(pressure value is only an example)</p>
<p>8. Aerating the system</p>	<p>After completion of the component tests, the display shows a prompt to remove the plug.</p>		
	<p>Remove the plug and press the 'Execute' button. The system is aerated and the indicator will show the current pressure.</p>		 <p>(pressure value is only an example)</p>
<p>9. Ending the self-test</p>	<p>After aeration of the system, the display changes automatically and indicates whether the unit passed the test. In the case of an error, the system stops the procedure and indicates an error code (see page 64), and the system is aerated.</p>		

Test preparation

Before the test, one must decide whether a dynamic or static test is to be performed.

Dynamic or static test?

With the **dynamic test**, in addition to the detection of static leaks, one can also determine whether a defective piston (contamination, scratches) has caused a leak.

During the measurement time, the pipette button must be pressed down completely and evenly 2-3 times. The associated piston movement allows errors to be recognized.

The button stroke is accomplished relatively slowly with 5 ml and 10 ml pipettes.

In the dynamic testing of pipettes, the piston should be in the upper position at the beginning, and should be brought back to that position at the end.

In the **static test**, by contrast, the pipette button is not pressed during the test procedure, i.e. the piston doesn't move. In power pipettes, one can position the piston in the uppermost location or at the 1st stop.

A leak that occurs during the piston stroke can thus not be determined at this point.

We recommend that the dynamic test be performed.

Test with or without pipette tip?

The test can take place with or without the pipette tip.

We recommend that the test be performed with a mounted, unused tip so that the entire pipette system is tested.

When a leak has been identified, the test can be repeated without a tip to determine whether the leak arises from the tip cone/tip coupling region. The adapter must be changed (see page 46) to be able to perform the test without the tip.

The pipette should be held firmly during the test.

Testing of single- or multi-channel pipettes

Changing the adapter

A total of 3 adapters are available

Adapter 1: already premounted, for testing single-channel pipettes up to 10 ml with a mounted pipette tip

Adapter 2: included in the delivery, for testing single-channel pipettes without a pipette tip

Adapter 3: available separately (optional), for testing of multi-channel pipettes with or without a pipette tip, wherein all 4 channels are tested simultaneously. The 4-channel adapter can be used for 4-, 8- and 12-channel pipettes.


When a leak is detected, three channels can be closed with plugs. The leaking channel can be determined through subsequent single-channel testing.

During multi-channel pipettes testing, the tips can be allowed to remain in the adapter when the instrument is disconnected. These can be removed simply by hand, or left inserted for the next test.

Important!

To protect the instrument from moisture and contamination, a readily changeable PE filter is included in the bottom the silicone adapter.



Test preparation

Changing the adapter and filter

1. Unscrewing the adapter lock nut.
2. Remove the pipette adapter.
3. If the PE filter is clogged with moisture or contaminants, it must be changed. In this case, carefully withdraw the PE filter from the bottom end of the silicone adapter.
4. Insert a new filter.
5. Insert the selected pipette adapter into the adapter housing (in the figure, Adapter 2 for pipettes without tips).
6. Screw the lock nut down firmly.
7. After changing the adapter, perform the Quick Test.

Note:

The 4-channel adapter consists of an adapter with an integrated lock nut. The entire ensemble is exchanged.

Testing single-channel pipettes









Example: 100 µl Air-displacement pipette, dynamic test

Before testing, mount a new, unused, disposable tip, and insert the pipette with the tip in the previously selected and attached Adapter 1 (see page 46).

What to do?	How to do it	Keys to press	Display readout												
1. Calling up the Main Menu	The 'Main Menu' is called up by pressing or turning the select knob.		<table border="1"> <tr><th colspan="2">Main Menu</th></tr> <tr><td>Pipette Type</td><td></td></tr> <tr><td>Self Test</td><td></td></tr> <tr><td>Switch OFF</td><td></td></tr> <tr><td>Settings...</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
2. Selecting the pipette type	Set to 'Pipette Type' by turning the select knob.		<table border="1"> <tr><th colspan="2">Main Menu</th></tr> <tr><td>Pipette Type</td><td></td></tr> <tr><td>Self Test</td><td></td></tr> <tr><td>Switch OFF</td><td></td></tr> <tr><td>Settings...</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
3. Confirming the pipette type	Press the 'Execute' button. The 'Channels' submenu is called up.		<table border="1"> <tr><th colspan="2">Channels</th></tr> <tr><td>Single Channel</td><td></td></tr> <tr><td>Multichannel</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Channels		Single Channel		Multichannel		select	back				
Channels															
Single Channel															
Multichannel															
select	back														
4. Selecting the single-channel pipette	Set to 'Single-Channel' by turning the select knob.		<table border="1"> <tr><th colspan="2">Channels</th></tr> <tr><td>Single Channel</td><td></td></tr> <tr><td>Multichannel</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Channels		Single Channel		Multichannel		select	back				
Channels															
Single Channel															
Multichannel															
select	back														
5. Confirming the single-channel pipette	Press the 'Execute' button. The 'Nominal Volume' submenu is called up.		<table border="1"> <tr><th colspan="2">Nominal Volume</th></tr> <tr><td>25 µl</td><td></td></tr> <tr><td>50 µl</td><td></td></tr> <tr><td>100 µl</td><td></td></tr> <tr><td>200 µl</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Nominal Volume		25 µl		50 µl		100 µl		200 µl		select	back
Nominal Volume															
25 µl															
50 µl															
100 µl															
200 µl															
select	back														
6. Selecting the pipette volume	'100 µl' is selected by turning the select knob.		<table border="1"> <tr><th colspan="2">Nominal Volume</th></tr> <tr><td>25 µl</td><td></td></tr> <tr><td>50 µl</td><td></td></tr> <tr><td>100 µl</td><td></td></tr> <tr><td>200 µl</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Nominal Volume		25 µl		50 µl		100 µl		200 µl		select	back
Nominal Volume															
25 µl															
50 µl															
100 µl															
200 µl															
select	back														

Continued on next page


Testing single-channel pipettes (cont.)

What to do?	How to do it	Keys to press	Display readout
7. Confirming the pipette volume	Press the 'Start/Execute' button. The 'Start' submenu is called up.		
8. Starting the test	Press the 'Start/Execute' button. The test is begun, and the test pressure is set (in progress...)		
Press the pipette button 2-3 times	Shortly before the end of the preparation period, the LED light shines 1 time red and 1 time green, then the progress bar appears and the system automatically starts the measurement. At this time, the pipette button is pressed down evenly 2-3 times. The piston must return the initial position before the measurement is finished! The progress bar graphic indicates the remaining time and the currently achieved vacuum. The LED indicator blinks orange.		 <p>(pressure value is only an example)</p>  <p>orange</p>
9. Ending the test	After the test is finished, the display indicates the leak rate that was determined and whether the system passed the test. The test can be repeated by pressing the "Start/Execute" button. Turning the select knob leads directly to the Main Menu.	 <p>green</p>  <p>red</p>	 <p>(leak rate is only an example)</p>  <p>(leak rate is only an example)</p>

Note:

The missing volume can be approximately determined using the correlation table on page 61.

What do the vertical progress bars in the display mean?

A mark is added to the progress bar. This indicates the threshold value. This corresponds to approx. 1/4 of the volume tolerance

- The progress bar is completely filled in to the bottom: The pipette is leak-tight.
- The progress bar is filled in from the top only down to the threshold value region: We recommend that a gravimetric test be performed.
- The progress bar is filled in to a point well above the threshold value: The pipette leaks. A gravimetric test should be performed, and if necessary the unit should be sent for repair.

What to do if there are leaks?

If a dynamic test indicates the presence of leaks, one can attempt to localize the source with further testing:

- A static test should be performed next.
The pipette button is not activated, and the piston is not moved. If the pipette is now leak-tight, this suggests there are scratches on the piston.

- The dynamic test should be performed without a tip.

Adapter 2 must be used in this case. If the pipette is now leak-proof, it is possible that an unsuitable tip was used.

What do the letters after some volumes mean?

The additional letters after some volumes during volume selection mean:

y = Yellow (instrument with yellow tips)

c = Clear tips

n = nano-cap tip

soft = During the testing of 5 ml and 10 ml pipettes, the piston might no longer move all the way back to top during the dynamic test, but might remain stuck in the instrument due to vacuum and friction. In this case, select the "5 ml soft" or "10 ml soft" test mode.

Note on duration of tests

In the first test after switching the instrument on, a relatively large volume must be pumped out, after which the instrument is given time to stabilize (relaxation). Regularly switching between large and small pipette types leads to longer test durations because the measurement pressure must be reset each time.



The measurement time for pipettes up to 1000 μl is generally 5 s, while it can be 20 s for 10 ml instruments.

Testing multichannel pipettes

Example: 200 µl 8-channel air-displacement pipette, dynamic test

Before testing, mount a new, unused, disposable tip, and insert four channels with tips in the previously selected and attached Adapter 3 (see page 46).

What to do?	How to do it	Keys to press	Display readout												
1. Calling up the Main Menu	The 'Main Menu' is called up by pressing or turning the select knob.		<table border="1"> <tr><th colspan="2">Main Menu</th></tr> <tr><td>Pipette Type</td><td></td></tr> <tr><td>Self Test</td><td></td></tr> <tr><td>Switch OFF</td><td></td></tr> <tr><td>Settings...</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
2. Selecting the pipette type	Set to 'Pipette Type' by turning the select knob.		<table border="1"> <tr><th colspan="2">Main Menu</th></tr> <tr><td>Pipette Type</td><td></td></tr> <tr><td>Self Test</td><td></td></tr> <tr><td>Switch OFF</td><td></td></tr> <tr><td>Settings...</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
3. Confirming the pipette type	Press the 'Execute' button. The 'Channels' submenu is called up.		<table border="1"> <tr><th colspan="2">Channels</th></tr> <tr><td>Single Channel</td><td></td></tr> <tr><td>Multichannel</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Channels		Single Channel		Multichannel		select	back				
Channels															
Single Channel															
Multichannel															
select	back														
4. Selecting the multi-channel pipette	Set to 'Multichannel' by turning the select knob.		<table border="1"> <tr><th colspan="2">Channels</th></tr> <tr><td>Single Channel</td><td></td></tr> <tr><td>Multichannel</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Channels		Single Channel		Multichannel		select	back				
Channels															
Single Channel															
Multichannel															
select	back														
5. Confirming the multichannel pipette	Press the 'Execute' button. The 'Nominal Volume' submenu is called up.		<table border="1"> <tr><th colspan="2">Nominal Volume</th></tr> <tr><td>50 µl</td><td></td></tr> <tr><td>100 µl</td><td></td></tr> <tr><td>200 µl</td><td></td></tr> <tr><td>250 µl</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Nominal Volume		50 µl		100 µl		200 µl		250 µl		select	back
Nominal Volume															
50 µl															
100 µl															
200 µl															
250 µl															
select	back														
6. Selecting the pipette volume	'200 µl' is selected by turning the select knob.		<table border="1"> <tr><th colspan="2">Nominal Volume</th></tr> <tr><td>50 µl</td><td></td></tr> <tr><td>100 µl</td><td></td></tr> <tr><td>200 µl</td><td></td></tr> <tr><td>250 µl</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Nominal Volume		50 µl		100 µl		200 µl		250 µl		select	back
Nominal Volume															
50 µl															
100 µl															
200 µl															
250 µl															
select	back														

What to do?	How to do it	Keys to press	Display readout
7. Confirming the pipette volume	Press the 'Start/Execute' button. The ' Start ' submenu is called up.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center; font-weight: bold;">ready</p> <p style="text-align: center;">MC – 200 µl</p> <hr/> <p style="text-align: left;">start</p> </div>
8. Starting the test	Press the 'Start/Execute' button. The test is begun, and the test pressure is set (in progress...)		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center; font-weight: bold;">Test</p> <p style="text-align: center;">preparing...</p> <hr/> <p style="text-align: left;">restart back</p> </div>
Press the pipette button 2-3 times	Shortly before the end of the preparation period, the LED light shines 1 time red and 1 time green, then the progress bar appears and the system automatically starts the measurement. At this time, the pipette button is pressed down evenly 2-3 times. The piston must return the initial position before the measurement is finished! The progress bar graphic indicates the remaining time and the currently achieved vacuum. The LED indicator blinks orange.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center; font-weight: bold;">Test</p> <div style="text-align: right;">  4s </div> <p style="text-align: center;">-121.21 hPa</p> <hr/> <p style="text-align: left;">restart back</p> <p style="font-size: small;">(pressure value is only an example)</p> <div style="text-align: center;">  <p style="font-weight: bold; color: orange;">orange</p> </div> </div>
9. Ending the test	After the test is finished, the display indicates the leak rate that was determined and whether the system passed the test. The 'Start/Execute' button can be used to repeat the test, or the test can be continued with the next 4 channels. Turning the select knob leads directly to the Main Menu.	 <p style="text-align: center; font-weight: bold; color: green;">green</p>  <p style="text-align: center; font-weight: bold; color: red;">red</p>	<div style="border: 1px solid black; padding: 5px; margin-bottom: 10px;"> <p style="text-align: center; font-weight: bold; font-size: 1.2em;">PASS</p> <p style="text-align: center;">MC – 200 µl</p> <p style="text-align: center;">Q_L: 0.04 hPa*ml/s</p> <hr/> <p style="text-align: left;">start</p> <p style="font-size: small;">(leak rate is only an example)</p> </div> <div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center; font-weight: bold; font-size: 1.2em;">FAIL</p> <p style="text-align: center;">MC – 200 µl</p> <p style="text-align: center;">Q_L: 4.6 hPa*ml/s</p> <hr/> <p style="text-align: left;">start</p> <p style="font-size: small;">(leak rate is only an example)</p> </div>

Note:


The missing volume can be approximately determined using the correlation table on page 61.

Testing multichannel pipettes

What do the vertical progress bars in the display mean?

A mark is added to the progress bar. This indicates the threshold value. This corresponds to approx. $\frac{1}{4}$ of the volume tolerance.

- The progress bar is completely filled in to the bottom: The instrument is leak-tight.
- The progress bar is filled in from the top only down to the threshold value region: We recommend that a gravimetric test be performed.
- The progress bar is filled in to a point well above the threshold value: The instrument leaks. A gravimetric test should be performed, and if necessary the unit should be sent for repair.

What to do if there are leaks?

If a dynamic test indicates the presence of leaks, one can attempt to localize the source with further testing:

- A static test should be performed next.
The pipette button is not activated, and the piston is not moved.
If the pipette is now leak-tight, this suggests there are scratches on the piston.
- The dynamic test should be performed without a tip.
Adapter 3 can be reused for this.
If the pipette is now leak-proof, it is possible that an unsuitable tip was used.
- Close 3 channels with the caps provided, and test the individual channels to localize the leak. Select the "single-channel" pipette type for this.





Note on duration of tests

In the first test after switching the instrument on, a relatively large volume must be pumped out, after which the instrument is given time to stabilize (relaxation). Regularly switching between large and small pipette types leads to longer test durations because the measurement pressure must be reset each time.

The measurement time for pipettes up to 500 μl is generally 5 s, while it can be 10 s for a volume of 1250 μl .

Switch the instrument off

The active PLT unit can be switched off at any time, independent of the shut-off time set (see page 57).


What to do?	How to do it	Keys to press	Display readout
1. Calling up the Main Menu	The ' Main Menu ' is called up by pressing the select knob.		
2. Selecting instrument OFF	Set to ' Instrument OFF ' by turning the select knob.		
3. Switch the instrument off	Press the left ' Execute ' button to switch the instrument off.		
Turning the instrument on again	Pressing any button turns the instrument on again. At first, the display shows the scheduled maintenance date for a few seconds Then it switches automatically to the 'Ready' indicator with the previously tested pipette.		 

Selecting the language


In the 'Settings' submenu, the following subitems can be selected: Language, contrast, maintenance date, shut-off time, pressure units, and factory defaults.

What to do?	How to do it	Keys to press	Display readout												
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing the select knob, and turn it to select 'Settings' .		<table border="1"><tr><td colspan="2">Main Menu</td></tr><tr><td colspan="2">Pipette Type</td></tr><tr><td colspan="2">Self Test</td></tr><tr><td colspan="2">Switch OFF</td></tr><tr><td colspan="2">Settings...</td></tr><tr><td>select</td><td>back</td></tr></table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
2. Confirming the settings	Press the 'Execute' button. The 'Settings' submenu is called up.		<table border="1"><tr><td colspan="2">Settings</td></tr><tr><td colspan="2">Language</td></tr><tr><td colspan="2">Contrast</td></tr><tr><td colspan="2">Inspection date</td></tr><tr><td colspan="2">Auto-Power-Off Time</td></tr><tr><td>select</td><td>back</td></tr></table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														
3. Selecting the language	Set to 'Language' by turning the select knob.		<table border="1"><tr><td colspan="2">Settings</td></tr><tr><td colspan="2">Language</td></tr><tr><td colspan="2">Contrast</td></tr><tr><td colspan="2">Inspection date</td></tr><tr><td colspan="2">Auto-Power-Off Time</td></tr><tr><td>Auswahl</td><td>zurück</td></tr></table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		Auswahl	zurück
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
Auswahl	zurück														
4. Calling up the Language Menu	Press the 'Execute' button. The 'Language' menu is called up. The language is selected by turning the select knob.		<table border="1"><tr><td colspan="2">Language</td></tr><tr><td colspan="2">Deutsch</td></tr><tr><td colspan="2">English</td></tr><tr><td colspan="2">Español</td></tr><tr><td>o.k.</td><td>back</td></tr></table>	Language		Deutsch		English		Español		o.k.	back		
Language															
Deutsch															
English															
Español															
o.k.	back														
5. Confirming the language	Press the 'Execute' button. The 'Settings' menu appears in the language selected. The 'Back' button leads to the Main Menu.		<table border="1"><tr><td colspan="2">Settings</td></tr><tr><td colspan="2">Language</td></tr><tr><td colspan="2">Contrast</td></tr><tr><td colspan="2">Inspection date</td></tr><tr><td colspan="2">Auto-Power-Off Time</td></tr><tr><td>select</td><td>back</td></tr></table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														


The color intensity of the indicator can be changed interactively.

What to do?	How to do it	Keys to press	Display readout												
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing the select knob, and turn it to select 'Settings'.		<table border="1"> <tr><th colspan="2">Main Menu</th></tr> <tr><td>Pipette Type</td><td></td></tr> <tr><td>Self Test</td><td></td></tr> <tr><td>Switch OFF</td><td></td></tr> <tr><td>Settings...</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		select	back
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
select	back														
2. Confirming the settings	Press the 'Execute' button. The 'Settings' submenu is called up.		<table border="1"> <tr><th colspan="2">Settings</th></tr> <tr><td>Language</td><td></td></tr> <tr><td>Contrast</td><td></td></tr> <tr><td>Inspection date</td><td></td></tr> <tr><td>Auto-Power-Off Time</td><td></td></tr> <tr><td>Auswahl</td><td>zurück</td></tr> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		Auswahl	zurück
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
Auswahl	zurück														
3. Selecting the contrast	Set to 'Contrast' by turning the select knob.		<table border="1"> <tr><th colspan="2">Settings</th></tr> <tr><td>Language</td><td></td></tr> <tr><td>Contrast</td><td></td></tr> <tr><td>Inspection date</td><td></td></tr> <tr><td>Auto-Power-Off Time</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														
4. Calling up the contrast	Press the 'Execute' button. The 'Contrast' menu is called up.		<table border="1"> <tr><th colspan="2">Contrast</th></tr> <tr><td></td><td>50</td></tr> <tr><td>o.k.</td><td>back</td></tr> </table>	Contrast			50	o.k.	back						
Contrast															
	50														
o.k.	back														
5. Setting the contrast	Set this from '10' to '100' by turning the select knob.		<table border="1"> <tr><th colspan="2">Contrast</th></tr> <tr><td></td><td>76</td></tr> <tr><td>o.k.</td><td>back</td></tr> </table>	Contrast			76	o.k.	back						
Contrast															
	76														
o.k.	back														
6. Confirming the contrast	Press the 'Execute' button. The 'Settings' menu reappears. The 'Back' button leads to the Main Menu.		<table border="1"> <tr><th colspan="2">Settings</th></tr> <tr><td>Language</td><td></td></tr> <tr><td>Contrast</td><td></td></tr> <tr><td>Inspection date</td><td></td></tr> <tr><td>Auto-Power-Off Time</td><td></td></tr> <tr><td>select</td><td>back</td></tr> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		select	back
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
select	back														

The interval for when the PLT unit should be sent to BRAND or an authorized repair shop for testing can be set individually. This interval is displayed each time the unit is switched on.

What to do?	How to do it	Keys to press	Display readout
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing the select knob, and turn it to select 'Settings' .		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings...</p> <p style="text-align: right;">select back</p> </div>
2. Confirming the settings	Press the 'Execute' button. The 'Settings' submenu is called up.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> <p style="text-align: right;">select back</p> </div>
3. Selecting the inspection date	Set to 'Inspection date' by turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Language Contrast Inspection date Auto-Power-Off Time</p> <p style="text-align: right;">select back</p> </div>
4. Calling up the inspection date	Press the 'Execute' button. The previously set maintenance date is displayed.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Inspection date</p> <p style="text-align: center;">1 - Apr - 2010</p> <p style="text-align: right;">o.k. back</p> </div>
5. Setting the inspection date	Turning the select knob changes the selection field in the series Year, Month, and Day. (Pressing today's date leads back to 'Settings'). Turning the select knob changes which position is enframed.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Inspection date</p> <p style="text-align: center;">1 - Apr - 2011</p> <p style="text-align: right;">o.k. select</p> </div>
6. Confirming the inspection date	Press the 'Execute' button. The indicator returns automatically to the 'Settings' menu. The 'Back' button leads to the Main Menu.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Pipette Type Self Test Switch OFF Settings...</p> <p style="text-align: right;">select back</p> </div>

The auto-power-off time of the PLT unit can be set individually by minutes from 'OFF' to 30 min.

What to do?	How to do it	Keys to press	Display readout												
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing or turning the select knob, and turn it to select 'Settings'.		<table border="1"> <tr><td colspan="2">Main Menu</td></tr> <tr><td colspan="2">Pipette Type</td></tr> <tr><td colspan="2">Self Test</td></tr> <tr><td colspan="2">Switch OFF</td></tr> <tr><td colspan="2">Settings...</td></tr> <tr><td>Auswahl</td><td>zurück</td></tr> </table>	Main Menu		Pipette Type		Self Test		Switch OFF		Settings...		Auswahl	zurück
Main Menu															
Pipette Type															
Self Test															
Switch OFF															
Settings...															
Auswahl	zurück														
2. Confirming the settings	The 'Settings' submenu is called up by pressing the 'Execute' button.		<table border="1"> <tr><td colspan="2">Settings</td></tr> <tr><td colspan="2">Language</td></tr> <tr><td colspan="2">Contrast</td></tr> <tr><td colspan="2">Inspection date</td></tr> <tr><td colspan="2">Auto-Power-Off Time</td></tr> <tr><td>Auswahl</td><td>zurück</td></tr> </table>	Settings		Language		Contrast		Inspection date		Auto-Power-Off Time		Auswahl	zurück
Settings															
Language															
Contrast															
Inspection date															
Auto-Power-Off Time															
Auswahl	zurück														
3. Selecting the auto-power-off time	Set to 'Auto-Power-Off time' by turning the select knob.		<table border="1"> <tr><td colspan="2">Settings</td></tr> <tr><td colspan="2">Contrast</td></tr> <tr><td colspan="2">Inspection date</td></tr> <tr><td colspan="2">Auto-Power-Off Time</td></tr> <tr><td colspan="2">Pressure Unit</td></tr> <tr><td>select</td><td>back</td></tr> </table>	Settings		Contrast		Inspection date		Auto-Power-Off Time		Pressure Unit		select	back
Settings															
Contrast															
Inspection date															
Auto-Power-Off Time															
Pressure Unit															
select	back														
4. Calling up the auto-power-off time	Press the 'Execute' button. The previously set auto-power-off time is displayed.		<table border="1"> <tr><td colspan="2">Auto-Power-Off Time</td></tr> <tr><td colspan="2">3 min</td></tr> <tr><td>o.k.</td><td>back</td></tr> </table>	Auto-Power-Off Time		3 min		o.k.	back						
Auto-Power-Off Time															
3 min															
o.k.	back														
5. Setting the auto-power-off time	The auto-power-off time can be set from 'OFF' to 30 min.		<table border="1"> <tr><td colspan="2">Auto-Power-Off Time</td></tr> <tr><td colspan="2">15 min</td></tr> <tr><td>o.k.</td><td>back</td></tr> </table>	Auto-Power-Off Time		15 min		o.k.	back						
Auto-Power-Off Time															
15 min															
o.k.	back														
6. Confirming the auto-power-off time	Press the 'Execute' button. The indicator returns automatically to the 'Settings' menu. The 'Back' button leads to the Main Menu.		<table border="1"> <tr><td colspan="2">Settings</td></tr> <tr><td colspan="2">Contrast</td></tr> <tr><td colspan="2">Inspection date</td></tr> <tr><td colspan="2">Auto-Power-Off Time</td></tr> <tr><td colspan="2">Pressure Unit</td></tr> <tr><td>select</td><td>back</td></tr> </table>	Settings		Contrast		Inspection date		Auto-Power-Off Time		Pressure Unit		select	back
Settings															
Contrast															
Inspection date															
Auto-Power-Off Time															
Pressure Unit															
select	back														

What to do?	How to do it	Keys to press	Display readout
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing or turning the select knob, and turn it to select 'Settings'.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings...</p> <hr/> <p>select back</p> </div>
2. Confirming the settings	Press the 'Execute' button. The 'Settings' submenu is called up.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Contrast Inspection date Auto-Power-Off Time Pressure Unit</p> <hr/> <p>select back</p> </div>
3. Selecting the pressure unit	Set to 'Pressure unit' by turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure Unit Factory Settings</p> <hr/> <p>select back</p> </div>
4. Calling up the pressure units	Press the left 'Execute' button. The currently set pressure units appear.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Pressure Unit</p> <p style="text-align: center;">hPa</p> <hr/> <p>o.k. back</p> </div>
5. Setting the pressure units	The pressure units can be set by turning the select knob to hPa, mbar, or Torr.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Pressure Unite</p> <p style="text-align: center;">mbar</p> <hr/> <p>o.k. back</p> </div>
6. Confirming the pressure units	Press the 'Execute' button. The indicator returns automatically to the 'Settings' menu. The 'Back' button leads to the Main Menu.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure Unit Factory Settings</p> <hr/> <p>select back</p> </div>

What to do?	How to do it	Keys to press	Display readout
1. Selecting 'Settings' in the Main Menu	Call up the 'Main Menu' by pressing the select knob, and turn it to select 'Settings'.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Main Menu</p> <p>Pipette Type Self Test Switch OFF Settings...</p> <hr/> <p>select back</p> </div>
2. Confirming the settings	Press the 'Execute' button. The 'Settings' submenu is called up.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure Unit Factory Settings</p> <hr/> <p>select back</p> </div>
3. Selecting the factory settings	Set to 'Factory Settings' by turning the select knob.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure Unit Factory Settings</p> <hr/> <p>select back</p> </div>
4. Calling up the factory settings	Press the 'Execute' button. The confirmation prompt appears.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Factory Settings</p> <p style="text-align: center;">please confirm!</p> <hr/> <p>o.k. back</p> </div>
5. Confirming the factory settings	Press the 'Execute' button. The indicator changes to the language selection of 'English'.		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Language</p> <p>Deutsch English Espanol</p> <hr/> <p>o.k. </p> </div>
	After selecting and confirming the language, the indicator changes back automatically to the 'Settings' menu. The 'Back' button leads to the Main Menu		<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p>Inspection date Auto-Power-Off Time Pressure Unit Factory Settings</p> <hr/> <p>select back</p> </div>

Firmware update

Firmware updates, e.g. when an additional pipette type (new volume) is to be included, are made available by BRAND on CD. Registered users are informed automatically. Unregistered users will receive this update upon request.

Installation

1. Place the CD-ROM in the drive (Drive "D" is used in this example).
2. Click on 'Start' on the desktop, and enter 'D:\setup.exe' in the 'Execute' command field to install the firmware update.
3. Follow the instructions on the screen.

Relationship of QL to volume loss (approx.)

With a leak of a given magnitude, the volume loss of the pipette during pipetting depends on various factors, e.g., the tip opening.

(Tolerance for the measured leak rate: $\pm 20\%$)

Pipette tested		Leak rate Q_L [hPa*ml/s]	Pipette volume loss (approx.) [μ l]
Channels	Nominal volume		
1-channel	1 μ l	1	0,02
1-channel	2 μ l	2	0,035
1-channel	2,5 μ l	2	0,04
1-channel	5 μ l y	5	0,05
1-channel	5 μ l	2	0,05
1-channel	10 μ l y	6	0,1
1-channel	10 μ l c/n	2	0,1
1-channel	20 μ l y	7	0,16
1-channel	20 μ l c	3	0,16
1-channel	25 μ l	7	0,2
1-channel	50 μ l	8	0,4
1-channel	100 μ l	5	0,6
1-channel	200 μ l	6	1,2
1-channel	250 μ l	6	1,5
1-channel	300 μ l	6	1,8
1-channel	500 μ l	8	3
1-channel	1000 μ l	10	6
1-channel	1250 μ l	15	7,5
1-channel	1500 μ l	32	9
1-channel	2000 μ l	59	15
1-channel	2500 μ l	39	15
1-channel	5 ml	31	30
1-channel	5 ml soft	15	30
1-channel	10 ml	38	60
1-channel	10 ml soft	19	60
4-channel	1 μ l	1	0,02
4-channel	2 μ l	2	0,035
4-channel	2,5 μ l	2	0,04
4-channel	5 μ l	2	0,05
4-channel	10 μ l	2	0,1
4-channel	20 μ l	3	0,16
4-channel	25 μ l	7	0,2
4-channel	50 μ l	8	0,4
4-channel	100 μ l	5	0,6
4-channel	200 μ l	6	1,2
4-channel	250 μ l	6	1,5
4-channel	300 μ l	6	1,8
4-channel	500 μ l	11	3
4-channel	1000 μ l	9	6
4-channel	1250 μ l	12	7,5

PLT unit (Pipet Leak Testing Unit)

Including one adapter each for testing air-displacement pipettes with tip (mounted) and without tip, 2 blind plugs, 3 replacement PE filters for the pipette adapters, a universal AC adapter, and a quality certificate.

Cat. No. 7039 70

Verp.-Einh. ???

Pipette adapter for testing with tips

for single-channel pipettes, including 1 blind plug

Cat. No. 7039 75

Pipette adapter for testing without tips

for single-channel pipettes, including 1 blind plug

Cat. No. 7039 76

4-Channel pipette adapter

for testing multi-channel pipettes with and without tips, including 4 blind plugs

Cat. No. 7039 77

Filters for the pipette adapter

PE, package of 10

Cat. No. 7039 78


Universal AC adapter

Input: AC 100 V - 240 V, 50/60 Hz
Output: DC 6.5 V, 800 mA

Cat. No. 7039 79


Errors during pipette tests

Problem	Error display	Possible cause	Corrective action
Test interruption		No pipette in adapter, or not inserted firmly enough	Insert pipette firmly enough
		Tested pipette is completely leaky	Pipette is defective, and cannot be tested with this testing instrument. Send the pipette out for repair.
		Tip does not seal onto shaft	Insert another tip Check shaft for damage
		Wrong pipette type adjusted	Select correct pipette type in the menu.
		Pipette adapter not leak-tight	Check adapter for lint or abrasion, and check filter for moisture. Clean adapter and exchange filter. Perform self-test.
		Pipette adapter is loose	Tighten the lock nut firmly. Perform self-test
		Moisture in the measurement room, adapter filter is clogged.	Exchange adapter filter, perform self-test. (Quick Test)
Measure leak rate implausible	Moisture or leaks in the PLT unit	Perform self-test (Full Test).	
PLT doesn't react at all, or behaves unusually.	Software error, electrical malfunction	Disconnect all cables from the instrument. Reconnect the power cable after approx. 5 seconds.	

Note: If the error cannot be resolved by the measures described, send the instrument out for repair.

Troubleshooting

Error during the self-test

Error display

Error number	Cause	Corrective action
1000 1010	Pressure sensor is defective	Send the instrument out for repair.
1020	Pump running continuously	Send the instrument out for repair.
1040 1050	Pumping down takes too long	Check whether the caps are firmly attached, and that the adapter is firmly screwed in. If they already are, then the pump is defective. Send the instrument out for repair.
1030 1060 1070 1080 1090	Defective valve	Send the instrument out for repair.
1020	Aeration through the aeration channel takes too long	Send the instrument out for repair.
1150 1160	Instrument leaks	Check the caps, adapter seating, and adapter filter, and replace if necessary. Repeat the self-test. If the error is not resolved, send the instrument out for repair.
1170	Aeration through the adapter takes too long	Remove the cap if this hasn't been done. If the adapter filter is contaminated, replace the filter.
1200	Electronic defect	Send the instrument out for repair.

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!

Return for Repair

- a) Clean and decontaminate the instrument carefully.
- b) Complete the 'Declaration on Absence of Health Hazards' (ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de).
- c) Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used.

The return transport of the instrument is at risk and cost of the sender.


Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

Disposal

The adjoining symbol means that storage batteries and electronic devices must be disposed of separately from household trash (mixed municipal waste) at the end of their service life.

- According to the Directive 2002/96/EC of the European Parliament and of the Council on Waste Electrical and Electronic Equipment (WEEE) of 27 January 2003, electronic equipment requires disposal according to the relevant national disposal regulations.

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