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| | Page |
|--------------------------------|------|
| Safety Instructions | 31 |
| Limitations of use | 32 |
| Operating and Control Elements | 33 |
| Pipetting | 34 |
| Checking the Volume | 37 |
| Accuracy Table | 39 |
| Adjustment – Easy Calibration | 40 |
| Autoclaving | 41 |
| UV sterilization | 41 |
| Servicing and Cleaning | 42 |
| Disassembling the Manifold | 44 |
| Ordering Information | 48 |
| Accessories | 49 |
| Spare Parts | 50 |
| Troubleshooting | 51 |
| Repairs | 52 |
| Contact addresses | 53 |
| Calibration Service | 54 |
| Warranty Information | 54 |
| Disposal | 54 |

Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

- 1.** Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
- 2.** Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
- 3.** Observe all specifications provided by reagent manufacturers.
- 4.** Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 32). If in doubt, contact the manufacturer or supplier.
- 5.** Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
- 6.** Avoid touching the tip orifices when working with hazardous samples.
- 7.** Never use force on the instrument!
- 8.** Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
- 9.** Before use check the instrument for visible damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 50), and contact the manufacturer if needed.

Limitations of Use

Purpose

The pipette is an air-displacement system for pipetting aqueous solutions with medium density and viscosity.

Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between +15 °C to +40 °C (59 °F to 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

Operating Limitations


Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than ± 1 °C/ ± 1.8 °F.

Operating Exclusions

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene and FKM
- for liquids attacking polycarbonate (viewing window)
- for liquids attacking polyphenyl sulfide (50 μ l, 100 μ l, 200 μ l and 300 μ l models)
- for liquids of a very high vapor pressure

Note:

Optimum analysis results can only be obtained with quality tips. We recommend pipette tips from BRAND. For further information, refer to the accuracy table on page 39.

1. Fitting the tips

Use the correct tips according to the volume range or the color code. Ensure that the tips are securely seated.

Pipette tips are disposables items!

2. Volume setting

- Push the volume-change protection upward to disengage (UNLOCK).
- Select the desired volume by rotating the volume-setting wheel. Avoid twisting and abrupt rotating motions during this adjustment.
- Push the volume-change protection down to re-engage (LOCK). Note: The volume-change protection tightens but does not lock volume-setting wheel.

3. Align the manifold

The manifold can turn freely in both directions.

4. Aspirate sample

- Press pipetting button to the first stop.
- Hold the pipette vertically and immerse the tips into the liquid.

| Volume range | Immersion depth in mm | Waiting time in s |
|------------------------------|-----------------------|-------------------|
| 0.1 μ l - 1 μ l | 1 - 2 | 1 |
| > 1 μ l - 100 μ l | 2 - 3 | 1 |
| > 100 μ l - 1000 μ l | 2 - 4 | 1 |
| > 1000 μ l | 3 - 6 | 3 |

- c) Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tips immersed for a few seconds.


5. Discharge sample

- a) Place the pipette tips against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.
- b) Press the pipetting button slowly to the first stop and hold it down. For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- c) The blow-out stroke empties the tips completely: Press the pipetting button down to the second stop.
- d) While doing this, wipe the pipette tips against the wall of the container.
- e) Remove the pipette tips from the container wall and let the pipetting button slide back.

6. Ejecting the tips

Hold the manifold over a suitable disposal container and press the tip ejection key to the stop.

Note:

ISO 8655 prescribes rinsing the pipette tips once with the sample liquid prior to the actual pipetting process.

Important!

Don't lay the instrument horizontal when the tips are filled. Liquid may enter and contaminate the instrument. The instrument should be stored without tips, placed upright in the supplied shelf/rack mount or bench top rack.

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 34 for procedure).

2. Condition the pipette

Condition the pipette before testing by using pipette tips to aspirate and discharge the test liquid (distilled H₂O) five times.

3. Carry out the test

Note:

According to DIN EN ISO 8655-2, it is recommended to replace the tips after each measurement. This rule may be waived in accordance with the DAkkS guideline DKD-R8-1.

- a) Aspirate liquid and pipette it into the weighing vessel.

Note:

Each individual channel must be tested separately.

- b) Weigh the pipetted quantity with an analytical balance. (Please follow the operating manual instructions from the balance manufacturer.)
- c) Calculate the volume, taking the temperature of the test liquid into account.
- d) 3-10 pipettings and weighings per channel in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

Checking the Volume

Calculation (for nominal volume)

x_i = Weighing results

n = Number of weighings

Z = Correction factor

(e.g., 1.0029 $\mu\text{l}/\text{mg}$
at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

Accuracy*

$$\text{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

V_0 = Nominal volume

Coefficient of Variation*

$$\text{CV\%} = \frac{100 s}{\bar{V}}$$

Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and variation coefficient (CV%):
A% and CV% are calculated according to the formulas for statistical control.

Note:

Testing instructions (SOPs) and a demo version of the EASYCAL™ 4.0 calibration software are available for download at www.brand.de.

Transferpette® S-8/-12

| Volume range μl | Volume step μl | A* ≤ ± % | CV* ≤ % | Increment μl | Rec. type of tips, μl |
|--------------------|-------------------|-------------|------------|-----------------|--------------------------|
| 0.5 - 10 | 10 | 1.6 | 1.0 | 0.01 | 0.5 - 20 |
| | 5 | 2 | 2 | | |
| | 1 | 8 | 6 | | |
| 5 - 50 | 50 | 0.8 | 0.4 | 0.1 | 2 - 200 |
| | 25 | 1.4 | 0.8 | | |
| | 5 | 6 | 3 | | |
| 10 - 100 | 100 | 0.8 | 0.3 | 0.2 | 2 - 200 |
| | 50 | 1.4 | 0.6 | | |
| | 10 | 4 | 2 | | |
| 20 - 200 | 200 | 0.8 | 0.3 | 0.2 | 2 - 200 |
| | 100 | 1.4 | 0.6 | | |
| | 20 | 4 | 1.5 | | |
| 30 - 300 | 300 | 0.6 | 0.3 | 0.5 | 5 - 300 |
| | 150 | 1.2 | 0.6 | | |
| | 30 | 3 | 1.5 | | |

* A = Accuracy, CV = Coefficient of Variation

DE-M  20 °C
Ex

Final test values related to the nominal capacity (maximum volume) or the indicated volume steps indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.


Adjustment – Easy Calibration

The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made using the Easy Calibration Technique.

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2. Remove the label window and the label. Push the hook forward, raise it slightly and then pull it back.


3. Using a paperclip or an unused pipette tip, remove the protective film (this protective film can be discarded).

4. Push the red adjustment slider completely back, raise the volume-setting wheel (decoupling) and release the adjustment slider.

5. Set the previously determined actual value (adjustment value) with the volume-setting wheel in the UNLOCK position.

A volume check is recommended after every adjustment.

6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Re-insert the label and the label window.

Note:

The change to the factory settings is indicated by the red adjustment slider now visible in the label window.

The Transferpette® S-8/-12 is completely autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi), with a holding time of at least 15 minutes according to DIN EN 285.

1. Eject the pipette tips.
2. Autoclave the complete pipette without any further disassembling.
3. Allow the Transferpette® S-8/-12 to completely cool and dry.

Note:

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

Attention:

Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between), with the volume-change protection set fully unlocked (UNLOCK).

If the pipette is autoclaved frequently, the pistons should be oiled with the supplied silicone oil in order to preserve smooth movement.

If necessary after autoclaving, tighten the connection between the hand grip and the manifold.

**UV sterilization**

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.


Servicing and Cleaning

The three main components of the manifold can be easily separated and disassembled for servicing, cleaning or replacing parts. The procedures are illustrated on the following pages.

Note: The changing of the V-rings/O-rings on the individual nose cones is described in detail in the instructions enclosed with the spare part.

The main components of the manifold

- A Piston unit** with piston support bar [I] and pistons inserted in this unit, which can be unscrewed individually for cleaning or replacement.
- B Nose cone assembly** with nose cone support bar [II] and central guide rod (Z), which is attached to this, and the nose cones and seals, which can be unscrewed individually for cleaning or replacement.
- C Manifold housing**, which is connected to the manifold housing cover [III] of the piston unit with two turn-lock fasteners.

Servicing

In order to assure proper functioning, the Transferpette® S-8/-12 should be serviced and cleaned at regular intervals.

What is to be inspected?

1. Inspect nose cones, pistons and seals for damage and contamination.
2. Test the sealing of the instrument. We recommend using the BRAND leak testing instrument PLT unit.

Alternatively: to do this aspirate a sample, and then hold the instrument in a vertical position for about 10 sec. If a drop forms at the tip orifices, see the troubleshooting guide, page 51.

Cleaning instructions

1. Clean single nose cones, pistons and nose cone support bar/piston support bar (**these components only**) with soap solution or isopropyl alcohol. Afterwards, rinse with distilled water.
2. Let these parts dry and cool down completely. Residual moisture in the nose cones may result in a loss of accuracy.
3. Lubricate the piston with a very thin coating of the silicone oil supplied. For the central guide rod (Z) only use the recommended fluorstatic grease!

Disconnecting handle from manifold

1. Eject the pipette tips.
2. To disconnect the manifold, pull it downward as far as possible, and **only then** turn it clockwise. After one rotation, it should no longer be pulled downward while it is being turned.

Attention:


Improper handling can damage the unit!


Disassembling the Manifold

Removing of nose cones and seals

for cleaning or replacing


1. Disconnect the manifold from the handle.
2. Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.

3. Push the mounting tool (M) (supplied with the pipette) on the nose cone and unscrew it.


Disassembling the Manifold

4. Push the piston unit down to the bottom. The seal will either remain inside the nose cone or will stay on the piston after the nose cone is removed. Remove the seal, inspect it and clean or replace if necessary. If required, lubricate the piston lightly with the supplied silicone oil.


(In the Transferpette® S-8/-12 30-300 µl unit, remove the additional pressure ring from the piston for cleaning.)

5. Push the seal on the piston with its flat side facing upward.

(In the Transferpette® S-8/-12 30-300 µl unit, first replace the additional pressure ring!)





Use the mounting tool to mount and tighten the cleaned or new nose cone.

6. Reassemble the manifold. The manifold must be screwed into the handle counter-clockwise, until it audibly snaps into place. Check the instrument for tightness (see page 43, "Servicing").

► Continued on page 46

Removing of pistons

for cleaning or replacing


1. Disconnect the manifold from the handle.
2. Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.
3. Remove both outer Phillips screws on the nose cone assembly.

Attention:


Do not loosen the central guide rod (Z)!

4. Pull the piston and nose cone assembly apart and separate. Remove stroke springs.
5. Place the mounting tool on the piston nut and unscrew the piston nut.


Disassembling the Manifold

- 6.** Remove the piston nut and pull out the piston with piston spring.

- 7.** Insert the piston spring (a) and cleaned or new piston. Screw on the piston nut and tighten with the mounting tool. Replace the stroke springs (b).

- 8.** Reassemble the manifold. The manifold must be screwed counter-clockwise into the handle so that it audibly snaps into place. Check the instrument for tightness, mobility and accuracy (see page 43, "Servicing").

Ordering Information

Items supplied

Confirm that your package includes: Transferpette® S -8 or -12 pipette, 2 x TipBox, filled with compatible pipette tips, container of silicone oil, 1 shelf/rack mount, 1 reagent reservoir, 1 mounting tool, 1 set of sealing rings made of FKM and this operating manual.

Transferpette® S -8

| Capacity | Description | Cat. No. |
|-------------|-------------|----------|
| 0.5 - 10 µl | M8-10 | 7037 00 |
| 5 - 50 µl | M8-50 | 7037 06 |
| 10 - 100 µl | M8-100 | 7037 08 |
| 20 - 200 µl | M8-200 | 7037 10 |
| 30 - 300 µl | M8-300 | 7037 12 |

Transferpette® S -12

| Capacity | Description | Cat. No. |
|-------------|-------------|----------|
| 0.5 - 10 µl | M12-10 | 7037 20 |
| 5 - 50 µl | M12-50 | 7037 26 |
| 10 - 100 µl | M12-100 | 7037 28 |
| 20 - 200 µl | M12-200 | 7037 30 |
| 30 - 300 µl | M12-300 | 7037 32 |

**Bench-top rack for 6 Transferpette® S
or 3 Transferpette® S-8/-12 pipettes**

Cat. No. 7048 05

**Bench-top rack for 1 Transferpette® S
-8/-12 pipette**

Cat. No. 7034 40

**Shelf/rack mount for 1 Transferpette® S
or 1 Transferpette® S-8/-12 pipette**

Cat. No. 7048 10

Additional accessories for Transferpette® S-8/-12

| Description | Cat. No. |
|---|----------|
| Label window , pack of 1 | 7046 50 |
| Blank labels , pack of 5 | 7046 51 |
| Silicone oil | 7036 77 |
| Fluorstatic grease | 7036 78 |
| Reagent reservoir , PP. Capacity 60 ml. Autoclavable (121 °C) | |
| non-sterile, with lid. Pack of 10. | 7034 59 |
| sterile, without lid. Packed individually. Pack of 100. | 7034 11 |
| sterile, without lid. 5 per bag. Pack of 200. | 7034 09 |
| PLT unit Pipette leak testing unit | 7039 70 |

Spare Parts

Parts will differ slightly depending on nominal volume of instrument.
 (Fig. shows spare parts for Transferpette® S-8/-12, 10-100 µl).

| Capacity | A | B* | C | D |
|-------------|---------|---------|---------|---------|
| 0.5 - 10 µl | 7056 59 | 7056 77 | 7033 80 | 7033 40 |
| 5 - 50 µl | 7056 66 | 7056 81 | 7034 91 | 7033 43 |
| 10 - 100 µl | 7056 62 | 7056 82 | 7034 91 | 7033 44 |
| 20 - 200 µl | 7056 63 | 7056 83 | 7034 91 | 7033 45 |
| 30 - 300 µl | 7056 64 | 7056 84 | 7034 91 | 7033 46 |

* incl. seal, V-ring and mounting tool.

Transferpette® S-8/-12 30-300 µl with additional pressure ring.

| Problem | Possible cause | Corrective action |
|---|--|---|
| Tip dripping (instrument leaks) | Unsuitable tip | Only use high-quality tips |
| | Tip not seated tightly | Press tip on firmly |
| The instrument does not aspirate or aspirates too little; the discharged volume is too low. | Seal contaminated | Clean seal |
| | The seal or cone is damaged | Replace seal or shaft |
| | The piston is contaminated or damaged | Clean or replace piston |
| Aspiration is too slow | Shaft clogged | Clean shaft |
| Discharged volume is too large | Pipetting button pressed too far into the blow-out position before sample uptake | Operate properly. See 'Pipetting', page 34. |
| Piston is difficult to move | The piston is contaminated or needs oil | Clean and oil the piston |

Return for repair

Important!

Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside the U.S. and Canada:

- Complete the "Declaration on Absence of Health Hazards" and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

In the U.S. and Canada:

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

BRAND GMBH + CO KG

Otto-Schott-Straße 25
97877 Wertheim (Germany)

Tel.: +49 9342 808-0
Fax: +49 9342 808-98000
E-Mail: info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)

Tel.: +1-860-767 2562
Fax: +1-860-767 2563
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park, Powai
Mumbai - 400 076 (India)

Tel.: +91 22 42957790
Fax: +91 22 42957791
E-Mail: info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)

Tel.: +86 21 6422 2318
Fax: +86 21 6422 2268
E-Mail: info@brand.cn.com
www.brand.cn.com

Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on www.brand.de or www.brandtech.com.

BRAND also offers you the possibility to have your instruments calibrated by the BRAND Calibration Service or the BRAND-owned DAkkS Calibration Service. Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND Calibration Service) or with a DAkkS Calibration Certificate. For further information, please contact your dealer or BRAND. Complete ordering information is available for download at www.brand.de (see Technical Documentation).

Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:

Information for warranty please see www.brandtech.com.

Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice.

We will not be held responsible for printing or typographical errors.

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