Liquid Handling · Easy Handling!

Transferpette® S -8/-12

F I R S T C L A S S · B R A N D

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25 Operating Manual

47 Mode d'emploi

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	Page
Safety Instructions	29
Limitations of use	30
Operating and Control Elements	31
Pipetting	32
Checking the Volume	34
Accuracy Table	36
Adjustment – Easy Calibration	37
Autoclaving	38
UV sterilization	38
Servicing and Cleaning	39
Disassembling the Manifold	41
Ordering Information	45
Accessories	46
Spare Parts	47
Troubleshooting	48
Repairs · Calibration Service	49
Warranty Information	50
Disposal	50

Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

- 1.** Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
- 2.** Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
- 3.** Observe all specifications provided by reagent manufacturers.
- 4.** Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 30).
If in doubt, contact the manufacturer or supplier.
- 5.** Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
- 6.** Avoid touching the tip orifices when working with hazardous samples.
- 7.** Never use force on the instrument!
- 8.** Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
- 9.** Before use check the instrument for visible damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 48), and contact the manufacturer if needed.

Limitations of Use

Purpose

The pipette is an air-displacement system for pipetting aqueous solutions with medium density and viscosity.

Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C to 40 °C (59 °F to 104 °F).
Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

Operating Limitations


Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than ± 1 °C/ ± 1.8 °F.

Operating Exclusions

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene and FKM
- for liquids attacking polycarbonate (viewing window)
- for liquids attacking polyphenyl sulfide
(50 µl, 100 µl, 200 µl and 300 µl models)
- for liquids of a very high vapor pressure

Note:

Optimum analysis results can only be obtained with quality tips. We recommend pipette tips from BRAND.

1. Fitting the tips

Use the correct tips according to the volume range or the color code. Ensure that the tips are securely seated.

Pipette tips are disposables items!

2. Volume setting

- Push the volume-change protection upward to disengage (UNLOCK).
- Select the desired volume by rotating the volume-setting wheel. Avoid twisting and abrupt rotating motions during this adjustment.
- Push the volume-change protection down to re-engage (LOCK). Note: The volume-change protection tightens but does not lock volume-setting wheel.


3. Align the manifold

The manifold can turn freely in both directions.

4. Aspirate sample

- Press pipetting button to the first stop.
- Hold the pipette vertically and immerse the tips 2 to 3 mm into the liquid.
- Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tips immersed for another 1-2 sec.
- Touch the tips against the container wall.


5. Discharge sample

- Place the pipette tips against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.
- Press the pipetting button slowly to the first stop and hold it down.

For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- The blow-out stroke empties the tips completely: Press the pipetting button down to the second stop.
- While doing this, wipe the pipette tips against the wall of the container.
- Remove the pipette tips from the container wall and let the pipetting button slide back.

6. Ejecting the tips

Hold the manifold over a suitable disposal container and press the tip ejection key to the stop.

Note:

ISO 8655 prescribes rinsing the pipette tips once with the sample liquid prior to the actual pipetting process.

Important!

Don't lay the instrument horizontal when the tips are filled. Liquid may enter and contaminate the instrument. The instrument should be stored without tips, placed upright in the supplied shelf/rack mount or bench top rack.

Checking the Volume

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 32 for procedure).

2. Condition the pipette

Condition the pipette before testing by using pipette tips to aspirate and discharge the test liquid (distilled H_2O) five times. After this, discard the pipette tips.

3. Carry out the test

- a) Attach new pipette tips and pre-rinse one time with test liquid.
- b) Aspirate liquid and pipette it into the weighing vessel.

Note: Each individual channel must be tested separately.

- c) Weigh the pipetted quantity with an analytical balance. (Please follow the operating manual instructions from the balance manufacturer.)
- d) Calculate the volume, taking the temperature into account.
- e) 3-10 pipettings and weighings per channel in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

Calculation (for nominal volume)

x_i = Weighing results
 n = Number of weighings

Z = Correction factor
 (e.g., 1.0029 $\mu\text{l}/\text{mg}$
 at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

Accuracy*

$$\text{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

V_0 = Nominal volume

Coefficient of Variation*

$$\text{CV\%} = \frac{100 \cdot s}{\bar{V}}$$

Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and variation coefficient (CV%):
 A% and CV% are calculated according to the formulas for
 statistical control.

Note:

Testing instructions (SOPs) and a demo version of the EASYCAL™ 4.0 calibration software are available for download at www.brand.de.

Accuracy Table

Final test values related to the nominal capacity (maximum volume) or the indicated volume steps indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/ 68 °F) and with smooth operation. According to DIN EN ISO 8655.

Transferpette® S-8/-12

Volume range µl	Volume step µl	A* ≤ ± %	CV* ≤ %	Increment µl	Type of tips µl
0.5 - 10	10	1.6	1.0	0.01	20
	5	2	2		
	1	8	6		
5 - 50	50	0.8	0.4	0.1	200
	25	1.4	0.8		
	5	6	3		
10 - 100	100	0.8	0.3	0.2	200 / 300
	50	1.4	0.6		
	10	4	2		
20 - 200	200	0.8	0.3	0.2	200 / 300
	100	1.4	0.6		
	20	4	1.5		
30 - 300	300	0.6	0.3	0.2	300
	150	1.2	0.6		
	30	3	1.5		

* A = Accuracy


* CV = Coefficient of Variation

The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made using the Easy Calibration Technique.

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2. Remove the label window and the label. Push the hook forward, raise it slightly and then pull it back.


3. Using a paperclip or an unused pipette tip, remove the protective film (this protective film can be discarded).

4. Push the red adjustment slider completely back, raise the volume-setting wheel (decoupling) and release the adjustment slider.

5. Set the previously determined actual value (adjustment value) with the volume-setting wheel in the UNLOCK position.

A volume check is recommended after every adjustment.

6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Re-insert the label and the label window.

Note:

The change to the factory settings is indicated by the red adjustment slider now visible in the label window.

Autoclaving

The Transferpette® S -8/-12 is completely autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi), with a holding time of at least 15 minutes according to DIN EN 285.

1. Eject the pipette tips.
2. Autoclave the complete pipette without any further disassembling.
3. Allow the Transferpette® S -8/-12 to completely cool and dry.

Note:

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

Attention:

Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between), with the volume-change protection set fully unlocked.

If the pipette is autoclaved frequently, the pistons should be oiled with the supplied silicone oil in order to preserve smooth movement.

If necessary after autoclaving, tighten the connection between the hand grip and the manifold.

UV sterilization


The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

The three main components of the manifold can be easily separated and disassembled for servicing, cleaning or replacing parts. The procedures are illustrated on the following pages.

Note: The changing of the V-rings/O-rings on the individual nose cones is described in detail in the instructions enclosed with the spare part.

The main components of the manifold

- A Piston unit** with piston support bar [I] and pistons inserted in this unit, which can be unscrewed individually for cleaning or replacement.
- B Nose cone assembly** with nose cone support bar [II] and central guide rod (Z), which is attached to this, and the nose cones and seals, which can be unscrewed individually for cleaning or replacement.
- C Manifold housing**, which is connected to the manifold housing cover [III] of the piston unit with two turn-lock fasteners.

➤ Continued on page 40

Servicing

In order to assure proper functioning, the Transferpette® S -8/-12 should be serviced and cleaned at regular intervals.

What is to be inspected?

1. Inspect nose cones, pistons and seals for damage and contamination.
2. Test the sealing of the instrument. We recommend using the BRAND leak testing instrument PLT unit.

Alternatively: to do this aspirate a sample, and then hold the instrument in a vertical position for about 10 sec. If a drop forms at the tip orifices, see the troubleshooting guide, page 48.

Cleaning instructions

1. Clean single nose cones, pistons and nose cone support bar/piston support bar (**these components only**) with soap solution or isopropyl alcohol. Afterwards, rinse with distilled water.
2. Let these parts dry and cool down completely. Residual moisture in the nose cones may result in a loss of accuracy.
3. Lubricate the piston with a very thin coating of the silicone oil supplied. For the central guide rod (Z) only use the recommended fluorstatic grease!

Disconnecting handle from manifold

1. Eject the pipette tips.
2. To disconnect the manifold, pull it downward as far as possible, and **only then** turn it clockwise. After one rotation, it should no longer be pulled downward while it is being turned.

Attention:




Improper handling can damage the unit!

Removing of nose cones and seals


for cleaning or replacing

1. Disconnect the manifold from the handle.
2. Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.
3. Push the mounting tool (M) (supplied with the pipette) on the nose cone and unscrew it.


► Continued on page 42

Disassembling the Manifold

4. Push the piston unit down to the bottom. The seal will either remain inside the nose cone or will stay on the piston after the nose cone is removed. Remove the seal, inspect it and clean or replace if necessary. If required, lubricate the piston lightly with the supplied silicone oil.


(In the Transferpette® S -8/-12 30-300 µl unit, remove the additional pressure ring from the piston for cleaning.)

5. Push the seal on the piston with its flat side facing upward.





(In the Transferpette® S -8/-12 30-300 µl unit, first replace the additional pressure ring!)

Use the mounting tool to mount and tighten the cleaned or new nose cone.

6. Reassemble the manifold. The manifold must be screwed into the handle counter-clockwise, until it audibly snaps into place. Check the instrument for tightness (see page 40, "Servicing").

Removing of pistons

for cleaning or replacing

1. Disconnect the manifold from the handle.
2. Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.
3. Remove both outer Phillips screws on the nose cone assembly.


Attention:

Do not loosen the central guide rod (Z)!


4. Pull the piston and nose cone assembly apart and separate. Remove stroke springs.
5. Place the mounting tool on the piston nut and unscrew the piston nut.

► Continued on page 44


Disassembling the Manifold

6. Remove the piston nut with piston spring and pull out the piston.

7. Insert the cleaned or new piston and piston spring (a). Screw on the piston nut and tighten with the mounting tool. Replace the stroke springs (b).

8. Reassemble the manifold. The manifold must be screwed counter-clockwise into the handle so that it audibly snaps into place. Check the instrument for tightness, mobility and accuracy (see page 40, "Servicing").

Items supplied

Confirm that your package includes: Transferpette®S -8 or -12 pipette, 1 Tip-Box N filled with compatible pipette tips, 1 Tip-Rack for refill, container of silicone oil, 1 shelf/rack mount, 1 reagent reservoir, 1 mounting tool, 1 set of sealing rings made of FKM and this operating manual.

Transferpette® S -8

Capacity	Description	Cat. No.
0.5 - 10 µl	M8-10	7037 00
5 - 50 µl	M8-50	7037 06
10 - 100 µl	M8-100	7037 08
20 - 200 µl	M8-200	7037 10
30 - 300 µl	M8-300	7037 12

Transferpette® S -12

Capacity	Description	Cat. No.
0.5 - 10 µl	M12-10	7037 20
5 - 50 µl	M12-50	7037 26
10 - 100 µl	M12-100	7037 28
20 - 200 µl	M12-200	7037 30
30 - 300 µl	M12-300	7037 32

Accessories

**Bench-top rack for 6 Transferpette® S
or 3 Transferpette® S-8/-12 pipettes**

Cat. No. 7048 05

**Bench-top rack for 1 Transferpette® S
-8/-12 pipette**

Cat. No. 7034 40

**Shelf/rack mount for 1 Transferpette® S
or 1 Transferpette® S-8/-12 pipette**

Cat. No. 7048 10

Additional accessories for Transferpette® S-8/-12

Description	Cat. No.
Label window , pack of 1	7046 50
Blank labels , pack of 5	7046 51
Silicone oil	7036 77
Fluorstatic grease	7036 78
Reagent reservoir , PP. Capacity 60 ml. Autoclavable (121 °C)	
non-sterile, with lid. Pack of 10.	7034 59
sterile, without lid. Packed individually. Pack of 100.	7034 11
sterile, without lid. 5 per bag. Pack of 200.	7034 09
PLT unit Pipette leak testing unit	7039 70

Parts will differ slightly depending on nominal volume of instrument.
 (Fig. shows spare parts for Transferpette® S -8/-12, 10-100 µl).

Piston Nose cone V-Ring/O-Ring

Capacity	A	B*	C	D
0.5 - 10 µl	7056 59	7056 77	7033 80	7033 40
5 - 50 µl	7056 66	7056 81	7034 91	7033 43
10 - 100 µl	7056 62	7056 82	7034 91	7033 44
20 - 200 µl	7056 63	7056 83	7034 91	7033 45
30 - 300 µl	7056 64	7056 84	7034 91	7033 46

* incl. seal, V-ring and mounting tool.

Transferpette® S -8/-12 30-300 µl with additional pressure ring.

Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (instrument leaks)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Press tip on firmly
The instrument does not aspirate or aspirates too little; the discharged volume is too low.	Seal contaminated	Clean seal
	The seal or cone is damaged	Replace seal or shaft
	The piston is contaminated or damaged	Clean or replace piston
Aspiration is too slow	Shaft clogged	Clean shaft
Discharged volume is too large	Pipetting button pressed too far into the blow-out position before sample uptake	Operate properly. See 'Pipetting', page 32.
Piston is difficult to move	The piston is contaminated or needs oil	Clean and oil the piston

Repairs · Calibration Service

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!

Return for repair

- a) Clean and decontaminate the instrument carefully.
- b) Complete the "Declaration on Absence of Health Hazards" (ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de).
- c) Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used.

The return transport of the instrument is at risk and cost of the sender.

Calibration Service

ISO 9001 and GLP guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on www.brand.de.

BRAND also offers you the possibility to have your instruments calibrated by the BRAND Calibration Service or the BRAND-owned DKD Calibration Service.

Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND Calibration Service) or with a DKD/DAkkS* Calibration Certificate. For further information, please contact your dealer or BRAND. Complete ordering information is available for download at www.brand.de.

* Based on the legal requirements the DKD Accreditation is successively transformed to the DAkkS Accreditation (Deutsche Akkreditierungsstelle GmbH), starting from January 1, 2010.

Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice.

We will not be held responsible for printing or typographical errors.

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