Gebrauchsanleitung | Operating manual | Mode
d'emploi | Instrucciones de manejo | Istruzione | 操作
手册

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Indien:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T: +91 22 42957790
F: +91 22 42957791
info@brand.co.in
www.brand.co.in

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BrandTech[®] Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T: +1-860-767 2562
F: +1-860-767 2563
info@brandtech.com
www.brandtech.com

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Table of Contents

1 Introduction.....	40	16 Calibration service.....	74
1.1 Scope of delivery.....	40	17 Warranty.....	74
1.2 Terms of use	40	18 Disposal	74
2 Safety Instructions	41		
2.1 General safety instructions.....	41		
2.2 Function.....	41		
2.3 Limitations of Use	43		
2.4 Operating Limitations.....	43		
2.5 Operating Exclusions	43		
2.6 Storage Conditions	44		
2.7 Recommended Application Range..	44		
2.8 Dispenser Selection Chart	45		
3 Functional and operating elements.....	46		
4 Assembly	48		
4.1 First Steps.....	48		
4.2 Priming	49		
5 Operation	51		
5.1 Dispensing.....	51		
5.2 Accessories	52		
6 Error limits	55		
7 Checking the Volume (Calibration).....	56		
8 Adjustment	57		
8.1 Digital model	57		
8.2 Analog-adjustable model	58		
9 Cleaning	59		
9.1 Cleaning, Analog-adjustable model, Fixed-volume model	59		
9.2 Cleaning Digital model.....	60		
9.3 Replacing discharge tube/valves	62		
10 Autoclaving.....	64		
11 Troubleshooting.....	64		
12 Product markings	65		
13 Ordering Information.....	66		
14 Accessories/spare parts	68		
15 Repairs	73		
15.1 Sending for repair	73		

1 Introduction

1.1 Scope of delivery

Bottle-top dispenser Dispensette® S or Dispensette® S Organic, discharge tube or discharge tube with recirculation valve, telescoping filling tube, recirculation tube (included only in recirculation valve models), assembly tool, various bottle adapters, a performance certificate and this operating manual.

Nominal volume [ml]	Adapters for bottle thread, PP	Filling tube Length [mm]
1, 2, 5, 10	GL 24-25, GL 28/S 28, GL 32-33, GL 38, S 40	125-240
25, 50, 100	GL 32-33, GL 38, S 40	170-330

1.2 Terms of use

- Please carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2.2 Symbols

Symbol	Meaning
	Hazardous area

1.2.3 Format

Layout	Meaning	Layout	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety Instructions

2.1 General safety instructions

Please read carefully!

The Dispensette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Comply with the operating exclusions (Operating Exclusions, p. 43)! If in doubt, contact the manufacturer or supplier.
6. Always use the instrument in such a way that neither the user nor any other person is endangered.
7. When dispensing, the discharge tube must always point away from you or any other person. Avoid splashes. Only use suitable vessels.
8. Never press down the piston when the discharge tube closure is attached.
9. Never remove the discharge tube while the dispensing cylinder is filled.
10. Reagents can accumulate in the screw cap of the discharge tube. Thus, the screw cap should be cleaned regularly.
11. For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
12. Never carry the mounted instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder may also lead to personal injury from chemicals (First Steps, p. 48, Fig. 3).
13. Never use force on the instrument. Use smooth gentle movements to operate the piston upwards and downwards.
14. Use only original manufacturer's accessories and spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
15. Always check the instrument for visible damage before use. The user can come into contact with media if the instrument has been insufficiently cleaned or inspected. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult the 'Troubleshooting' section of this manual (Troubleshooting, p. 64). Contact the manufacturer if required.

2.2 Function

With the Dispensette® S and Dispensette® S Organic bottle-top dispensers, liquids can be dispensed directly from the supply bottle. Available in Digital, Analog and Fixed models. The instruments are marked DE-M and optionally equipped with recirculation valve.

2.2.1 Dispensette S (red color code)

Digital

Analog

Fixed

2.2.2 Dispensette S Organic (yellow color code)

Digital

Analog

Fixed

2.2.3 Operation

NOTICE

Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

Dispensette® S

Borosilicate glass, Al₂O₃-ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium, PP (screw cap).

Dispensette® S Organic

Borosilicate glass, Al₂O₃-ceramic, ETFE, FEP, PFA, PTFE, tantalum, PP (screw cap).

Alternatively ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 68). The suitability of ETFE/PTFE bottle adapters must be checked by the user.

2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- Operating temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- Vapor pressure up to max. 600 mbar. Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm²/s (dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Density up to 2.2 g/cm³

2.4 Operating Limitations

- Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions). If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately (Cleaning, p. 59).
- When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
- The instrument is designed for general laboratory applications and complies with the relevant standards, e.g. DIN EN ISO 8655. Compatibility of the instrument for a specific application (e.g., trace material analysis, food sector etc.) must be checked by the user. Approvals for specific applications, e.g. for production and administration of food, pharmaceuticals or cosmetics are not available.

2.5 Operating Exclusions

2.5.1 Dispensette S

Never use Dispensette® S for:

- Liquids attacking Al₂O₃-ceramic, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide*)
- Liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Liquids which are decomposed catalytically by platinum-iridium (e.g., H₂O₂)
- Hydrochloric acid > 20 % and nitric acid > 30 %
- Tetrahydrofuran
- Trifluoroacetic acid
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument
- Liquids attacking PP (screw cap and bottle adapter)**

* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 68). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

2.5.2 Dispensette S Organic

Never use Dispensette® S Organic for:

- Liquids attacking Al_2O_3 -ceramic, tantalum, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide*)
- Liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Bases and saline solutions
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument
- Liquids attacking PP (screw cap and bottle adapter)**

* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 68). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

2.6 Storage Conditions

Store the instrument and accessories only in cleaned condition in a cool and dry place.

Storage temperature: from -20 °C to 50 °C (from -4 °F to 122 °F).

2.7 Recommended Application Range

Dispensette® S offers a broad range of applications for dispensing aggressive reagents, including concentrated acids (e.g., H_3PO_4 , H_2SO_4), bases (e.g., NaOH, KOH), saline solutions, as well as many polar solvents.

Dispensette® S Organic is ideal for dispensing organic solvents including chlorinated and fluorinated hydrocarbons (e.g., trichlorotrifluoroethane and dichloromethane), concentrated acids (e.g., HCl and HNO_3), trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides.

NOTICE

Instrument selection

For guidelines on selecting the right dispenser observe the corresponding Operating Exclusions and the 'Dispenser selection chart' on the next page.

NOTICE

Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

2.8 Dispenser Selection Chart

Reagent	Disp. S	Disp. S Organ	Reagent	Disp. S	Disp. S Organ	Reagent	Disp. S	Disp. S Organ
Acetaldehyde	+	+	Cyclohexane			Methylene chloride		+
Acetic acid (glacial), 100%	+	+	Cyclohexanone	+	+	Mineral oil (Engine oil)	+	+
Acetic acid, ≤ 96%	+	+	Cyclopentane		+	Monochloroacetic acid		+
Acetic anhydride		+	Decane	+	+	Nitric acid, ≤ 30%	+	+
Acetone	+	+	1-Decanol	+	+	Nitric acid, 30-70% ***		+
Acetonitrile	+	+	Dibenzyl ether	+	+	Nitrobenzene	+	+
Acetophenone		+	Dichloroacetic acid		+	Oleic acid	+	+
Acetyl chloride		+	Dichlorobenzene	+	+	Oxalic acid	+	+
Acetylacetone	+	+	Dichloroethane		+	n-Pentane		+
Acrylic acid	+	+	Dichloroethylene		+	Peracetic acid		+
Acrylonitrile	+	+	Dichloromethane		+	Perchloric acid	+	+
Adipic acid	+		Diesel oil (Heating oil), bp 250-350 °C		+	Perchloroethylene		+
Allyl alcohol	+	+	Diethanolamine	+	+	Petroleum, bp 180-220 °C		+
Aluminium chloride	+		Diethyl ether		+	Petroleum ether, bp 40-70 °C		+
Amino acids	+		Diethylamine		+	Phenol	+	+
Ammonia, ≤ 20%	+	+	Diethyleneglycol	+	+	Phenylethanol	+	+
Ammonia, 20-30%		+	1,2-Diethylbenzene	+	+	Phenyldiazine	+	+
Ammonium chloride	+		Diethylene glycol	+	+	Phosphoric acid, ≤ 85%	+	+
Ammonium fluoride	+		Dimethyl sulfoxide (DMSO)	+	+	Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1	+	+
Ammonium sulfate	+		Dimethylaniline		+	Piperidine	+	+
n-Amyl acetate	+	+	Dimethylformamide (DMF)	+	+	Potassium chloride	+	
Amyl alcohol (Pentanol)	+	+	1,4-Dioxane		+	Potassium dichromate	+	
Amyl chloride (Chloropentane)	+		Diphenyl ether	+	+	Potassium hydroxide	+	
Aniline	+	+	Essential oil		+	Potassium permanganate	+	
Barium chloride	+		Ethanol	+	+	Propionic acid	+	+
Benzaldehyde	+	+	Ethanolamine		+	Propylene glycol (Propanediol)	+	+
Benzene (Benzol)	+	+	Ethyl acetate		+	Pyridine	+	+
Benzene (Petroleum benzene) bp 70-180 °C		+	Ethylbenzene		+	Pyruvic acid	+	+
Benzoyl chloride	+	+	Ethylene chloride	+	+	Salicylaldehyde	+	+
Benzyl alcohol	+	+	Fluoroacetic acid		+	Scintillation fluid	+	+
Benzylamine	+	+	Formaldehyde, ≤ 40%	+	+	Silver acetate	+	
Benzylchloride	+	+	Formamide	+	+	Silver nitrate	+	
Boric acid, ≤ 10%	+	+	Formic acid, ≤ 100%		+	Sodium acetate	+	
Bromobenzene	+	+	Glycerol	+	+	Sodium chloride	+	
Bromonaphthalene	+	+	Glycol (Ethylene glycol)	+	+	Sodium dichromate	+	
Butanediol	+	+	Glycolic acid, ≤ 50%	+	+	Sodium fluoride	+	
1-Butanol	+	+	Heating oil (Diesel oil), bp 250-350 °C		+	Sodium hydroxide, ≤ 30%	+	
n-Butyl acetate	+	+	Heptane		+	Sodium hypochlorite	+	
Butyl methyl ether	+	+	Hexane		+	Sulfuric acid, ≤ 98%	+	+
Butylamine	+	+	Hexanoic acid	+	+	Tartaric acid	+	
Butyric acid	+	+	Hexanol	+	+	Tetrachloroethylene		+
Calcium carbonate	+		Hydriodic acid, ≤ 57% **	+	+	Tetrahydrofuran (THF) ***		+
Calcium chloride	+		Hydrobromic acid		+	Tetramethylammonium hydroxide	+	
Calcium hydroxide	+		Hydrochloric acid, ≤ 20%	+	+	Toluene		+
Calcium hypochlorite	+		Hydrochloric acid, 20-37% **		+	Trichloroacetic acid	+	
Carbon tetrachloride		+	Hydrogen peroxide, ≤ 35%		+	Trichlorobenzene	+	
Chloro naphthalene	+	+	Isomyl alcohol	+	+	Trichloroethane		+
Chloroacetaldehyde, ≤ 45%	+	+	Isobutanol	+	+	Trichloroethylene		+
Chloroacetic acid	+	+	Isocetane		+	Trichlorotrifluoro ethane		+
Chloroacetone	+	+	Isopropanol (2-Propanol)	+	+	Triethanolamine	+	+
Chlorobenzene	+	+	Isopropyl ether	+	+	Triethylene glycol	+	+
Chlorobutane	+	+	Lactic acid	+	+	Trifluoro ethane		+
Chloroform		+	Methanol	+	+	Trifluoroacetic acid (TFA)		+
Chlorosulfonic acid	+	+	Methoxybenzene	+	+	Turpentine		+
Chromic acid, ≤ 50%	+	+	Methyl benzoate	+	+	Urea	+	
Chromosulfuric acid	+	+	Methyl butyl ether	+	+	Xylene		+
Copper sulfate	+		Methyl ethyl ketone	+	+	Zinc chloride, ≤ 10%	+	
Cresol		+	Methyl formate	+	+	Zinc sulfate, ≤ 10%	+	
Cumene (Isopropyl benzene)	+	+	Methyl propyl ketone	+	+			

* use ETFE/PTFE bottle adapter

** use PTFE seal

For dispensing hydrofluoric acid, we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BRAND. Status as of: 02/19/13


3 Functional and operating elements

Dispensette® S Digital

- 1 Volume adjustment
- 2 Cover plate
- 3 Housing
- 4 Protective sleeve/dispensing cylinder
- 5 Piston
- 6 Valve block
- 7 Valve block adapter (GL 45 bottle thread)
- 8 Recirculation tube
- 9 Telescoping filling tube
- 10 Latch
- 11 Safety lock
- 12 Piston mounting nut
- 13 Knob
- 14 Screw cap
- 15 Discharge tube
- 16 Valve lever, recirculation valve

Dispensette® S Analog

- 1 Cover plate
- 2 Piston seat
- 3 Housing
- 4 Pointer
- 5 Volume adjustment
- 6 Protective sleeve/dispensing cylinder
- 7 Valve block
- 8 Valve block adapter (GL 45 bottle thread)
- 9 Recirculation tube
- 10 Telescoping filling tube
- 11 Screw cap
- 12 Discharge tube
- 13 Valve lever, recirculation valve
- 14 Piston

Telescoping filling tube and recirculation tube

Assembly tool

4 Assembly

4.1 First Steps


⚠ WARNING

Follow the safety instructions

- > Wear protective clothing, eye protection and protective gloves!
- > Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- > Follow all safety instructions and comply with the limitations of use, see Limitations of Use, p. 43.
- > Comply with the operating limitations, see Operating limitations, p. 43.

1. Mounting the filling tube/recirculation tube

- a. Adjust the length of the telescoping filling tube to the bottle height and attach it.
- b. Center and attach the filling tube (part with smaller diameter) carefully to avoid damaging the nozzle.
- ⇒ If a discharge tube with a recirculation valve is used, the recirculation tube must also be installed.
- c. Insert the recirculation tube with the opening pointing outward.

2. Mount the instrument on a bottle and align

NOTICE

To avoid tipping over, use a bottle stand for small bottles.

- a. Screw the instrument (GL 45 threads) onto the reagent bottle, and then align the discharge tube with the bottle label. This is done by rotating the valve block with the discharge tube.

3. Transporting the instrument

NOTICE

Polypropylene (PP) bottle adapters are supplied with the instrument. These may be used only for media that do not attack the PP. Alternatively ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 68). The suitability of ETFE/PTFE bottle adapters must be checked by the user.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- c. When mounted to a reagent bottle, always carry the instrument as shown in the figure!

4.2 Priming

⚠ WARNING

Follow during every use, especially with hazardous media

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 43 and Operating exclusions, p. 43.

NOTICE

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splashes, as ethanol and glycerin residues may be present in the instrument.

Instruments with recirculation valve

- a. Open the screw cap of the dispensing tube.

- b. Turn valve to 'Recirculate'.

- c. For deaerating gently pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this process at least 5 times.

- d. Turn valve to 'Dispense'.

- e. To avoid splashes when deaerating hold the discharge tube on the inner wall of a suitable receiving vessel and dispense liquid to vent deaerate the discharge tube until it is bubble-free. Wipe away any remaining drops from the discharge tube.

Instruments without recirculation valve

- Open the screw cap of the discharge tube (see 'instrument with recirculation valve', Fig. a.). To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel.
- For deaerating pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until the discharge tube is bubble-free.

5 Operation

5.1 Dispensing

1. Selecting the volume

Digital: Rotate the volume-setting wheel until the desired volume is indicated (mechanical counter).

Analog: Loosen the volume selector thumb screw $\frac{3}{4}$ turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

Fixed: The volume is non-adjustable and cannot be changed.

2. Dispensing

⚠ WARNING

Follow during every use, especially with hazardous media

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 43 and Operating exclusions, p. 43.

- a. Unscrew the screw cap of the dispensing tube.
- b. When using instruments equipped with a recirculation valve, turn the valve to 'Dispensing'.
- c. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- d. Gently lift the piston until the upper stop and then depress piston slowly and steadily with minimal force until the lower stop.
- e. Wipe off the discharge tube against the inner wall of the receiving vessel.
- f. Seal the dispensing tube with the screw cap.

NOTICE

After using the piston, always press it down to the lower stop (parking position). If the piston is not pressed down to the lower stop, unintentional media leaks may occur.

NOTICE

The filled status of the instrument must be specially marked during cleaning.

5.2 Accessories

5.2.1 Flexible discharge tube with recirculation valve

For serial dispensing the flexible discharge tube can be used for the bottle-top dispenser Dispensette® S and Dispensette® S Organic (Accessories/spare parts, p. 68).

The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, the tube must lie in regular loops and must not be twisted. The operating exclusions for the respective instrument in use apply.

Assembly

⚠ WARNING

Use only undamaged tubing


There should be no visible damage to the discharge tube (e.g. kinks, etc.). Each time you are going to use the tubing, examine it carefully.

- To dispense aggressive liquids, you should take safety measures in addition to the normal precautions. We recommend use of a protective shield.
- The bottle must be supported using a bottle stand.
- To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace into the holder after use.
- For cleaning rinse the tube carefully.
- Do not dismantle!

Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, Cleaning, p. 59).
- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.

- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions (Fig. 1).
- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it (Fig. 2). In doing so, the instrument must not be mounted to a bottle. Install the receiver tube.
- d. Press the plug of the recirculation valve downwards.
- e. Slide the flexible discharge tube housing into the valve block up to the stop (Fig. 3).
- f. Slide the discharge tube housing all the way down (Fig. 4).
- g. Attach the valve handle that fits the discharge valve and press it in firmly. Note the color coding and marking (see assembly instruction 'Flexible discharge tube for Dispensette® S').

NOTICE

Use a bottle holder (Accessories/spare parts, p. 68).

5.2.2 Drying tube

Use of a drying tube, filled with a suitable absorbent (purchased separately), might be necessary for moisture- and CO_2 -sensitive media.

(Accessories/spare parts, p. 68)

Assembly

- a. Use the assembly tool to unscrew the air vent cap.
- b. Screw in the filled drying tube.
- c. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

NOTICE

If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

5.2.3 Sealing ring for valve block

For highly volatile media we recommend to seal the connection from valve block to bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts, p. 68).

Assembly

- a.** Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

5.2.4 Air vent cap for micro filter with Luer-cone

For sterile media we recommend the air vent cap with Luer-cone to attach a micro filter. This provides increased protection against contamination by displacement air (Accessories/spare parts, p. 68).


Assembly

- a.** Unscrew the air vent cap (see 'Assembly Drying tube', Fig. a.).
- b.** Screw in the air vent cap with a Luer cone.
- c.** Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

- d.** Insert a commercially available sterile filter into the Luer cone.

Type Digital • Easy Calibration is manufactured under U.S. Patent 5,957,330.

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. Testing takes place according DIN EN ISO 8655-6 with a completely filled instrument and with uniform and smooth dispensing.

Error limits

Nominal volume ml	A* ≤ ± %	µl	CV* ≤ %	µl
1	0.6	6	0.2	2
2	0.5	10	0.1	2
5	0.5	25	0.1	5
10	0.5	50	0.1	10
25	0.5	125	0.1	25
50	0.5	250	0.1	50
100	0.5	500	0.1	100

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

Partial volumes

The percentage values for A and CV are relative to the nominal volume (V_N) and must be converted for partial volumes (V_p).

e.g.	Volume	A* ≤ ± %	µl	CV* ≤ %	µl
V_N	25.0	0.5	125	0.1	25
$V_T = 50\% N$	12.5	1.0	125	0.2	25
$V_T = 10\% N$	2.5	5.0	125	1.0	25

*A = Accuracy, CV = Coefficient of Variation

NOTICE

The error limits in DIN EN ISO 8655-5 are significantly lower. The maximum error for a single measurement can be approximated from the sum of error limits $EL = A + 2 \times CV$ (for 25 ml: $125 \mu l + 2 \times 25 \mu l = 175 \mu l$).

7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from www.brand.de.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits, p. 55) is performed as follows:

1. Prepare the instrument

Clean the instrument (Cleaning, p. 59), fill it with distilled H₂O and then prime it carefully.

2. Check the volume

- 10 dispensing operations with distilled H₂O in 3 volume ranges (100 %, 50 %, 10 %) are recommended
- For discharge, depress the piston slowly and steadily without force until the lower stop
- Wipe off the tip of discharge tube.
- Weigh the dispensed quantity on an analytical balance. (Please refer to the user manual of the scale manufacturer.)
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

3. Calculation

x_i = Weighing results n = Number of weighings

Z = Correction factor
(e.g. 1.0029 µl/mg at 20 °C, 1013 hPa)

$$\text{Mean} \quad \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume} \quad \bar{V} = \bar{x} \cdot z$$

$$\text{Accuracy*} \quad A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

$$V_0 = \text{Nominal volume}$$

$$CV\% = \frac{\bar{V}}{\bar{V}}$$

$$\text{Standard deviation*} \quad s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

8 Adjustment

After a long period of usage, an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (Checking the Volume (Calibration), p. 56).
- Calculate mean volume (actual value) (Checking the Volume (Calibration), p. 56).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm the appropriate adjustment.

Example

The gravimetric check gives an actual value of 9.90 ml for a 10 ml instrument set for a nominal volume of 10.00 ml.

8.1 Digital model

- a.** Slide the latch to the left and remove the front housing shell.

- b.** Pull out the safety lock. The adjustment cover will then come off. Discard the adjustment cover.

- c.** Pull the red knob to disengage the gears. Set the display to actual delivered volume (e.g., 9.90 ml).

- d.** Press in the red knob and then the safety lock again.


- e.** Replace housing and slide the latch to the right.
 ⇒ Adjustment is completed.
 ⇒ Alteration of the factory setting is indicated by a red recalibration flag (circle).

Adjustment range

Nominal volume [ml]	Digital max. +/- [μ l]
1	12
2	24
5	60
10	120
25	300
50	600
100	–

8.2 Analog-adjustable model

- a. Insert the pin of the assembly tool into the cover plate (Pos. 1), and break it off with a rotating motion. Discard the adjustment cover.
 - b. Insert the pin of the assembly tool into the adjustment screw and rotate to the left in order to increase the dispensing volume, or rotate to the right to decrease the dispensing volume (e.g. for an actual value of 9.97 ml, rotate approx. 1/2 turn to the left).
- ⇒ Adjustment is completed
- ⇒ The change in the adjustment is indicated by a red disk (circle in figure).

Adjustment range

Nominal volume	Analog-adjustable/Fixed-volume max. +/- [μ l]	One rotation corresponds to [μ l]
1 ml	6	~ 8
2 ml	12	~ 16
5 ml	30	~ 40
10 ml	60	~ 80
25 ml	150	~ 130
50 ml	300	~ 265
100 ml	600	~ 400

9 Cleaning

⚠ WARNING

Components filled with reagent

The cylinder, valves, telescoping filling tube and discharge tube contain reagent!

- > Never remove the discharge tube while the dispensing cylinder is filled.
- > Point the valves and tube openings away from your body.
- > Wear protective clothing, eye protection and appropriate hand protection!

The instrument must be cleaned in the following situations to assure correct operation:

- immediately when the piston is difficult to move
- prior to a reagent change
- prior to long term storage
- prior to dismantling the instrument
- prior to autoclaving
- prior to changing the valve
- regularly when using liquids which form deposits (e.g., crystallizing liquids)
- regularly when liquids accumulate in the screw cap.

For proper cleaning and removal of any deposits in the parts through which liquids pass, also always completely withdraw the piston from the cylinder after rinsing with a suitable cleaning solution. If necessary, the parts can also be cleaned in an ultrasonic bath.

9.1 Cleaning, Analog-adjustable model, Fixed-volume model

1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

2. Rinsing the instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

3. Removing pistons

NOTICE

Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- a. Hold the housing securely and unscrew the piston seat completely by turning it to the left.

- b. Pull out the piston carefully.

4. Cleaning the piston and cylinder and reassembling the instrument

NOTICE

Special instructions for Dispensette® S Organic

For the Dispensette® S Organic, slide the piston into the cylinder vertically with a rotating motion.

- a. Clean the piston and cylinder and reassemble the instrument. If necessary carefully remove deposits at the edge of the glass cylinder.
- b. Rinse the piston and cylinder with deionized water, and dry them carefully.
- c. Insert the piston completely into the cylinder and then reassemble the instrument. The piston can only be inserted into the cylinder when the discharge tube is attached.

9.2 Cleaning Digital model

1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

2. Rinsing the instrument


- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

3. Removing pistons

NOTICE

Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- a. Carry out assembly and disassembly at the maximum volume setting only. To do this, turn the volume adjustment wheel as far as it will go in + direction (= maximum volume).
- b. Slide the latch to the left.
- c. Remove the front housing shell.
- d. Use the assembly tool to loosen the piston mounting nut and pull the rear housing section backwards.
- e. Then carefully pull the piston out of the cylinder.

4. Cleaning the piston and cylinder and reassembling the instrument



NOTICE

Special instructions for Dispensette® S Organic

For the Dispensette® S Organic, slide the piston into the cylinder vertically with a rotating motion.

- a. Clean the piston and cylinder and reassemble the instrument. If necessary carefully remove deposits at the edge of the glass cylinder.
- b. Rinse the piston and cylinder with deionized water, and dry them carefully.
- c. Insert the piston completely into the cylinder and then reassemble the instrument. The piston can only be inserted into the cylinder when the discharge tube is attached.

NOTICE**Cleaning Digital model**

The stop segment must engage underneath the cylinder's stop ring. When tightening the piston mounting nut using the assembly tool, use your thumbs to press the piston/cylinder unit back towards the back side of the housing.

9.3 Replacing discharge tube/valves

NOTICE

A functional test must always be carried out after exchanging components.

9.3.1 Replacing discharge tube

- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.
- c. Hold coupling piece of the new discharge tube and pull housing up. Push housing into the valve block until it meets the stop.
- d. Slide the discharge tube housing all the way down.
- e. For instruments with a recirculation valve, pull up the valve lever to the 'Recirculate' position, and press it in tightly.

9.3.2 Replacing valves

**NOTICE**

Always install the valve intended for the particular instrument model and size! (See ordering information, Ordering Information, p. 66) Identical filling valves are used in the Dispensette® S and Dispensette® S Organic, but the discharge valves are different. For easy identification, the discharge valves for the Dispensette® S Organic are marked with 'ORG'.

9.3.2.1 Replacing the discharge valve

NOTICE

The discharge valve has a safety ball that closes automatically when the discharge tube is removed. Attaching the discharge tube reopens the safety ball.

- a. After disassembling the discharge tube (see Replacing discharge tube, p. 62), use the assembly tool to unscrew the discharge valve.
- b. Firmly screw in the new discharge valve first by hand and then tighten it with the assembly tool. The thread should not be visible.

9.3.2.2 Replacing the filling valve

- a. Pull out the recirculation tube and the telescoping filling tube.
- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

9.3.2.3 Dislodging a stick valve ball

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 µl plastic pipette tip (see the figure at the side).

10 Autoclaving

The instrument can be autoclaved at 121 °C (250 °F), 2 bar with a hold time of at least 15 minutes according to DIN EN 285.

- a. The instrument must be carefully cleaned prior to autoclaving (Cleaning, p. 59).
- b. Open the closure cap on the discharge tube, and for instruments with a recirculation valve, set the valve to 'Dispense'.
- c. Check that the filling valve is securely seated (Fig. 1a). In the Digital model, additionally check that the piston mounting is securely seated (Fig. 1b).
- d. To ensure unhindered access for the steam and to prevent the valve ball in the filling valve from possibly becoming stuck, hold the instrument with the discharge piston pressed vertically downward, and gently tap against the casing with your hand (Fig. 2). Then lay it horizontally in the autoclave. Be sure to avoid the instrument coming into contact with metal surfaces in the autoclave!

NOTICE


Do not reassemble the instrument until it has cooled down to room temperature (Cooling time approx. 2 hours). After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. It is the user's responsibility to ensure effective autoclaving.

11 Troubleshooting

Problem	Possible cause	Corrective action
Piston sluggish or stuck	Crystal deposits, Contaminants	Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Carry out cleaning (Cleaning, p. 59).
Filling not possible	Volume adjusted to minimum setting	Set the desired volume (Dispensing, p. 51).
	Filling valve stuck	Unscrew the filling valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Dislodging a stick valve ball, p. 63). If necessary, replace the filling valve.
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block, clean it and replace if necessary. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it.
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently	Discharge valve is not screwed in deeply enough	Tighten the discharge valve with the assembly tool until it meets the stop so that the threads are no longer visible.

Problem	Possible cause	Corrective action
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Valve screw connections loose	Tighten the valves firmly with the assembly tool.
	The instrument has not been primed	Prime the instrument (Priming, p. 49).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Carry out cleaning (Cleaning, p. 59). Tighten the valves using the assembly tool.
Dispensed volume too low	Filling tube is loose or damaged	Carry out cleaning (Cleaning, p. 59). Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube (Replacing discharge tube, p. 62).
	Filling valve is loose or damaged	Carry out cleaning (Cleaning, p. 59). Tighten the valves using the assembly tool. If necessary, replace filling valves.
Leaking liquid between instrument and bottle	Recirculation tube not connected	Connect recirculation tube (First Steps, p. 48, Fig. 3).
	Volatile reagent dispensed without sealing ring	Install sealing ring for valve block (First Steps, p. 48).

12 Product markings

Symbol or number	Meaning
	General warning symbol
	Refer to the operating manual
	Use eye protection
	Use hand protection
	Wear protective clothing
XXZXXXXX	Serial number
 18	The device is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the last two digits of the year in which the marking was affixed (here: 2018).
www.brand.de/ip	Patent information

13 Ordering Information

Dispensette® S, Digital

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
0.1 – 1	0.005	4600 310	4600 311
0.2 – 2	0.01	4600 320	4600 321
0.5 – 5	0.02	4600 330	4600 331
1 – 10	0.05	4600 340	4600 341
2.5 - 25	0.1	4600 350	4600 351
5 – 50	0.2	4600 360	4600 361

Dispensette® S, Analog-adjustable

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
0.1 - 1	0.02	4600 100	4600 101
0.2 - 2	0.05	4600 120	4600 121
0.5 - 5	0.1	4600 130	4600 131
1 - 10	0.2	4600 140	4600 141
2.5 - 25	0.5	4600 150	4600 151
5 - 50	1.0	4600 160	4600 161
10 - 100	1.0	4600 170	4600 171

Dispensette® S, Fixed-volume

Volume [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
1	4600 210	4600 211
2	4600 220	4600 221
5	4600 230	4600 231
10	4600 240	4600 241
Special fixed volumes: 0.5-100 ml	4600 290	4600 291

Volume [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
(Please specify when ordering!)		

NOTICE

Scope of supply, see Scope of delivery, p. 40

Dispensette® S Organic, Digital

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
0.5 - 5	0.02	4630 330	4630 331
1 - 10	0.05	4630 340	4630 341
2.5 - 25	0.1	4630 350	4630 351
5 - 50	0.2	4630 360	4630 361

Dispensette® S Organic, Analog-adjustable

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
0.5 - 5	0.1	4630 130	4630 131
1 - 10	0.2	4630 140	4630 141
2.5 - 25	0.5	4630 150	4630 151
5 - 50	1.0	4630 160	4630 161
10 - 100	1.0	4630 170	4630 171

Dispensette® S Organic, Fixed-volume

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
5		4630 230	4630 231
10		4630 240	4630 241
Special fixed volumes: 2-100 ml		4630 290	4630 291

Volume [ml]	Graduation [ml]	Discharge tube without recirculation valve Order No.	Discharge tube with recirculation valve Order No.
(Please specify when ordering!)			

NOTICE

Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

14 Accessories/spare parts

Bottle adapters

External thread	for bottle thread/ground joint	Material	Order No.
GL 32	GL 25	PP	7043 25
GL 32	GL 28/ S 28	PP	7043 28
GL 32	GL 30	PP	7043 30
GL 32	GL 45	PP	7043 45
GL 45	GL 32	PP	7043 96
GL 45	GL 35	PP	7044 31
GL 45	GL 38	PP	7043 97
GL 45	S* 40	PP	7043 43
GL 45	S* 50	PP	7043 50
GL 45	S* 54	PP	7044 30
GL 45	S* 60	PP	7043 48
GL 32	GL 25	ETFE	7043 75
GL 32	GL 28/ S 28	ETFE	7043 78
GL 32	GL 45	ETFE	7043 95
GL 45	GL 32	ETFE	7043 98
GL 45	GL 38	ETFE	7043 99
GL 45	S* 40	PTFE	7043 91
GL 32	NS 19/26	PP	7044 19
GL 32	NS 24/29	PP	7044 24
GL 32	NS 29/32	PP	7044 29

* buttress thread

Discharge tubes without recirculation valve

Packaged unit 1 pc

Description	Nominal volume ml	Shape	Length mm	Order No.
for Dispensette® S	1, 2, 5, 10	fine tip	108	7080 02
	5, 10	standard	108	7080 05
	25, 50, 100	fine tip	135	7080 06
	25, 50, 100	standard	135	7080 08
for Dispensette® S Organic	5, 10	fine tip	108	7080 12
	5, 10	standard	108	7080 14
	25, 50, 100	fine tip	135	7080 16
	25, 50, 100	standard	135	7080 19

Discharge tubes with recirculation valve

Packaged unit 1 pc

Description	Nominal volume ml	Shape	Length mm	Order No.
for Dispensette® S	1, 2, 5, 10	fine tip	108	7081 02
	5, 10	standard	108	7081 04
	25, 50, 100	fine tip	135	7081 06
	25, 50, 100	standard	135	7081 09
for Dispensette® S Organic	5, 10	fine tip	108	7081 12
	5, 10	standard	108	7081 14
	25, 50, 100	fine tip	135	7081 16
	25, 50, 100	standard	135	7081 19

Flexible discharge tube with recirculation valve

for Dispensette® S and Dispensette® S Organic:
PTFE, spiraled, approx. 800 mm long, with safety handle.
Packaged unit 1 pc
Not suitable for hydrofluoric acid (HF)

Nominal volume ml	Discharge tube Outer diameter	Discharge tube Inner diameter	Order No.
1, 2, 5, 10	3	2	7081 32
25, 50, 100	4.5	3	7081 34

Dispensette® S discharge valve

PFA/Boro 3.3/ceramic/platinum-iridium.
No valve marking.
Packaged unit 1 pc

For nominal volume ml	Order No.
1, 2*	6749
5, 10	6727
25, 50, 100	6728

* with valve marking '1 + 2'

Dispensette® S Organic discharge valve

PFA/Boro 3.3/ceramic/tantalum.
Valve marked with 'ORG'.
Packaged unit 1 pc

For nominal volume ml	Order No.
5, 10	6729
25, 50, 100	6730

Dispensette® S and Dispensette® S Organic filling valve


PFA/ETFE/Boro 3.3/ceramic.

No valve marking.

Packaged unit 1 pc

For nominal volume ml	Order No.
1, 2, 5, 10	6734
25, 50, 100	6735

Telescoping filling tube for Dispensette® S and Dispensette® S Organic

FEP.

Individually adjustable lengths.

Packaged unit 1 pc

For nominal volume ml	Outer diameter mm	Length mm	Order No.
1, 2, 5, 10	6	70-140	7082 10
		125-240	7082 12
		195-350	7082 14
		250-480	7082 16
25, 50, 100	7.6	170-330	7082 18
		250-480	7082 20

Air vent cap for micro filter with Luer-cone

PP. Air vent cap and PTFE-sealing ring. Packaged unit 1 pc
Order No. 7044 95

Recirculation tube

FEP
Packaged unit 1 pc
Order No. 6747

Adjustment/assembly tool


Packaged unit 1 pc
Order No. 6748

Bottle stand

PP. Support rod, 325 mm, base plate 220 x 160 mm.
Packaged unit 1 pc
Order No. 7042 75

Sealing ring for valve block

PTFE, for highly volatile media.
Packaged unit 1 pc
Order No. 7044 86

Drying tube


Incl. PTFE sealing ring (without drying agent). Packaged unit 1 pc
Order No. 7079 30

Screw cap with fastener

Packaged unit 1 pc

Description	Nominal volume ml	Order No.
PP, red, for Dispensette® S	1, 2, 5, 10	7060 18
	25, 50, 100	7060 19
PP, yellow, for Dispensette® S Organic	5, 10	7060 25
	25, 50, 100	7060 27
ETFE*, for Dispensette® S and Dispensette® S Organic	1, 2, 5, 10	7060 29
PTFE*, for Dispensette® S and Dispensette® S Organic	25, 50, 100	7060 31

Dispensette® cask withdrawal system

Packaged unit 1 pc

Description	Order No.
Dispensette® cask withdrawal system	7042 61
Tripod connection for wall mount, stainless steel, Dispensette® withdrawal system	7042 68
Table/rack clamp for wall mount for Dispensette® withdrawal system	7042 72

15 Repairs

15.1 Sending for repair

NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

Clean the device thoroughly and decontaminate!

- When returning products, please include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside the U.S. and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within the U.S. and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc. **before** sending the device in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

BRAND GMBH + CO KG
Otto-Schott-Strasse 25
97877 Wertheim (Germany)
T: +49 9342 808-0
F: +49 9342 808-98000
info@brand.de
www.brand.de

India:
BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T: +91 22 42957790
F: +91 22 42957791
info@brand.co.in
www.brand.co.in

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T: +1-860-767 2562
F: +1-860-767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)
T: +86 21 6422 2318
F: +86 21 6422 2268
info@brand.cn.com
www.brand.cn.com

16 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring instruments. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media. The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or by the BRAND DAkkS laboratory.

Simply send us the instrument to be calibrated, accompanied by details about which type of calibration you would like. The instrument will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from BRAND.

The order document is available for download on www.brand.de (see Technical Documents).

17 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

18 Disposal

Prior to disposal, please observe the corresponding national disposal regulations and ensure that the product is disposed of properly.

Subject to technical changes, errors, and misprints.

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