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Table of contents

1	Introduction.....	31
1.1	Scope of delivery	31
1.2	Terms of use	31
2	Safety regulations	33
2.1	Safety regulations	33
2.2	Purpose.....	34
2.3	Limitations of use.....	34
2.4	Application restrictions.....	34
2.5	Operating exclusions	35
3	Functional and operational components	36
4	Pipetting.....	37
5	Cleaning	39
6	Maintenance	40
6.1	Change capillary	40
6.2	Change piston rod	42
6.3	Change seal	45
7	Checking the volume	46
8	Troubleshooting.....	48
9	Accuracy table.....	49
10	Ordering data for instruments and spare parts	50
10.1	Accessories	50
11	Product markings	51
12	Repairs	51
12.1	Sending for repair	51
13	Calibration service.....	53

14 Warranty	54
15 Disposal.....	54

1 Introduction

1.1 Scope of delivery

Transferpettor, quality certificate and replacement capillaries.

1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2.2 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety regulations

2.1 Safety regulations

Please read carefully!

The instrument micro can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and observe this operating manual before using the instrument.
2. Follow the general hazard instructions and safety regulations (e.g., wear protective clothing, eye protection and protective gloves). When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
3. Follow the instructions given by the reagent manufacturer.
4. Use the instrument only for pipetting liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see Operating exclusions. In case of doubt, contact the manufacturer or dealer.
5. Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels.
6. Avoid touching the tip opening when working with aggressive media.
7. Never use force.
8. Use only original spare parts. Do not make any technical modifications. Do not disassemble the instrument further than described in the user manual.

9. Always check that the instrument is in proper working condition before use. If the instrument appears to have a malfunction (e.g., sluggish piston, leakage), stop pipetting immediately and observe Troubleshooting. Contact the manufacturer, if necessary.

2.2 Purpose

The instrument is used for pipetting liquids. It works on the principle of direct displacement.

2.3 Limitations of use

Suitable for media with:

- densities of up to 13.6 g/cm^3
- vapor pressure up to 500 mbar
- viscosity (see table)

Nominal volume [μl]	suitable for viscosities of up to (kinematic viscosity) [mm^2/s]
1, 2, 5	6,000
10	20,000
20	50,000
25	50,000
50	100,000
100	140,000

Working temperature range: $15\text{--}40^\circ\text{C}$ ($59\text{--}104^\circ\text{F}$)

2.4 Application restrictions

Crystallization of concentrated salt solutions and other crystallizing liquids can damage the seal.

2.5 Operating exclusions

NOTICE


Do not autoclave

The instrument is not autoclavable.

The user is responsible for checking the compatibility of the instrument with the intended application. The instrument cannot be used:

- For liquids that attack glass, stainless steel (instruments up to 10 µl), or PTFE (instruments from 20 µl).

3 Functional and operational components

1 Pipetting button

3 Volume display (digital instrument)

5 Allen screws (under the housing)

7 Union nut

9 Slotted screw for housing mounting

11 Piston rod with seal (from 20 µl)

2 Lever for locking the volume setting (digital instrument)

4 Housing (here: fixed volume)

6 Shaft with inspection window
Union nut


8 End plate

10 Capillary

12 Piston rod without seal (up to 10 µl)

4 Pipetting

Setting the volume

- a. Set lever (1) to 
- b. Turn the pipetting button (2) until the desired volume appears in the display (3).
- c. Set lever (1) to 

Recommended working range:


20–100% of the nominal volume. For volumes below 20 %, accuracy is no longer guaranteed.

Aspirating liquid

- a. Press the pipetting button as far as it will go.

b.

b. Immerse the capillary in the liquid (2–3 mm).


c. Allow the pipetting button to slowly slide back.

d.

d. Gently wipe capillary along vessel wall.

e.

e. Carefully wipe the outside of the capillary with a lint-free cloth or swab. Do not touch the opening; otherwise, liquid will be sucked out.


Eject liquid

a.


a. Place capillary against vessel wall.

b. Press the pipetting button as far as it will go.

- c. Gently wipe capillary along vessel wall.
- d. Allow the pipetting button to slide back.

Check for leaks daily

- a. Set the nominal volume.
- b. Aspirate the liquid.
- c. Hold the instrument vertically for 10 s (capillary downwards).
- d. If a drop forms (see Troubleshooting, p. 48)

5 Cleaning

- a. Disassemble the instrument (see Maintenance, p. 40).
- b. Clean the piston rod.
- c. Allow to dry completely.
- d. Assemble the instrument.

6 Maintenance

6.1 Change capillary


Remove old capillary

CAUTION

Risk of contamination

Do not touch the contaminated area of the capillaries. Wear protective gloves.

- Loosen the union nut (half turn).
- Pull the capillary off the piston rod.

Remove damaged capillaries (instruments from 20 µl)


CAUTION


Risk of injury

Pull only front fragment over the seal. Wear protective gloves.


- a. Loosen the union nut (half turn).
- b. Pull off the front fragment (1) from the seal.

- c. Pull rear fragments from holder, wrap in cloth, and crush with pliers.
- d. Insert new capillary.

Insert new capillary

- a. Push the capillary up to the stop (1).
- b. Tighten the union nut by hand.
- c. Check correct seating/attachment of the capillary through the inspection window.

NOTICE

The color marking (color code) of the pipette button, the capillaries, and the accessories must match.

6.2 Change piston rod

Remove capillary

CAUTION

Risk of contamination

Do not touch the contaminated area of the capillaries. Wear protective gloves.

- Loosen the union nut.
- Pull the capillary off the piston rod.
- Unscrew the union nut.

Remove piston rod

(Tools in repair set, see Ordering data for instruments and spare parts, p. 50)

a.

b.

- Unscrew the slotted screw.
- Pull off the housing.
- Press the pipetting button.


d.

- d. Allen screws appear in the recess (1). With digital instruments, the pipetting button may have to be turned slightly for this purpose.
- e. Loosen Allen screws (half turn).
- f. Pull out the piston rod.

Install piston rod

b.

- a. Place the setting gauge from the repair set vertically on a hard, level surface. Frosted end upwards.
- b. Insert the piston rod into the shaft, and place the shaft on the setting gauge.

c.

- c. Press the pipetting button as far as it will go. The piston rod and shaft must be in a vertical position (do not tilt).
- d. Tighten the Allen screws (with the pipetting button pressed).


e.

e. Place the housing, and screw in the slotted screw.

f. Screw the union nut onto the shaft (do not tighten).

h.


g. Push the capillary up to the stop (1) (color code upwards).

h. Tighten the union nut by hand.


i. Check correct seating/attachment of the capillary through the inspection window.

6.3 Change seal

Remove old seal

- a. Remove capillary (see Change capillary, p. 40).
- b. Cut the old seal with a sharp knife (scalpel or similar), and remove it.

- c. Insert a new seal (1) into the attachment block (2) (included in the repair set, see page 22).
- d. Press the piston rod firmly into the seal.
- e. Pull off the attachment block.
- f. Check that the seal is firmly seated.
- g. Slide on the capillary (see Change capillary, p. 40).

Check instrument

- a. Check the volume after each change of capillary, piston rod, or seal.

See also

→ Ordering data for instruments and spare parts (p. 50)

7 Checking the volume

We recommend testing the instrument every 3 to 12 months depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. For gravimetric volume testing, we recommend the test procedure as described (e.g., in DIN EN ISO 8655 Part 6).

1. Setting the nominal volume

- a. Set the maximum specified instrument volume (for procedure, see).

2. Conditioning the pipette

- a. Condition pipette before testing by pipetting once with test liquid (distilled water).

3. Performing the test

- a. Aspirate the test liquid and pipette into the weighing vessel.
- b. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- c. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- d. At least 10 pipettings and weighings are recommended. With a digital pipette, perform at least 10 weighings in three volume ranges (100%, 50%, 20%).

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPA)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and coefficient of variation (CV%): A% and CV% are calculated according to the formulas of statistical quality control.

NOTICE


Test Instructions (SOPs) are available for download from www.brand.de.

8 Troubleshooting

Problem	Possible cause	Corrective action
Pipetted volume too large/ air bubbles in liquid taken up.	Capillary is not seated correctly.	Loosen the union nut, push the capillary as far as it will go, and tighten the union nut. See Change capillary, p. 40
Pipette is dripping.	Capillary damaged.	Change capillary (see Change capillary, p. 40)
	Piston rod damaged (instruments up to 10 μ l).	Repair instrument (see Repairs, p. 51) Change piston rod (see Change piston rod, p. 42) Clean instrument (see Cleaning, p. 39)
	Seal damaged (instru- ments from 20 μ l).	Change seal (see Change seal, p. 45) Repair instrument (see Repairs, p. 51)
	Instrument soiled.	Clean instrument (see Cleaning, p. 39)

9 Accuracy table

Volume	Accuracy $R \leq \pm$		Coefficient of variation $CV \leq$	
	%	μl	%	μl
Variable				
2.5–10	3.0	0.3	0.8	0.08
5–25	2.4	0.6	0.5	0.125
10–50	1.8	0.9	0.4	0.2
20–100	1.5	1.5	0.4	0.4
Fixed				
1	12.0	0.12	4.0	0.04
2	7.5	0.15	2.0	0.04
5	3.0	0.15	0.8	0.04
10	3.0	0.3	0.8	0.08
20	2.4	0.48	0.5	0.1
25	2.4	0.6	0.4	0.1
50	1.8	0.9	0.4	0.2
100	1.5	1.5	0.4	0.4
200	1.5	3	0.2	0.4

Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

10 Ordering data for instruments and spare parts

Volume	Color code	Instru- ment	Capillaries ¹	Seals ²	Piston rod ³	Repair set ⁴
variable						
2.5–10	Orange	701807	701902	—	701930	701965
5–25	2 × white	701812	701906	701920	701932	701966
10 - 50	Green	701817	701908	701922	701934	701967
20 - 100	Blue	701822	701910	701924	701936	701968
Fixed						
1	White	701842	701900	—	701928	701964
2	White	701844	701900	—	701928	701964
5	White	701853	701900	—	701928	701964
10	Orange	701858	701902	—	701930	701965
20	Black	701863	701904	701920	701932	701966
25	2 × white	701864	701906	701920	701932	701966
50	Green	701868	701908	701922	701934	701967
100	Blue	701873	701910	701924	701936	701968
200	Red	701878	701910	701924	701938	701968

1) Packaging unit 100 pieces; for 100/200 µl: 50 pieces.

2) Packaging unit 3 pieces.

3) Packaging unit 3 pieces; from 20 µl, already equipped with seal.

4) Repair set consisting of: 1 Allen key, 1 screwdriver, 1 setting gauge, 1 attachment block (from 20 µl), 1 piston rod (from 20 µl with seal attached), 3 seals, 3 sealing rings, 1 slotted screw.

10.1 Accessories

Description	Pkg unit	Order No.
Transferpettor station For storage of four instruments up to 200 µl with accessories.	1 pc.	701960

11 Product markings

Symbol or number	Meaning
	Read the user manual.
	Read the user manual.
XXZXXXXX	Serial number

12 Repairs

12.1 Sending for repair

NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Only send the device without a battery installed.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

13 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

14 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

15 Disposal

Before disposal, observe the relevant national disposal regulations, and ensure that the product is disposed of properly.

Subject to technical changes, errors, and misprints.

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