Liquid Handling · Easy Handling!

# Transferpette®

## 2 ml, 0,5 - 5 ml

FIRST CLASS · BRAND

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### Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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## Contents

Seite		Page
4	Safety Instructions	4
6	Purpose	6
6	Operating Limitationes	6
7	Operating Elements	7
8	Handling	8
9	Pipetting	9
11	Cleaning	11
12	Maintenance	12
13	UV sterilization	13
13	Autoclaving	13
15	Adjustment	15
18	Checking the Volume	18
	Calculation of Accuracy/ Coefficient of Variation	
19		19
20	Troubleshooting	20
22	Specifications/Spares	22
23	Accessories	23
24	Return for repair	24
26	Contact addresses	26
27	Calibration Service	27
28	Warranty Information	28
30	Disposal	30

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## Safety Instructions

### Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe all specifications provided by reagent manufacturers.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations.

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Observe operating exclusions (see page 7).  
If in doubt, contact the manufacturer or supplier.

5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visible damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 21), and contact the manufacturer if needed.

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## Purpose

Air-displacement pipette for pipetting aqueous solutions of medium density and viscosity.

## Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C to 40 °C (59 °F to 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

## Operating Limitations

Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than  $\pm 5$  °C/ $\pm 41$  °F.

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## Operating Exclusions

The user has to ensure in each case the compatibility of the instrument with the intended application. This instrument cannot be used:


- for liquids incompatible with polypropylene or FKM
- for liquids of a very high vapor pressure
- in Reverse Mode (ISO 8655/2)

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

## Handling

### Use proper pipette tips

Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend to use BRAND pipette tips.

Pipette tips are disposables.

### **(A) Volume setting**

### **(B) Hand posture while pipetting**

Put your thumb across the pipetting key – unlike with conventional pipettes.

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

## Pipetting

### Aspirating the sample

Rinse the tip once with the sample liquid.

1. Press the pipetting key to the first stop (A).
2. Immerse the pipette tip 2 - 3 millimeters into the sample.
3. Let the pipetting key slide back slowly.

#### Note:

Leave the tip immersed in the liquid for about another second to avoid the intake of air.

4. Lightly wipe the tip against the wall of the vessel.

#### Note:

Don't lay the instrument horizontal when the tip is full, or liquid may penetrate inside the instrument.

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

## Dispensing the sample

Place the pipette tip against the wall of the vessel.

1. Press the pipetting key down to the stop (A) and hold it down.

### Note:

For serums and liquids of high viscosity or low surface tension, observe adequate waiting time to improve the accuracy.

2. The blow-out stroke empties the tip completely: Press down to the stop (B).

### Note:

When working with serum and liquids of high viscosity, low surface tension or at low volumes: Rinse the tip with the reagent contained in the vessel to improve the accuracy.



3. Wipe the pipette tip against the wall of the vessel.  
Let the pipetting key slide back.

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## Cleaning

### Note:

The Transferpette® requires regular maintenance, depending on frequency of use and pipetted media. After cleaning, always check volume (see page 18).

- Eject the tip
- Press the two lateral closures (A) and pull off the lower part of the tip ejector (B).
- Unscrew the piston assembly with lower shaft end (C) from the handle.
- Unscrew lower shaft end (E) from the piston assembly (D).
- Pull off the O-ring on the piston assembly toward the front, and clean it.
- Clean the inside of the piston assembly (D) and the lower shaft end (E) and clean them with isopropanol.
- Lubricate the O-ring inside and outside, put it back in place ('Maintenance', page 12).
- Reassemble parts.

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
## Maintenance

### At regular intervals:

- Check the piston assembly and the seal for contamination. Clean them if necessary.
- Check volume (see page 18).
- Check the pipetting system - consisting of Transferpette® and pipette tip - for tightness.

We recommend using the BRAND leak testing instrument PLT unit.

### Alternatively:

After drawing in the sample:  
Hold down the instrument vertically for about 10 seconds. If a drop forms at the tip, "Trouble Shooting", page 21.

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### Lubricating the seal (O-ring)

Lubricating is required:

- after cleaning,
- if motion is uneven
- after autoclaving.

#### Note:

Use only silicone grease (see page 23).

## UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

## Autoclaving

For sterile work or for decontamination, the pipette parts shown here are autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285.

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
- continued -

### Attention:

It is the user's responsibility to ensure effective autoclaving.

**Note:** After autoclaving, always check volume (see page 18). Only autoclave parts that have been cleaned (see 'Cleaning', page 11).

- Eject the tip
- Press the two lateral closures and pull off the lower part of the tip ejector.
- Unscrew the entire lower part of the pipette from the handle.
- Unscrew the pipetting assembly.
- Autoclave the piston assembly and the lower shaft end.
- Allow the parts to cool completely and assemble in the reverse order.
- Check volume (see page 18).

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

1/250 vom Nennvolumen

### Note:

The handle is not autoclavable. For disinfection you may use commercially available disinfecting solutions as long as they do not attack PP.

## Adjustment

The instrument is permanently adjusted for aqueous solutions. It can be adjusted for solutions of other density and viscosity.

### Fixed-volume models

- Check the volume, determine actual value (see page 18).
- Remove ejector cap (see page 16).
- Adjust by means of the calibration key.

### Note:

Maximum permissible adjustment is 3 full turns.

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


## Adjustment

– continued –

- Use a paper clip to extract the red distance plate (B) from the top.
- Pull out volume setting knob mechanism (C) to the stop in the direction of its axis (approx. 2 mm).

- Set to the previously measured actual value. Push back volume setting mechanism to the stop (approx. 2 mm), without turning it.
- Replace the red security plate.
- Mount ejector cap (A).
- Secure ejector cap by turning the lateral closures into position .

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## Checking the volume

The gravimetric volume test of the pipette requires the following steps:

(Testing procedures are described e.g., in ISO 8655 Part 6, Testing Instruction see page 25).

Adjust volume to nominal capacity.

### **Conditioning the pipette:**

- To condition the pipette before testing, mount a new pipette tip and pipette five times with testing liquid ( $\text{H}_2\text{O}$  dist.). Discard the pipette tip. Carrying out the test (observe instructions by the balance manufacturer):
- Mount a new pipette tip and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature into account.
- Carry out at least 10 pipetting and weighing operations.

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## Calculation of accuracy and coefficient of variation

Accuracy (A%) and coefficient of variation (CV%) are calculated according to the equations used in statistical quality control.

**Calculations** (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad \begin{array}{l} x_i = \text{results of weighings} \\ n = \text{number of weighings} \end{array}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

$$Z = \text{correction factor} \\ (\text{e.g., } 1,0029 \text{ } \mu\text{l/mg at } 20 \text{ } ^\circ\text{C, } 1013 \text{ hPa})$$

### Accuracy

$$\mathbf{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100 \quad V_0 = \text{nominal volume}$$

### Standard deviation Coefficient of variation

$$\mathbf{s} = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}} \quad \mathbf{CV\%} = \frac{100 \text{ s}}{\bar{V}}$$

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## Troubleshooting

Trouble	Possible Causes	Action to be taken
Tip dripping (instrument leaks) or volume too low	<ul style="list-style-type: none"><li>- Unsuitable tip</li><li>- Tip not seated tightly</li><li>- Seal contaminated</li><li>- Seal damaged</li> <li>- Non-aqueous solutions pipetted</li></ul>	<ul style="list-style-type: none"><li>- Only use high-quality tips</li><li>- Press tip on firmly</li><li>⇒ 'Cleaning', page 11</li><li>⇒ 'Specifications, Spares', page 22</li> <li>⇒ 'Adjustment', page 15</li></ul>
Volume too high	<ul style="list-style-type: none"><li>- Pipetting key pressed too far (into blow out position) before sample take up</li></ul>	<ul style="list-style-type: none"><li>- Operate properly</li><li>⇒ 'Pipetting', page 9</li></ul>

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
**20 °C  
EX**

## Specifications/Spares

(Spare parts can be fitted by the user)

Volumen (ml)	Best.-Nr.	Richtigkeit* (R)		Variations- koeffizient* (VK)		Abwerfer- haube (A) Best.-Nr.	Kolbeneinheit, kpl. mit O-Ring (B) Best.-Nr.	O-Ring für Kolbeneinheit (C) Best.-Nr.	Schaftunterteil (D) Best.-Nr.	Abwerferunter- teil, kpl. (E) Best.-Nr.
Capacity (ml)	Cat.No.	Accuracy* (A) ± %    ± µl		Coefficient of variation (CV) %        µl		Ejector cap (A) Cat. No.	Piston assem- bly, complete, with O-ring (B) Cat. No.	O-ring for piston assem- bly (C) Cat. No.	Lower shaft end (D) Cat. No.	Tip ejector, sleeve (E) Cat. No.
2	7041 64	0,6	12	0,2	4	9398	7289	7287	7032 46	7299
0,5 - 5	7041 82	0,6	30	0,2	10	9397	7290	7288	7032 47	7299

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\* Final test values relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist. H<sub>2</sub>O, and with smooth, jerk-free operation.

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## Accessories

<b>Rack</b> for 3 instruments	Cat. No. 7032 03
<b>Rack</b> for 2 x 3 instruments	Cat. No. 7032 08
<b>Silicon grease</b>	Cat. No. 7036 77
<b>Filter</b> , pack of 25 pc.	Cat. No. 7046 52
<b>PLT unit</b>	Cat. No. 7039 70

### **BRAND pipette tips**, 5 ml, PP

1 bag of 200 tips, bulk	Cat. No. 7025 95
5 bags of 200 tips, bulk	Cat. No. 7026 00
1 box, PP, with 28 tips	Cat. No. 7026 05

Following tolerances are typically (manufacturer's experience) obtained by new instruments under optimized testing conditions (qualified operators and standardized ambience conditions):

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## Return for Repair

### **Attention:**

Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

### **Outside the U.S. and Canada:**

- Complete the "Declaration on Absence of Health Hazards" and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

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### **In the U.S. and Canada:**

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

## **Kontaktadressen/Contact addresses**

### **BRAND GMBH + CO KG**

Otto-Schott-Straße 25  
97877 Wertheim (Germany)

Tel.: +49 9342 808-0

Fax: +49 9342 808-98000

E-Mail: [info@brand.de](mailto:info@brand.de)

[www.brand.de](http://www.brand.de)

### **USA and Canada:**

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)

Tel.: +1-860-767 2562

Fax: +1-860-767 2563

[www.brandtech.com](http://www.brandtech.com)

### **India:**

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park, Powai  
Mumbai - 400 076 (India)

Tel.: +91 22 42957790

Fax: +91 22 42957791

E-Mail: [info@brand.co.in](mailto:info@brand.co.in)

[www.brand.co.in](http://www.brand.co.in)

### **China:**

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B  
No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)

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## Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend to check the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers you the possibility, to have your instruments calibrated by the BRAND calibrating service or the BRAND-owned DAkkS calibration service. Just send in the instruments to be calibrated accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND calibrating service) or with a DAkkS calibration certificate. For further information please contact your dealer or BRAND.

Complete ordering information is available for download at [www.brand.de](http://www.brand.de) (see technical documentation).

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## **Warranty**

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual.

We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:  
Information for warranty please see [www.brandtech.com](http://www.brandtech.com).

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## **Disposal**

For the disposal of instruments and tips, please observe the relevant national disposal regulations valid in your area.

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Subject to technical modification without notice.  
Errors excepted.