Liquid Handling · Easy Handling!

# Transferpette -8/-12

FIRST CLASS · BRAND

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## Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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## **Contents**

	Page
Safety Instructions	4
Purpose, Limitations of use	6
Operating Limitations, Operating exclusions	7
Operating Elements	8
Handling	8
Pipetting	10
Maintenance	13
Adjustment	14
Checking the volume	15
Components	18
Cleaning	19
Repairs	22
UV sterilization	29
Autoclaving	30
Trouble Shooting	32
Specifications	33
Ordering Data	34
Spare Parts and Accessories	34
Return for Repair	37
Contact addresses	39
Calibration Service	40
Warranty	41
Disposal	42

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## Safety Instructions

### Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe all specifications provided by reagent manufacturers.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations.

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Observe operating exclusions (see page 7).  
If in doubt, contact the manufacturer or supplier.

5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visual damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the ‘Troubleshooting’ section of this manual (see page 32), and contact the manufacturer if needed.

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## Purpose

The multichannel piston-operated pipette is an air-displacement system for pipetting aqueous solutions with medium density and viscosity into microtiter plates.

## Limitations of use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C to 40 °C (59 °F to 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

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## Operating limitations

Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than  $\pm 5\text{ °C}/\pm 41\text{ °F}$ .

## Operating exclusions

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:


- for liquids incompatible with polypropylene or FKM
- for liquids of a very high vapour pressure
- in Reverse Mode (ISO 8655/2)

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## Handling

### Use proper pipette tips

Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend the use of BRAND pipette tips. Pipette tips are disposables.

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## Adjusting the manifold


The **manifold** can be rotated counterclockwise to the most convenient pipetting position.

Counterclockwise rotation  
for ease of handling

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### Note:

The manifold must be completely mounted to the handle for accurate pipetting. During the counterclockwise rotation, a clicking noise must be heard. Turning in the opposite direction will loosen the manifold.


## Hand posture while pipetting

Securely attach tips to nose cone assembly. Put your thumb **across** the pipetting key – unlike with conventional pipettes.

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
## Pipetting

### Aspirating the sample

- Rinse the tips once with the sample liquid (at normal pipetting speed).
1. Press the lateral pipetting key to **the first stop (A)**.
  2. Hold the pipette vertically and immerse the tip into the liquid.

#### Note:

Leave the tips immersed in the liquid for a few seconds to avoid the intake of air.

3. Let the pipetting key slide back **slowly**.
4. Lightly wipe the tips against the wall of the vessel.


#### Note:

Don't lay the instrument horizontally when the tips are filled, or liquid may flow inside the instrument.

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
## Discharging the sample

Lightly place the pipette tips against the wall of the vessel.

- Press the pipetting key down to the **stop (A)** and hold it down.

### Note:

For serums and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- The blow-out stroke empties the tips completely: Press down to the **stop (B)**.

### Note:

To improve accuracy when working at low volumes: rinse the tips with the reagent contained in the vessel

- Lightly wipe the pipette tips against the wall of the vessel.
- Let the pipetting key slide back.

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

### Ejecting the tips (A)

Press the ejector cap down firmly in the region of the color code (see arrow).

Tips that are jammed on too tightly may need to be pulled off separately by hand.

### Volume setting (B)

Set the desired volume by turning the volume setting knob (► "Volume setting knob", page 8).

### After use:

Always store the instrument upright on the **stand** supplied.

#### Note:

To reduce strain on the sealing rings (► "V-rings", page 26), always store the instrument without the tips attached.

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
## Maintenance

### Check for tightness

At regular intervals, check the pipetting system – consisting of Transferpette®-8/-12 and pipette tip – for tightness.

We recommend using the BRAND leak testing instrument PLT unit.

Alternatively: after drawing in the sample: Hold down the instrument vertically for about 10 seconds. If a drop forms at one of the tips, ➤ “Trouble Shooting”, page 32.

### Lubricating the plunger

After autoclaving the manifold several times, it is necessary to lubricate the pistons to maintain smooth operation:

Unscrew nose cones, push pistons downward, and lubricate (➤ “Repairs”, page 22, 23).



To lubricate the pistons, use only the silicone grease supplied with the instrument.

The silicone grease is also available separately as accessory (➤ “Ordering Data”, page 34).

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
## Adjustment

The instrument is permanently adjusted for aqueous solutions. It can be adjusted for solutions of other density and viscosity.

1. Check the volume, determine actual value (► page 17).
2. Release lateral closures by turning into position ①.
3. Remove ejector cap (A).
4. Use a paper clip to extract the red distance plate (B) from the top.
5. Pull out volume setting knob (C) to the stop in the direction of its axis (approx. 2 mm).

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6. Set to the previously measured actual value. Push back volume setting knob to the stop (approx. 2 mm), without turning it.
7. Replace the red distance plate across the axis.
8. Mount ejector cap (A).
9. Secure ejector cap by turning the lateral closures into position .

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## Checking the volume

The instrument should be calibrated at least once a year. A more frequent interval may be dictated by individual lab rules, applications and conditions.

The gravimetric volume test of the pipette requires the following steps (page 16):

(Testing procedures are described e.g., in ISO 8655, Part 6, ► “Testing Instruction”, page 37).

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Adjust volume to nominal capacity.

### Conditioning the pipette:

To condition the pipette before testing, mount new pipette tips and pipette five times with testing liquid ( $\text{H}_2\text{O}$  dist.). Discard the pipette tips.

### Carrying out the test (observe instructions by the balance manufacturer):

- Mount a new pipette tip and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature of the test liquid into account.
- The number of tests per channel should be selected according to the required accuracy. We recommend 5-10 tests per channel in 3 volume ranges (100 %, 50 %, 10 %).

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## Calculation of accuracy and coefficient of variation

Accuracy (A%) and coefficient of variation (CV%) are calculated according to the equations used in statistical quality control.

**Calculations** (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad \begin{array}{l} x_i = \text{results of weighings} \\ n = \text{number of weighings} \end{array}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

$$Z = \text{correction factor} \\ (\text{e.g., } 1,0029 \text{ } \mu\text{l/mg at } 20 \text{ } ^\circ\text{C, } 1013 \text{ hPa})$$

### Accuracy


$$\mathbf{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100 \quad V_0 = \text{nominal volume}$$

### Standard deviation Coefficient of variation

$$\mathbf{s} = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}} \quad \mathbf{CV\%} = \frac{100 \text{ s}}{\bar{V}}$$

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## Components

- ① Handle
- ② Manifold
- ③ Manifold housing
- ④ Nose cone assembly
- ⑤ Individual nose cone
- ⑥ Piston spring
- ⑦ Piston
- ⑧ Piston support bar
- ⑨ Manifold housing cover
- ⑩ Seal
- ⑪ Mounting tool
- ⑫ Central guide rod

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
## Cleaning

**Note:** The Transferpette®-8/-12 requires regular maintenance, depending on frequency of use and pipetted media.

## Dismantling the instrument

**Be careful:** Work over a tray to collect any part that may fall out during dismantling.

- Eject the tips.
- Set the instrument to nominal volume.
- Carefully remove the manifold (②) from the handle (①).
- Turn both closures of the manifold housing cover (⑨) with a coin (rotation through 90°) and slide off the housing.
- **Remove both outer screws (Phillips). Do not loosen the central guide rod (⑫)!**
- Remove the nose cone assembly (④), piston springs (⑥), and piston support bar (⑧).
- Unscrew single nose cones with seals from the nose cone support bar (► "Repairs", page 22).

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## Cleaning the components

- Clean single nose cones (⑤), pistons (⑦), piston support bar (⑧) and piston springs (⑥) (**these components only**) with an appropriate cleaning agent (e.g., water, alcohol. Do not use organic solvents!).
- Let these parts dry completely. Residual moisture in the nose cone assembly may result in a loss of accuracy.
- Lubricate the pistons lightly (only with the supplied silicone grease).

### Important:

To lubricate the central guide rod, use only the recommended fluorstatic grease (► page 34)!

### Note:

To protect against contamination, filter tips can also be used.

## Assembling the instrument

- Reassemble the nose cone assembly without cones (④), piston springs (⑥), piston support bar (⑧) and manifold housing cover (⑨). Tighten all screws.

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- Set volume to minimum.
- Push the seals on the pistons with the flat side facing upward (direction of the handle) (► Fig. 4, page 23); use the mounting tool to mount and tighten the single nose cones clockwise.
- Slide on the manifold housing (③). Turn the closures into place by rotation through 90°.
- Push the handle and the manifold together and turn until they click into place.

#### **Important:**

Do not press on the ejector cap during this procedure.

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#### **After cleaning:**

- Check if pipetting and ejector functions are correct. If this is not the case, dismantle the instrument and carefully reassemble it. If necessary, lubricate the pistons (► "Maintenance", page 13).
- Check for proper fit of tips.
- Check for correct volume delivery (► page 15).

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

## Repairs

► “Accessories and spare parts”,  
pages 34-36.

### Replacing of single nose cones and seals



#### Removing a single nose cone

- Set the instrument to nominal volume.
- Turn both closures of the manifold housing cover (⑨) with a coin (rotation through 90°) and slide off the housing.
- Push the mounting tool (⑪) on a single nose cone (Fig. 1).
- Turn counterclockwise to loosen and remove the nose cone (Fig. 2).

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### Note:

The seal will either remain inside the nose cone, or will stay on the piston after the nose cone is removed. By pressing the pipetting key it can easily be removed and cleaned (Fig. 3).

- Set volume to minimum.
- Push down the pipetting key completely.
- If required, lubricate the piston lightly with the **supplied** silicone grease.
- Push the seal on the piston with its flat side facing upward (direction of the handle) (Fig. 4).
- Use the mounting tool to mount and tighten the new nose cone clockwise.
- Reassemble the manifold (► page 20).
- Check instrument for tightness (► “Check for tightness”, page 13).

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

### Replacing the pistons

- Unscrew the manifold from the handle.
- Turn both closures of the manifold housing cover (Ⓢ) with a coin (rotation through 90°) and slide off the housing.
- **Remove both outer screws (Phillips)** which hold the nose cone assembly (④) to the metal support rods of the manifold housing cover (Ⓢ). **Do not loosen the central guide rod (► page 18)!**

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- Push the mounting tool (⑩) on a piston (⑦).
- Turn counterclockwise to loosen and remove the piston (⑦) (along with the spring ⑬).
- Insert a new piston.
- Use the mounting tool to mount and tighten the new piston.
- Reassemble the instrument (▶ page 20) and check for tightness (▶ page 13).

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## **Replacement of nose cone V-rings (O-rings)**




### **Required materials:**

- 1 V-Ring removal tool (square)
- 1 V-Ring mounting tool (cone shaped)\*
- Replacement V-rings (O-rings)

\* The mounting tool is only necessary for replacing the V-rings on instruments of the volume ranges 2.5-25 µl, 5-50 µl, 10-100 µl, 20-200 µl and 30-300 µl.

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
### To remove an old V-ring (O-ring):

1. Position the square removal tool ① under one of the nose cones.
2. Move the tool up against the nose as far as possible to push the V-ring (O-ring) up, away from the nose cone.
3. If necessary, cut the dislodged portion of the V-ring (O-ring) with a scissor and discard. Do not damage the nose cone.

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
### Mounting a new V-ring (O-ring):

4. Only for instruments of size 2.5-25  $\mu\text{l}$ , 5-50  $\mu\text{l}$ , 10-100  $\mu\text{l}$ , 20-200  $\mu\text{l}$  and 30-300  $\mu\text{l}$ : Slide the cone-shaped mounting tool ① over the nose cone.


For instruments of size 0.5-10  $\mu\text{l}$  and 2-20  $\mu\text{l}$ , no mounting tool is required; the O-ring can be placed directly onto the nose cone.

### Note:

Moisten mounting tool, nose cone and V-ring (O-ring) with dist.  $\text{H}_2\text{O}$  in order to facilitate mounting.

5. Place the new V-ring (O-ring) onto the mounting tool. The widest side of the V-ring should face towards nose cone manifold. Push the V-ring (O-ring) over the mounting tool into the groove on the nose cone.

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6. Remove the mounting tool.  
Check correct seat of the V-ring (O-ring). If twisted, turn it several times around the nose cone to restore the V-ring (O-ring) to its correct shape.

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### **After any repair:**

- Check if pipetting and ejector functions are correct.
- Check for proper fit of the tips.
- Check for correct volume delivery (► page 15).

### **After every change of the complete manifold:**

- Check for correct volume delivery (► page 15).  
Read just instrument if necessary.

## **UV sterilization**

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

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## Autoclaving

### Caution:

Only the manifold may be autoclaved. Eject the tips before autoclaving! For disinfection the handle you may use commercially available disinfecting solutions as long as they do not attack PP.

- Remove the manifold from the handle.
- For sterile work or for decontamination, the complete manifold is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285. It is the user's responsibility to ensure effective sterilization and correct autoclave temperature.

### After autoclaving

- Allow the manifold to cool down completely.
- Push the manifold and the handle together and turn until they click into place. **Important:** Do not press the ejector cap during this procedure.
- If necessary, tighten the nose cones with a sterile pipette tip of proper size.
- Check for smooth function of pipetting and ejector keys.
- Check for proper fit of pipette tips.
- Check for correct volume delivery (► page 15).

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## Trouble Shooting

Problem	Possible Cause	Corrective Action
Pipette tip leaks or volume too low	<ul style="list-style-type: none"> <li>- Tip does not sit properly</li> <li>- Damaged V-ring (O-ring) on single nose cone</li> <li>- Piston damaged</li> <li>- Seal dirty, damaged</li> <li>- Instrument is dirty</li> <li>- Non-aqueous solutions pipetted</li> </ul>	<ul style="list-style-type: none"> <li>- Press tip on firmly. Use only suitable high-quality tips, page 8.</li> <li>- Replace V-ring (O-ring) on nose cones. "Repairs", pages 26-28.</li> <li>- Replace piston. "Repairs", pages 24, 25.</li> <li>- "Repairs", pages 22, 23.</li> <li>- Clean instrument. "Cleaning", page 19.</li> <li>- "Adjustment", page 14.</li> </ul>
Poor sample aspiration	<ul style="list-style-type: none"> <li>- Single nose cone clogged</li> </ul>	<ul style="list-style-type: none"> <li>- Clean nose cone. Remove cones only; do not disassemble instrument! "Cleaning" page 19.</li> </ul>
Volume too large	<ul style="list-style-type: none"> <li>- Pipetting key pressed too far down into the blow-out stroke before aspirating liquid.</li> </ul>	<ul style="list-style-type: none"> <li>- Operate pipette properly. "Operating the pipette", page 10.</li> </ul>
Pipetting key moving stiffly	<ul style="list-style-type: none"> <li>- Key movement reduced</li> </ul>	<ul style="list-style-type: none"> <li>- Lubricate piston and centrale guide rod with the appropriate grease, pages 13, 34.</li> </ul>

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The error limits of ISO 8655 are not exceeded.

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\* Test values determined at the same temperature of instrument, liquid (H<sub>2</sub>O dist.) and ambience (20 °C) and smooth, jerk-free operation.

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When pipetting liquids with different viscosity and/or flow properties than H<sub>2</sub>O dist., deviations from the values stated can occur.

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Accuracy and coefficient of variation refer to each single channel.

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**BRAND TipBox**, empty Cat. No. 7329 92

**Reservoirs** for reagent aspiration  
(set of 10), autoclavable at 121 °C Cat. No. 7034 59

**Stand** for one instrument Cat. No. 7034 40

**Silicone grease** for piston  
(tube, approx. 1 g) Cat. No. 7036 77

**Fluorstatic grease** for central  
guide rod (tube, approx. 1 g) Cat. No. 7036 78

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**Individual nose cone** (2 pieces, incl. seals and mounting tool)

For capacity, $\mu\text{l}$	Cat. No.
0.5 - 10	7056 77
2 - 20	7056 78
2.5 - 25	7033 33
5 - 50	7056 81
10 - 100	7056 82
20 - 200	7056 83
30 - 300	7056 84

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Für Volumen, $\mu\text{l}$	Best.-Nr.
0,5 - 10	7033 40
2 - 20	7033 41
2,5 - 25	7033 42
5 - 50	7033 43
10 - 100	7033 44
20 - 200	7033 45
30 - 300	7033 46

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#### PLT unit

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Best.-Nr. 7039 70

**V-Rings** (Set of 24), incl. removing and mounting device  
(except 10 µl and 20 µl capacity: no mounting tool needed)

For capacity, µl	Cat. No.
10, 20 (O-ring)	7033 80
25, 50, 100, 200, 300	7034 91

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## Return for Repair

### **Attention:**

Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

### **Outside the U.S. and Canada:**

- Complete the "Declaration on Absence of Health Hazards" and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

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### **In the U.S. and Canada:**

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

## Kontaktadressen/Contact addresses

### **BRAND GMBH + CO KG**

Otto-Schott-Straße 25  
97877 Wertheim (Germany)

Tel.: +49 9342 808-0

Fax: +49 9342 808-98000

E-Mail: [info@brand.de](mailto:info@brand.de)

[www.brand.de](http://www.brand.de)

### **India:**

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park, Powai  
Mumbai - 400 076 (India)

Tel.: +91 22 42957790

Fax: +91 22 42957791

E-Mail: [info@brand.co.in](mailto:info@brand.co.in)

[www.brand.co.in](http://www.brand.co.in)

### **USA and Canada:**

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)

Tel.: +1-860-767 2562

Fax: +1-860-767 2563

[www.brandtech.com](http://www.brandtech.com)

### **China:**

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B  
No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)

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## Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend to check the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers you the possibility, to have your instruments calibrated by the BRAND calibrating service or the BRAND-owned DAkkS calibration service. Just send in the instruments to be calibrated accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND calibrating service) or with a DAkkS calibration certificate. For further information please contact your dealer or BRAND.

Complete ordering information is available for download at [www.brand.de](http://www.brand.de) (see technical documentation).

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## Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual.

We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

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## **Disposal**

For the disposal of instruments and tips, please observe the relevant national disposal regulations valid in your area.

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Subject to technical modification without notice.  
Errors excepted.