Liquid Handling · Easy Handling!

Transferpette®

FIRST CLASS · BRAND

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Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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Contents

Page

4	Safety Instructions	4
6	Purpose, Limitations of use	6
	Operating limitations, Operating exclusions	7
7	Operating Elements	8
8	Handling	8
8	- Use proper pipette tips	8
8	- Volume Setting	9
9	- Hand posture while pipetting	9
9	- Ejecting the tip	10
10	Pipetting	11
11	- Aspirating the sample	11
11	- Dispensing the sample	12
12	Maintenance	13
13	- Lubricating the plunger	13
13	Adjustment	14
14	- Fixed-volume models	14
14	- Variable-volume models	15
15	Checking the volume	17
17	Calculation of accuracy/coefficient of variation	18
18	Cleaning	19
19	UV sterilization / Autoclaving	21
21	Troubleshooting	25
24	Specifications / Spares	26
26	Accessories	29
29	Return for repair	30
30	Contact addresses	32
32	Calibration Service	33
33	Warranty	34
34	Disposal	35
35		

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Safety Instructions

Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe all specifications provided by reagent manufacturers.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations.

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Observe operating exclusions (see page 7).

If in doubt, contact the manufacturer or supplier.

5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visual damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the ‘Troubleshooting’ section of this manual (see page 25), and contact the manufacturer if needed.

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Purpose

Air-displacement pipette for pipetting aqueous solutions of medium density and viscosity.

Limitations of use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C and 40 °C (59 °F and 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

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Operating limitations

Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than $\pm 5\text{ }^{\circ}\text{C}/\pm 41\text{ }^{\circ}\text{F}$.

Operating exclusions

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene
- for liquids of a very high vapour pressure
- in Reverse Mode (ISO 8655/2)

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
Handling

Use proper pipette tips


Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend the use of BRAND pipette tips. Pipette tips are disposables

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- continued -

Fitting the tip

- in regular intervals, check the pipetting system-consisting of Transferpette® and pipette tip for tightness. We recommend using the BRAND leak testing instrument PLT unit.
- turn them lightly in the direction of the arrow

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
Note:

If the tip is turned in the other direction, the screw connection of the shaft tip may be loosened, with air leaks resulting. The shaft tip must be tightened at the latest when a red ring becomes visible on it.

Ⓐ Volume setting

Ⓑ Hand posture while pipetting


Put your thumb across the pipetting key - unlike with conventional pipettes.

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

Pipetting

Aspirating the sample

- Rinse the tip once with the sample liquid.

1. Press the lateral pipetting key to **the first stop (A)**.
2. Immerse the pipette tip 2-3 millimeters into the sample.
3. Let the pipetting key slide back **slowly**.

Note:

Leave the tip immersed in the liquid for about another second to avoid the intake of air.

4. Lightly wipe the tip against the wall of the vessel.

Note:

Don't lay the instrument horizontal when the tip is filled, or liquid may penetrate inside the instrument.

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

Dispensing the sample

Place the pipette tip against the wall of the vessel.

- Press the pipetting key down to the stop (A) and hold it down.

Note:

For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- The blow-out stroke empties the tip completely: Press the pipetting key down to the stop (B).

Note:

When working with serum and liquids of high viscosity, low surface tension or at low volumes: Rinse the tip with the reagent contained in the vessel to improve the accuracy.

- Wipe the pipette tip against the wall of the vessel.
- Let the pipetting key slide back.

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
Maintenance

In regular intervals:

1. Check firm seat of pipette shaft cone, tighten if necessary.
2. Examine piston and seal for contamination (see page 19 “Cleaning”).
3. Check volume.

Lubricating the plunger

Very thinly lubricate the plunger:

- after cleaning,
- if motion is uneven,
- after autoclaving.

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
Adjustment

The instrument is permanently adjusted for aqueous solutions. It can be adjusted for solutions of other density and viscosity.

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
1/250 do volume nominal

Note:


Maximum permissible adjustment is 3 full turns.

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
Variable-volume models


1. Check the volume, determine actual value (see page 17, "Checking the volume").
2. Release lateral closures by turning into position .
3. Pull off ejector cap (A).

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4. Use a paper clip to extract the red distance plate (B) from the top.
5. Pull out volume setting knob (C) to the stop in the direction of its axis (approx. 2 mm).

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6. Set to the previously measured actual value. Push back volume setting knob to the stop (approx. 2 mm), without turning it.
7. Replace the red distance plate.
8. Mount ejector cap (A).
9. Secure ejector cap by turning the lateral closures into position .

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Checking the volume

The gravimetric volume test of the pipette requires the following steps:

(Testing procedures are described e.g. in ISO 8655 Part 6).

Adjust volume to nominal capacity.

Conditioning the pipette:

To condition the pipette before testing, mount a new pipette tip and pipette five times with testing liquid (H₂O dist.). Discard the pipette tip.

Carrying out the test (observe instructions by the balance manufacturer):

- Mount a new pipette tip and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature into account.
- Carry out at least 10 pipetting and weighing operations.

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$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

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$$\text{CV\%} = \frac{100 s}{\bar{V}}$$

Calculation of accuracy (A%) and coefficient of variation (CV%):

Accuracy (A%) and coefficient (CV%) are calculated according to the equations used in statistical quality control.

Calculations (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad \begin{array}{l} x_i = \text{results of weighings} \\ n = \text{number of weighings} \end{array}$$

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
Cleaning

Note:

The Transferpette® requires regular maintenance, depending on frequency of use and pipetted media.

After cleaning, always check volume (see page 17).

Pipettes up to 200 µl (Grey and yellow colour code)

1. Eject the tip.
2. Unscrew pipette shaft (C) from the handle.
3. Remove plunger (A) and plunger spring (B) from the handle.
4. Unscrew pipette shaft cone (D) from the pipette shaft (C).
5. Clean parts shown on drawing with soap solution or isopropyl alcohol. Afterwards rinse with distilled water.
6. Dry the parts (max. 120 °C/248 °F).
7. Lubricate plunger very thinly ("Maintenance", page 13).
8. Reassemble parts after they have cooled down.

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
Pipettes of 200 µl upwards

1. Eject the tip.
2. Unscrew pipette shaft (C) from handle.
3. Unscrew retention bushing (F) from pipette shaft (C).
4. Remove plunger (A) and plunger spring with seal (B), (E).

Note:

Do not disassemble parts (B) and (E) any further. Take care not to deform the seal; clean it only with a soft lint-free cloth.

5. Unscrew pipette shaft cone (D) from the pipette shaft (C).
6. Clean parts shown on drawing with soap solution or isopropyl alcohol. Afterwards rinse with distilled water.
7. Dry the parts (max. 120 °C).
8. Lubricate plunger ("Maintenance", page 13).
9. Reassemble parts after they have cooled down.

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– continued –

Note:

Do not overtighten the retention bushing (F). Tighten only by hand.

10. After assembling, operate the pipetting lever five times to centre the seal.

UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

Autoclaving

For sterile work or for decontamination, the pipette parts shown here are autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285.

Note:

To protect against contamination, filter tips can also be used.

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
Attention:

It is user's responsibility to ensure effective autoclaving.

After each autoclaving:

- Check volume (see page 17).

Note:

Only autoclave parts that have been cleaned.

Pipettes up to 200 µl

(Grey and yellow colour code)

1. Eject the tip.
2. Cautiously unscrew pipette shaft (C) from handle.
3. Remove plunger (A) and spring (B) from the handle.

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
Pipettes up to 200 µl

- continued -

4. Sterilize the complete shaft, plunger and plunger spring in the autoclave.
5. Allow the parts to cool completely and assemble in the reverse order.

Pipettes of 200 µl upwards

1. Eject the tip.
2. Cautiously unscrew pipette shaft (C) from handle.
3. Autoclave the complete shaft (C) in the autoclave without having it disassembled.
4. Allow to dry and to cool completely.
5. Screw pipette shaft into handle.

Note:

The handle is not autoclavable. For disinfection you may use commercially available disinfecting solutions as long as they do not attack PP.

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Troubleshooting

Trouble	Possible Causes	Action to be taken
<ul style="list-style-type: none"> - Tip dripping (instrument leaks) - Aspirated volume too low - Delivered volume too low 	<ul style="list-style-type: none"> - Unsuitable tip - Tip not seated tightly - Pipette shaft cone loose - Seal contaminated - Seal damaged - Non-aqueous solutions pipetted - Piston damaged 	<ul style="list-style-type: none"> - Only use high-quality tips - Press tip on firmly - Tighten clockwise see "Handling", page 8 - see "Cleaning", page 19 - see "Specifications/Spares", page 26 - see "Adjustment", page 14 - see "Spares", page 26
<ul style="list-style-type: none"> - Slow suction or none 	<ul style="list-style-type: none"> - Shaft clogged 	<ul style="list-style-type: none"> - Clean shaft see "Cleaning", page 19
<ul style="list-style-type: none"> - Volume too high 	<ul style="list-style-type: none"> - Pipetting key pressed too far (into the blow-out position) before sample take up. 	<ul style="list-style-type: none"> - Operate properly see "Pipetting", page 11

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Specifications/Spares

(Spare parts can be fitted by the user)

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** Corresponde à cor das ponteiros

* Final test values relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C/68 °F) of instrument, ambience and dist H₂O, and with smooth, jerk-free operation. The absolute error margins (µl) are not exceeded at any partial volume.


** Corresponding to the color of the tips.

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Following tolerances are typically (manufacturer's experience) obtained by new instruments under optimized testing conditions (qualified operators and standardized ambience conditions):

Pipetas até 200 µl/Pipettes up to 200 µl

Pipetas acima de 200 µl/Pipettes of 200 µl upwards

Volume (µl)	Ref.	Haste de Pipeta*** Ref.	Haste de Pipeta (E) Ref.	Haste com selo E) Ref.	Cone de haste (F) Ref.	Capa ejetora Ref.	Pistão (B) Ref.	Selo com mola (D) Ref.
Volume (µl)	Cat. No.	Pipette shaft*** Cat. No.	Pipette shaft (E) Cat. No.	Pipette shaft with seal (E) Cat. No.	Pipette shaft cone (F) Cat. No.	Ejector cap Cat. No.	Piston (B) Cat. No.	Seal with spring (D) Cat. No.
5	7041 06	7032 74	-	7032 40	7033 60	10066	7275	-
10	7041 08	7032 74	-	7032 40	7033 60	10066	7275	-
20	7041 16	7032 75	-	7032 41	7033 60	10066	7276	-
25	7041 20	7032 75	-	7032 41	7033 60	10066	7276	-
50	7041 28	7032 76	-	7032 41	7033 60	10066	7276	-
100	7041 38	7032 77	-	7032 43	7033 60	10066	7277	-
200	7041 44	7032 68	-	7032 42	7033 47	10066	7267	-
200	7041 46	7032 69	7032 44	-	7033 64	10067	7278	6530 57
250	7041 48	7032 69	7032 44	-	7033 64	10067	7278	6530 57
500	7041 54	7032 71	7032 45	-	7033 66	10067	7279	6530 58
1000	7041 62	7032 73	7032 45	-	7033 66	10067	7279	6530 58
0,5 - 10	7041 02	7032 78	-	7032 38	7032 48	10048	7274	-
2 - 20	7041 03	7032 79	-	7032 39	7032 48	10048	7271	-
2 - 20	7041 04	7032 80	-	7032 39	7033 60	10068	7273	-
5 - 50	7041 72	7032 76	-	7032 41	7033 60	10068	7276	-
10 - 100	7041 74	7032 77	-	7032 43	7033 60	10068	7277	-
20 - 200	7041 78	7032 68	-	7032 42	7033 47	10068	7267	-
25 - 250	7041 76	7032 69	7032 44	-	7033 64	10069	7278	6530 57
100 - 1000	7041 80	7032 73	7032 45	-	7033 66	10069	7279	6530 58

*** haste da pipeta completa com partes (A) a (F)

*** Pipette shaft complete with parts (A) to (F)

Acessórios

Suporte para 3 instrumentos	Ref. 7032 03
Suporte para 2 x 3 instrumentos	Ref. 7032 08
Suporte de parede para 3 instrumentos	Ref. 7032 10
Óleo de silicone V 140 000 mm 2/s	Ref. 7055 02
Mola (C)	Ref. 9884
Anel de retenção (A)	Ref. 6530 92
PLT unit	Ref. 7039 70

Ponteiras **BRAND**, soltas, PP

Volume	Embalagem	Ref.
0,1 - 20 µl	2000	7320 02
0,5 - 20 µl	2000	7320 04
1 - 50 µl	2000	7320 06
2 - 200 µl	1000	7320 08
5 - 300 µl	1000	7320 10
50 - 1000 µl	1000	7320 12
0,5 - 5 ml	200	7025 95

Accessories

Rack for 3 instruments	Cat. No. 7032 03
Rack for 2 x 3 instruments	Cat. No. 7032 08
Wall support for 3 instruments	Cat. No. 7032 10
Silicon oil V 140 000 mm ² /s	Cat. No. 7055 02
Spring (C)	Cat. No. 9884
Retention bushing (A)	Cat. No. 6530 92
PLT unit	Cat. No. 7039 70

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Return for Repair

Attention:

Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside the U.S. and Canada:

- Complete the "Declaration on Absence of Health Hazards" and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

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In the U.S. and Canada:

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

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India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park, Powai
Mumbai - 400 076 (India)
Tel.: +91 22 42957790
Fax: +91 22 42957791
E-Mail: info@brand.co.in
www.brand.co.in

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Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend to check the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on www.brand.de or www.brandtech.com.

BRAND also offers you the possibility, to have your instruments calibrated by the BRAND calibrating service or the BRAND-owned DAkkS calibration service. Just send in the instruments to be calibrated accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND calibrating service) or with a DAkkS calibration certificate. For further information please contact your dealer or BRAND.

Complete ordering information is available for download at www.brand.de (see technical documentation).

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Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual.

We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:

Information for warranty please see
www.brandtech.com.

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Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations valid in your area.

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Subject to technical modification without notice.
Errors excepted.

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© Expulsão da ponteira