# Transferpette® S

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	Page
Safety Instructions	29
Limitations of Use	30
Operating and Control Elements	31
Pipetting	32
Checking the Volume	35
Accuracy Table	37
Adjustment – Easy Calibration	39
Autoclaving	40
UV sterilization	40
2 ml, 5 ml and 10 ml Filter	40
Servicing and Cleaning	41
Ordering Information · Accessories	43
Spare Parts · Accessories	44
Troubleshooting	46
Repairs	47
Contact addresses	48
Calibration Service	49
Warranty Information	50
Disposal	50

### Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe the reagent manufacturers' information.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 30). If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Only use original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visible damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 46), and contact the manufacturer if needed.

## Limitations of Use

---

### Purpose

The pipette is an air-displacement system for pipetting aqueous solutions with medium density and low to medium viscosity.

### Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- use temperature from +15 °C to +40 °C (59 °F to 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

### Operating Limitations

Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than  $\pm 1$  °C/ $\pm 1.8$  °F.


### Operating Exclusions

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene
- for liquids incompatible with polycarbonate (inspection window)
- for liquids of a very high vapor pressure
- for liquids incompatible with FKM and polyetheretherketone
- for liquids attacking polyphenyl sulfide (PPS) (50  $\mu$ l pipet, adjustable volume)

## Operating and Control Elements

### Note:

Optimum analysis results can only be obtained with quality tips. We recommend pipette tips from BRAND. For further information, refer to the accuracy table on pages 37 and 38.

## Pipetting

- 2 ml, 5 ml and 10 ml instruments should only be used with the PE filter installed (see page 40).
- Pipette tips are disposables items!

### 1. Fitting the tip

Use the correct tip according to the volume range or the color code. Ensure that the tip is securely seated.

### 2. Volume setting

- Push the volume-change protection upward to disengage (UNLOCK).
- Select the desired volume by rotating the volume-setting wheel. Avoid twisting and abrupt rotating motions during this adjustment.
- Push the volume-change protection down to re-engage (LOCK). Note: The volume-change protection tightens but does not lock volume-setting wheel.


### 3. Aspirate sample

- Press pipetting button to the first stop.
- Hold the pipette vertically and immerse the tip into the liquid.

Volume range	Immersion depth in mm	Waiting time in s
0.1 $\mu$ l - 1 $\mu$ l	1 - 2	1
> 1 $\mu$ l - 100 $\mu$ l	2 - 3	1
> 100 $\mu$ l - 1000 $\mu$ l	2 - 4	1
> 1000 $\mu$ l	3 - 6	3

- c) Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tip immersed for a few seconds.

- a) Place the pipette tip against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.

- b) Press the pipetting button slowly to the first stop and hold it down. For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- c) The blow-out stroke empties the tip completely: Press the pipetting button down to the second stop.
- d) While doing this, wipe the pipette tip against the wall of the container.
- e) Remove the pipette tip from the container wall and let the pipetting button slide back.

### 5. Ejecting the tip

Hold the pipette shaft over a suitable disposal container and press the tip ejection key to the stop.

#### Note:

ISO 8655 prescribes rinsing the pipette tip once with the sample liquid prior to the actual pipetting process.

#### Important!

Don't lay the instrument horizontal when the tip is filled. Liquid may enter and contaminate the instrument. The instrument should be stored without tip, placed upright in the supplied shelf/rack mount or bench top rack.

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

### 1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 32 for procedure).

### 2. Condition the pipette

Condition the pipette before testing by using a pipette tip to aspirate and discharge the test liquid (distilled H<sub>2</sub>O) five times. After this, discard the pipette tip.

### 3. Carry out the test

**Note:**

According to DIN EN ISO 8655-2, it is recommended to replace the tip after each measurement. This rule may be waived in accordance with the DAkkS guideline DKD-R8-1.

- a) Aspirate liquid and pipette it into the weighing vessel.
- b) Weigh the pipetted quantity with an analytical balance. (Please follow the operating manual instructions from the balance manufacturer.)
- c) Calculate the volume, taking the temperature of the test liquid into account.
- d) At least 10 pipettings and weighings in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

## Checking the Volume

### Calculation (for nominal volume)

$x_i$  = Weighing results

$n$  = Number of weighings

$Z$  = Correction factor

(e.g., 1.0029  $\mu\text{l}/\text{mg}$   
at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

### Accuracy\*

$$\text{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

$V_0$  = Nominal volume

### Coefficient of Variation\*

$$\text{CV\%} = \frac{100 s}{\bar{V}}$$

### Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

\*) = Calculation of accuracy (A%) and variation coefficient (CV%):  
A% and CV% are calculated according to the formulas for  
statistical control.

### Note:

Testing instructions (SOPs) and a demo version of the EASYCAL™ 4.0 calibration software are available for download at [www.brand.de](http://www.brand.de).

## Transferpette® S, adjustable volume

Volume range $\mu\text{l}$	Volume step $\mu\text{l}$	A* $\leq \pm \%$	CV* $\leq \%$	Increment $\mu\text{l}$	Recommended type of tip, $\mu\text{l}$
0.1 - 1	1	2	1.2	0.001	0.1 - 20
	0.5	4	2.4		
	0.1	20	12		
0.1 - 2.5	2.5	1.4	0.7	0.002	0.5 - 20
	1.25	2.5	1.5		
	0.25	12	6		
0.5 - 10	10	1	0.5	0.01	0.5 - 20
	5	1.6	1		
	1	7	4		
2 - 20	20	0.8	0.4	0.02	2 - 200
	10	1.2	0.7		
	2	5	2		
5 - 50	50	0.8	0.3	0.05	2 - 200
	25	1.2	0.5		
	5	4	2		
10 - 100	100	0.6	0.2	0.1	2 - 200
	50	0.8	0.4		
	10	3	1		
20 - 200	200	0.6	0.2	0.2	2 - 200
	100	0.8	0.3		
	20	3	0.6		
100 - 1000	1000	0.6	0.2	1	50 - 1000
	500	0.8	0.3		
	100	3	0.6		
500 - 5000	5000	0.6	0.2	5	500 - 5000
	2500	0.8	0.3		
	500	3	0.6		
1000 - 10000	10000	0.6	0.2	10	1000 - 10000
	5000	0.8	0.3		
	1000	3	0.6		

\* A = Accuracy, CV = Coefficient of Variation

Final test values related to the nominal capacity (maximum volume) indicated on the instrument, and the indicated volume steps obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.

## Accuracy Table

### Transferpette® S, Fixed volume

Volume range $\mu\text{l}$	A* $\leq \pm \%$	CV* $\leq \%$	Recommended type of tip, $\mu\text{l}$
10	1	0.5	0.5 - 20
20	0.8	0.4	2 - 200
25	0.8	0.4	2 - 200
50	0.8	0.4	2 - 200
100	0.6	0.2	2 - 200
200	0.6	0.2	2 - 200
500	0.6	0.2	50 - 1000
1000	0.6	0.2	50 - 1000
2000	0.8	0.3	500 - 5000

\* A = Accuracy, CV = Coefficient of Variation

Final test values related to the nominal capacity (maximum volume) indicated on the instrument obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.

#### Note:



The device is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation:

**DE-M** 18

Character string: DE-M (DE for Germany), framed by a rectangle, as well as the last two digits of the year in which the marking was affixed (here: 2018).

## Adjustment – Easy Calibration

The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made using the Easy Calibration Technique.

1. Check the volume, determine actual value (see page 35).
2. Remove the label window and the label. Push the hook forward, raise it slightly and then pull it back.
3. Using a paperclip or an unused pipette tip, remove the protective film (this protective film can be discarded).
4. Push the red adjustment slider completely back, raise the volume-setting wheel (decoupling) and release the adjustment slider.
5. Set the adjustment value:
  - Transferpette® S, adjustable volume: with the volume-setting wheel in the UNLOCK position, set to the previously determined actual value.
  - Transferpette® S, fixed volume: set the volume by rotating in the +/- direction. A volume check is recommended after every adjustment.
6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Re-insert the label and the label window.

### Note:

The change to the factory settings is indicated by the red adjustment slider now visible in the label window.

## Autoclaving

---

The Transferpette® S is completely autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes, according to DIN EN 285.

1. Eject the pipette tip.
2. Autoclave the complete pipette without any further disassembling.
3. Allow the Transferpette® S to completely cool and dry.

### Note:

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

### Attention:

**Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between), with the volume-change protection set fully unlocked (UNLOCK).**

If the pipette is autoclaved frequently, piston and seal should be greased with the supplied silicone grease in order to preserve smooth movement. If necessary after sterilization, tighten the connection between the hand grip and the pipette shaft.

### UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

### 2 ml, 5 ml + 10 ml Filter Transferpette® S

A hydrophobic PE filter is used as a safeguard against liquid entering the pipette.

Change the filter if it becomes wet or contaminated.

- Use a flat object such as a screwdriver.
- Remove the filter without damaging the tip cone.

Remove the filter before autoclaving!

The instrument can be operated without a filter.

## Servicing

Inspect the pipette tip cone for damage.

Inspect the piston and seal for contamination.

Test the instrument's piston seal. To do this, affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 46.

## Disassembly and cleaning

1. Unscrew the pipette shaft (S) from the hand grip.
2. Unscrew the upper part of the ejector (A) from the pipette shaft.
3. Pull the shaft (B, C and D) out of the lower part (E) of the ejector.
4. Unscrew the piston unit (B).

**Note:** Piston remains connected with piston unit (B)!

5. Remove the seal with spring (C) (this is non-removable on 1 µl, 2.5 µl and 10 µl Transferpette® S models)
6. Clean the parts shown with a mild soap solution or isopropanol and then rinse with distilled water.
7. Allow the parts to dry (max. 120 °C/248 °F).
8. Grease piston and seal with a very thin layer of supplied silicone grease.
9. Assemble the ambient temperature parts in reverse order from above. Piston unit and upper part of the ejector (A, B) should only be hand-tight.

### Servicing

Inspect the pipette tip cone for damage.

Inspect the piston and O-Ring-seal for contamination.

Test the instrument's piston seal. We recommend using the BRAND leak testing instrument PLT unit. Alternatively: to do this, affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 46.

### Disassembly and cleaning

1. Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
2. Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
3. Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
4. Remove the O-Ring-seal from the piston unit and clean it.

**Note:** Do not disassemble piston unit (G) any further!

5. Clean piston unit (G) and lower part of pipette shaft (H) with a soap solution or isopropanol and then rinse with distilled water.
6. Allow the parts to dry (max. 120 °C/ 248 °F) and to cool down.
7. Carefully lubricate the inside and outside of the O-ring and mount it on the piston.
8. Assemble the individual components in the reverse order from that described above.

(For illustration purpose only)

### Transferpette® S, fixed volume

Capacity	Description	Cat. No.
10 µl	F-10	7047 08
20 µl	F-20	7047 16
25 µl	F-25	7047 20
50 µl	F-50	7047 28
100 µl	F-100	7047 38
200 µl	F-200	7047 44
500 µl	F-500	7047 54
1000 µl	F-1000	7047 62
2000 µl	F-2000	7047 64

### Transferpette® S, adjustable volume

Capacity	Description	Cat. No.
0.1 - 1 µl	D-1	7047 68
0.1 - 2.5 µl	D-2.5	7047 69
0.5 - 10 µl	D-10	7047 70
2 - 20 µl	D-20	7047 72
5 - 50 µl	D-50	7047 73
10 - 100 µl	D-100	7047 74
20 - 200 µl	D-200	7047 78
100 - 1000 µl	D-1000	7047 80
0.5 - 5 ml	D-5000	7047 82
1 - 10 ml	D-10000	7047 84

### Bench-top rack for 6 Transferpette® S pipettes

Cat. No. 7048 05

### Shelf/rack mount for 1 Transferpette® S pipette

Cat. No. 7048 10


## Spare Parts

### Transferpette® S up to 1000 µl

Before reordering a piston unit or seal and spring for your Transferpette® S with glass piston (up to serial number 08N), please contact [info@brand.de](mailto:info@brand.de).

Parts will differ slightly depending on nominal volume of instrument. (Fig. shows spare parts for Transferpette® S 20-200 µl)

### Transferpette® S, fixed volume

Capacity	A	B	C	D	E
10 µl	7055 10	7046 01	–	7046 21*	7046 38
20 µl	7055 10	7046 02	7046 10	7046 22	7046 39
25 µl	7055 10	7046 08	7046 14	7046 22	7046 40
50 µl	7055 10	7046 54	7046 61	7046 23	7046 41
100 µl	7055 10	7046 54	7046 61	7046 23	7046 42
200 µl	7055 10	7046 55	7046 62	7046 24	7046 43
500 µl	7055 10	7046 56	7046 63	7046 25	7046 44
1000 µl	7055 10	7046 56	7046 63	7046 25	7046 45

\* The seal is permanently built into the shaft – it cannot be removed.


### Transferpette® S, adjustable volume

Capacity	A	B	C	D	E
0.1 - 1 µl	7055 10	7046 00	–	7046 20*	7046 30
0.1 - 2.5 µl	7055 10	7046 16	–	7046 18*	7046 49
0.5 - 10 µl	7055 10	7046 01	–	7046 21*	7046 31
2 - 20 µl	7055 10	7046 02	7046 10	7046 22	7046 32
5 - 50 µl	7055 10	7046 15	7046 17	7046 59	7046 65
10 - 100 µl	7055 10	7046 54	7046 61	7046 23	7046 33
20 - 200 µl	7055 10	7046 55	7046 62	7046 24	7046 34
100 - 1000 µl	7055 10	7046 56	7046 63	7046 25	7046 35

\* 0.1-1 µl / 0.1 - 2.5 µl / 0.5-10 µl including seal

### Transferpette® S, 2 ml, 5 ml und 10 ml

Parts will differ slightly depending on nominal volume of instrument.  
(Fig. shows spare parts for Transferpette® S 5 ml).

### Transferpette® S, fixed volume and adjustable volume

Capacity	F + F'	G	H	I
2 ml	7046 66	7046 06	7032 47	7046 26
0.5 - 5 ml	7046 36	7046 06	7032 47	7046 26
1 - 10 ml	7046 37	7046 07	7046 28	7046 26

### Additional accessories for Transferpette® S

Description	Cat. No.
<b>Filter</b> for Transferpette® S 2 ml + 5 ml, pack of 25.	7046 52
<b>Filter</b> for Transferpette® S 10 ml, pack of 25.	7046 53
<b>Silicone grease</b> for Transferpette® S up to 1000 µl	7055 02
<b>Silicone grease</b> for Transferpette® S 2 ml/5 ml/10 ml	7036 77
<b>Label window</b> , pack of 1	7046 50
<b>Blank labels</b> , pack of 5	7046 51
<b>PLT unit</b> Pipette leak testing unit	7039 70

## Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (instrument leaks)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Press tip on firmly
The instrument does not aspirate or aspirates too little; the discharged volume is too low.	Seal contaminated	Clean seal
	The seal or cone is damaged	Replace seal or shaft
	The piston is contaminated or damaged	Clean or replace piston
Aspiration is too slow	Shaft clogged	Clean shaft
	The filter in the 2 ml, 5 ml and 10 ml models is contaminated	Change the filter
Discharged volume is too large	Pipetting button pressed too far into the blow-out position before sample uptake	Operate properly. See 'Pipetting', page 32.
Piston is difficult to move	The piston is contaminated or needs grease	Clean and grease the piston

## Return for repair

**Important!** Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

### Outside the U.S. and Canada:

- Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

### In the U.S. and Canada:

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

### **BRAND GMBH + CO KG**

Otto-Schott-Straße 25  
97877 Wertheim (Germany)

Tel.: +49 9342 808-0

Fax: +49 9342 808-98000

E-Mail: [info@brand.de](mailto:info@brand.de)

[www.brand.de](http://www.brand.de)

### **USA and Canada:**

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)

Tel.: +1-860-767 2562

Fax: +1-860-767 2563

[www.brandtech.com](http://www.brandtech.com)

### **India:**

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park, Powai  
Mumbai - 400 076 (India)

Tel.: +91 22 42957790

Fax: +91 22 42957791

E-Mail: [info@brand.co.in](mailto:info@brand.co.in)

[www.brand.co.in](http://www.brand.co.in)

### **China:**

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B  
No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)

Tel.: +86 21 6422 2318

Fax: +86 21 6422 2268

E-Mail: [info@brand.cn.com](mailto:info@brand.cn.com)

[www.brand.cn.com](http://www.brand.cn.com)

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers you the possibility to have your instruments calibrated by the BRAND Calibration Service or the BRAND-owned DAkkS Calibration Service. Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND Calibration Service) or with a DAkkS Calibration Certificate. For further information, please contact your dealer or BRAND. Complete ordering information is available for download at [www.brand.de](http://www.brand.de) (see Technical Documentation).

## **Warranty**

---

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:

Information for warranty please see [www.brandtech.com](http://www.brandtech.com).

## **Disposal**

---

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice.

We will not be held responsible for printing or typographical errors.