

Thermocycler-Kompatibilität

Stand: 0909

Mit einem Blick auf die Tabelle sehen Sie, welche BRAND PCR-Platten mit Ihrem Thermocycler kompatibel sind. Durch Herstellerinformationen und Feedback unserer Kunden können wir die Tabelle laufend aktualisieren. Fordern Sie unverbindlich kostenlose Muster unserer PCR-Platten an (www.brand.de), um die Kompatibilität zu Ihrem Gerät sicherstellen zu können. Mit einer kurzen Information über das Ergebnis helfen Sie uns, die Tabelle zu vervollständigen.

	24-, 48-, 96-well, o. Rahmen 7814 11, 7814 15, 7813 50	96-well, erhöhter Rahmen 7813 52	96-well, ganzer Rahmen 7813 53	96-well, halber Rahmen 7814 00	384-well, ganzer Rahmen 7813 45	384-well, ganzer Rahmen 7813 47	384-well, ganzer R., starr 7813 48
Applied Biosystems							
2700	●	●		●	●	●	●
3100	●	●				●	●
3130	-	-				●	●
3700	●	●			●	●	●
3730/3730x	●	●				●	●
qPCR 5700		●					
qPCR 7000	●	●					
qPCR 7300	●	●					
qPCR 7500	●	●					
qPCR 7700	●	●					
qPCR 7900 HT		●				●	
9600	●	●		●			
9700	●	●		●	●	●	●
Amersham Biosciences							
MegaBACE 500			●				
MegaBACE 1000			●				
MegaBACE 4000					●	●	●
Biometra							
Uno	●	●	●	●			
Uno II	●	●		●	●	●	●
T1 Thermal Cycler	●	●	●	●	●	●	●
T3 Thermal Cycler	-			-			
Tgradient	●	●	●	●			
Trobot	●		●	●	●	●	●
BioRad							
iCycler	●	●	●				
MyCycler	●						
qPCR MyiQ	●		●				
qPCR iQ5	●		●				
Corbett Research							
PalmCycler 96			●				
PalmCycler 384						●	●
Eppendorf							
Mastercycler Gradient	●	●	●	●			
Mastercycler ep	●	●	●				
Mastercycler M 384					●	●	●
Mastercycler ep			●				
qPCR Realplex			●				

	24-, 48-, 96-well, o. Rahmen 7814 11, 7814 15, 7813 50	96-well, erhöhter Rahmen 7813 52	96-well, ganzer Rahmen 7813 53	96-well, halber Rahmen 7814 00	384-well, ganzer Rahmen 7813 45	384-well, ganzer Rahmen 7813 47	384-well, ganzer R., starr 7813 48
Ericomp							
Single Block	●	●					
Twin Block	●	●					
Delta Cycler	●	●					
Hybaid							
Multiblock System MBS	●		●		●	●	●
Omnigene	●	●	●	●	●	●	●
Omn-E	●	●	●	●			
PCR Express	●	●	●	●	●	●	●
PCR Sprint	-		●	-			
pxe	●		●		●	●	●
px2	●		●		●	●	●
Touchdown	●	●	●	●	●	●	●
MJ Research							
BaseStation			●				
qPCR Chromo 4			●				
Dyad/Disciple	●		●			●	●
qPCR Opticon			●				
qPCR Opticon 2			●				
PTC-100	●	●	●	●		●	●
PTC-200	●	●	●	●	●	●	●
PTC-225 Tetrad	●	●	●	●	●	●	●
MWG							
Primus 96	●	●	●	●			
Primus 384					●	●	●
Stratagene							
qPCR Mx4000	●	●					
qPCR Mx3000	●		●				
Robocycler	●	●	●	-	●	●	●
TaKaRa							
TP240			●				
TP3000	●		●				
Techne							
TC-412/Flexigene	●	-	●	●	●	●	●
Genius	●	-	●	●	●	●	●
TC-512/Touchgene	●	-	●	●	●	●	●
Gradient							
TC-3000X	●*	-	-	-	-	-	-
Transgenomic							
Wave System			●				

* kompatibel mit 7814 11 und 7814 15

● = Real Time PCR geeignet ● = geeignet - = nicht geeignet □ = noch keine Info

