



QuikSip™

Safe and fast extraction of non-critical media

BRAND. For lab. For life.®



- + Safe removal of supernatants.
- + Fingertip vacuum control using the cell-culture™ unit.
- + Suitable for use with disposable pipette tips, micropipettes and Pasteur pipettes.

Get to know the QuikSip™

The handy QuikSip™ bottle-top aspirator from BRAND is designed for the safe and quick aspiration of non-critical media in biology. Use it to easily remove residual liquids such as culture media from dishes, tissue culture flasks, and multiwell plates or supernatants after precipitation and centrifugation of proteins and nucleic acids.

The instrument can be dismantled within a few seconds for cleaning or maintenance without the need for any tools. This means that wearing parts can be replaced quickly, easily, and cost-effectively.



Easy maintenance

all wearing parts can be replaced quickly and easily

Efficient work

up to max. 25 ml with one plunger stroke

Many possible applications

use with cell culture or 8-channel pipetting unit



QuikSip™ aspirator with single channel device for removal of liquids out of reaction vessels, tubes, etc.



QuikSip™ aspirator with 8 channel manifold for removal of liquids from multiwell plates.

- + Up to max. 25 ml per plunger stroke
- + No vacuum pump required
- + Vacuum can be easily controlled with the cell-culture™ unit
- + 8-channel pipetting unit available as an accessory
- + For use with disposable pipette tips, micropipettes and glass pasteur pipettes.
- + Adapter and suction tube of the cell-culture™ are autoclavable

Ordering information



QuikSip™ BT-Aspirator

Items supplied:

1 QuikSip™ BT-Aspirator, spare dispensing cartridge and 2 PP adapters (GL 45/32 and GL 45/S 40).

Pack of	Cat. No.
1	4723180



cell-culture™-unit

Single channel device. Complete with suction tube and 3 adapters.

Pack of	Cat. No.
1	25950

Accessories



Dispensing cartridge

Piston (PE), cylinder (PP).

Pack of	Cat. No.
3	704504



8-channel manifold

PP. Autoclavable (121 °C).

Pack of	Cat. No.
1	704526



Actuator

for QuikSip™ BT-Aspirator.

Pack of	Cat. No.
1	704570



Filling tube

PP, suitable for QuikSip™ BT-Aspirator, with filling valve, PP/EPDM.

Pack of	Cat. No.
1	704575

BRAND GMBH + CO KG

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



BRAND. For lab. For life.®

BRAND®, BRAND. For lab. For life.®, cell-culture™, QuikSip™ as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. The BRANDGROUP figurative mark is registered trademark or trademark of Brand Group SE & Co. KG, Germany All other trademarks mentioned or depicted here are the property of the respective owners.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65.

Subject to technical modification without notice. Errors excepted.

994382 © 2024 BRAND GMBH + CO KG | Printed in Germany | 0724



Find accessories and replacement parts, user manuals, test instructions (SOP) and product videos at shop.brand.de



Further information on products and applications can be found on our YouTube channel: [mylabBRAND](https://www.youtube.com/mylabBRAND)



BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com