

总览

应用	材质	最低耐受温度 °C	最高耐受温度 °C	可热封	可移除	可穿刺	无DNA, DNA酶, RNA酶	灭菌	包装规格	货号
	qPCR, ELISA, 比色实验	聚酯	-40	110		✓	✓		100片	781391
 	qPCR, PCR和储存 (封膜条)	聚酯	-40	120			✓		400条 (8条/片, 50片)	781383
	PCR, ELISA, EIA 和其他光学应用	聚酯	-40	120			✓		100片	781390
 	PCR, ELISA 其他光学应用和储存	聚丙烯	-80	120					100片	701367
 	PCR 和冷冻储存	聚酯	-80	120		✓	✓		100片	701376
 	PCR 和冷冻储存	铝箔	-80	120		✓	✓		100片	701377
 	PCR 和储存 (封膜条)	铝箔	-80	120			✓	✓	300片 (6条/片, 50片)	781382
 	qPCR, PCR 存储和处理	聚乙烯	-80	110	✓		✓		100片	701431
 	PCR, 存储和运输	铝箔	-20	120	✓		✓	✓	100片	701430
 	低温储存, 避光和自动化	铝箔	-80	120			✓	✓	100片	781381
	存储和运输	铝箔	-200	120	✓	✓	✓	✓	100片	701432
	自动化	聚乙烯/ 聚丙烯	-40	90			✓		50片	701370
	自动化	聚乙烯	-40	90			✓	✓	100片	701374
	细胞/组织培养	人造丝	-20	80					100片	701364
	细胞/组织培养	人造丝	-20	80				✓	50片	701365
 	荧光检测, 储存	聚乙烯, 黑色	-40	80					50片	701371
	发光检测, 显微镜观察	聚乙烯, 白色	-40	80					50片	701372

 PCR/qPCR

 自动化

 储存和运输

 细胞培养

 荧光和发光检测