


# Cuvette, macro

## 759105

PRODUCT  
INFORMATION




Dimensions in mm

Filling volume min.	2.5 ml
Light path	10 mm
Material	PMMA

Subject to technical modification without notice. Errors excepted.

# Cuvette, semi-micro 759115

PRODUCT  
INFORMATION



Dimensions in mm

<b>Filling volume min.</b>	1.5 ml
<b>Light path</b>	10 mm
<b>Material</b>	PMMA

Subject to technical modification without notice. Errors excepted.