


50 – 1250 µl Filter tips XL

732618 | 732718



Dimensions in mm


Volume	50 – 1250 µl XL
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

50 – 1250 µl Filter tips XL

732638 | 732738



Dimensions in mm

Volume	50 – 1250 µl XL
Material	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml