


1 – 50 µl Pipette tips

732006 | 732026

PRODUCT
INFORMATION




Dimensions in mm

Volume	1 – 50 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	bulk
Quality	Produced under cleanroom conditions acc. to ISO 14 644-1 class 8, CE-marked acc. to EC-Directive 98/79 EC (IVD).

1 – 50 µl Pipette tips

732106 | 732206



Dimensions in mm


Volume	1 – 50 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

1 – 50 µl Pipette tips

732126 | 732226



Dimensions in mm

Volume	1 – 50 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP ,CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml