


# 1 – 10 ml Pipette tips

## 702603 | 702604

PRODUCT  
INFORMATION




Dimensions in mm

<b>Volume</b>	1 – 10 ml
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	bulk
<b>Quality</b>	Produced under cleanroom conditions. CE-marked acc. to EC-Directive 98/79 EC (IVD).

# 1 – 10 ml Pipette tips

## 702606

PRODUCT  
INFORMATION



Dimensions in mm

<b>Volume</b>	1 – 10 ml
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Produced under cleanroom conditions. CE-marked acc. to EC-Directive 98/79 EC (IVD).