


2 – 200 μ l Pipette tips

732008 | 732028

PRODUCT
INFORMATION




Dimensions in mm

Volume	2 – 200 μ l
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	yellow
Packaging	bulk
Quality	Produced under cleanroom conditions acc. to ISO 14 644-1 class 8, CE-marked acc. to EC-Directive 98/79 EC (IVD).

2 – 200 µl Pipette tips

732108 | 732208 | 732248



Dimensions in mm


Volume	2 – 200 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

2 – 200 µl Pipette tips

732128 | 732228 | 732268



Dimensions in mm

Volume	2 – 200 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP ,CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml