


0.5 – 5 ml Pipette tips

702595 | 702600

PRODUCT
INFORMATION




Dimensions in mm

Volume	0.5 – 5 ml
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	bulk
Quality	Produced under cleanroom conditions. CE-marked acc. to EC-Directive 98/79 EC (IVD).

0.5 – 5 ml Pipette tips

702605


 PRODUCT
INFORMATION


Dimensions in mm

Volume	0.5 – 5 ml
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Produced under cleanroom conditions. CE-marked acc. to EC-Directive 98/79 EC (IVD).