


0.1 – 1 µl Filter tips ULR 732802



PRODUCT
INFORMATION




Dimensions in mm

Volume	0.1 – 1 µl
Material	PP, Filter PE; Ultra Low Retention surface; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

0.1 – 1 µl Filter tips ULR 732822

Dimensions in mm

Volume	0.1 – 1 µl
Material	PP, Filter PE; Ultra Low Retention surface; autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP ,CE-marked acc. to EC-Directive 98/79 EC (IVD).



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml