


# 2 – 200 µl Pipette tips ULR

## 732308 | 732348



Dimensions in mm


<b>Volume</b>	2 – 200 µl
<b>Material</b>	PP; Ultra Low Retention surface; autoclavable at 121°C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, CE-marked acc. to EC-Directive 98/79 EC (IVD).



\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

# 2 – 200 µl Pipette tips ULR

## 732328 | 732368



Dimensions in mm

<b>Volume</b>	2 – 200 µl
<b>Material</b>	PP; Ultra Low Retention surface; autoclavable at 121°C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP ,CE-marked acc. to EC-Directive 98/79 EC (IVD).



\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml