

A strong team for better lab hygiene



Fleece Wipes

Pursept® Wipes XL

F I R S T C L A S S · B R A N D

Pursept® Wipes XL and Pursept®-A Xpress, to disinfect the workplace, lab table or clean bench quickly and reliably. Wipes can be stored in the convenient dispenser box.

■ **Advanced microstructure**

The lint-free and tear-proof wipes are designed for both optimal disinfectant delivery and excellent contaminant uptake.

■ **Refillable dispenser box**

Wipes are dispensed individually from the space-saving, reclosable box.

■ **Pursept®-A Xpress**

Neutralizes surface bacteria, viruses and fungi - Effective in only 15 seconds - DGHM/VAH compliant.



Hygienic safety in only 15 seconds with this strong team: Pursept® Wipes XL, convenient dispenser box and Pursept®-A Xpress disinfectant

The Pursept® Wipes XL roll can be stored in the dispenser box. For high-volume use, wipes can be soaked with 1.5 liters of liquid disinfectant directly in the handy dispenser box (lasts up to 28 days). For less frequent use, individual wipes can be soaked immediately before use.

Ordering Information

Technical Information

Pursept® Wipes XL

Wipes are made of 100% polyester, specially designed to be tear-proof and lint-free. They offer optimal disinfectant delivery together with excellent contaminant uptake.

Description	Pack of	Cat. No.
70 wipes (230 x 280 mm) per roll	6	448 70



Dispenser Box (empty) for Pursept® Wipes XL

Refillable dispenser box for practical and space-saving storage of Pursept® Wipes XL, a perfect way to keep pre-soaked wipes from drying out.

Description	Pack of	Cat. No.
Dispenser box for 70 wipes (230 x 280 mm)	3	448 75



Pursept®-A Xpress

Alcoholic surface disinfectant
Ready to use, fresh-scented solution.
DGHM/VAH listed and RKI compliant, works in only 15 seconds.

Description	Pack of	Cat. No.
1 liter bottle	10	448 30
5 liter can	3	448 35



Microbiological efficiency of the undiluted working solution:

Surface disinfection: Bacteria and fungi with little organic burden with mechanics acc. to DGHM/VAH directives: 15 seconds
Limited virucidal efficacy*
(incl. HIV, HBV, HCV): 15 seconds
Effectiveness against TbB: 30 seconds

* As recommended by the RKI, Bundesgesundheitsblatt 01/2004

Ingredients:

100 g Pursept®-A Xpress contain:
55 g ethanol,
0.03 g n-alkyl-aminopropyl-glycine,
auxiliary agents

Chemical and physical specifications:

Clear, colorless liquid
Density: approx. 0.9 g/ml
pH value: approx. 7.5

BRAND is a trademarks of BRAND GMBH + CO KG, Germany.
Other reproduced brands are the property of the respective owner.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

Subject to technical modification without notice. Errors excepted.

