

**BRAND GMBH + CO KG**

Postfach 11 55  
97861 Wertheim  
Germany  
Tel: +49 9342 808-0  
Fax: +49 9342 808-98000  
info@brand.de  
www.brand.de



BRAND GMBH + CO KG · Otto-Schott-Straße 25 · 97877 Wertheim · Germany

## Product information: BLAUBRAND® Volumetric Instruments

Dear valued customer,

Since 1995 BLAUBRAND® Volumetric Instruments by BRAND are marked with a batch number and come with a corresponding certificate. The batch number is shown on the certificate and also appears on the instrument, for instance **DE-M** 18.01.

Now we have improved the batch numbering:

**DE-M** 18.40 (manufacturing year 2018/batch 40) now reads 40 **DE-M** 18 (batch 40/ manufacturing year 2018).

With the improved identification the respective instrument complies with the **German Weights and Measures Act** and the regulations based on it. We are audited regularly by a German Authority so that we are always in compliance with this act. In the course of 2018 all BLAUBRAND® Volumetric Instruments will be changed to the new system of batch numbering. Order processes at BRAND will remain the same.

Our **high product quality** regarding

- robust blanks,
- quality printing inks blue, white and amber as well as Color-Codes,
- accuracy with AQL-level 0,4 (certainty of at least 99,6 %) and
- compliance with the respective standards for volumetric flasks, graduated & mixing cylinders, bulb pipettes, graduated pipettes, burets and automated burets

is hence secured in a regulatory manner. In addition all these features are approved in a declaration of conformity including a certificate of performance (one per pack). The improved numbering assures the use of an accurate, durable and regulatorily optimized **BLAUBRAND® Volumetric Instrument**.

Kind regards

Peter Schütte  
- Managing Director Marketing & Sales -

i. A.

Dr. Christian Schurz  
- Product Manager Glass & Plastic Ware -