

BRAND Repair and Calibration Service



BRAND GMBH + CO KG · Otto-Schott-Straße 25 · 97877 Wertheim · Germany

To:

BRAND GMBH + CO KG

– Repair and Calibration Service –

Otto-Schott-Str. 25
97877 Wertheim

Fax: +49 9342 808-91290

ORDER

Calibration is offered for instruments listed below, made by BRAND or from other manufacturers.
Repair and maintenance is provided only for BRAND instruments.

Calibration costs for sent in Liquid Handling Instruments

Prices as of 01.01.2019 All prices are recommended sales prices in € without VAT	works calibration price per unit	DAkkS-Calibration price per unit
Single-channel piston-operated pipettes, variable Dispensers (bottle-top), variable Repetitive dispensers (Stepper), Piston-operated burettes, Titration systems	60.40 €	106.60 €
Maintenance and cleaning (BRAND instruments only)	24.60 €	24.60 €
Additional certificate documenting status before maintenance	26.50 €	47.90 €
Single-channel piston-operated pipettes, fixed volume Dispensers, fixed volume	43.35 €	76.40 €
Maintenance and cleaning (BRAND instruments only)	24.60 €	24.60 €
Additional certificate documenting status before maintenance	19.90 €	36.00 €
Multichannel piston-operated pipettes, 4 or 8 channels	83.15 €	195.70 €
Maintenance and cleaning (BRAND instruments only)	53.75 €	53.75 €
Additional certificate documenting status before maintenance	51.25 €	89.60 €
Multichannel piston-operated pipettes, 12 channels	113.45 €	217.60 €
Maintenance and cleaning (BRAND instruments only)	72.05 €	72.05 €
Additional certificate documenting status before maintenance	69.50 €	119.50 €
With repetitive dispensers (Stepper) and titration systems Calibration Certificate for each additional dispensing cartridge	26.50 €	47.90 €
Pipette Leak Testing unit – PLT unit	97.50 €	
Additional certificate documenting status before maintenance	35.60 €	

With repetitive dispensers, extended operational testing is done instead of maintenance (BRAND instruments only).

Calibration costs for sent in Volumetric Instruments made of glass or plastic

Prices as of 01.01.2019 All prices are recommended sales prices in € without VAT	DAkkS-Calibration price per unit
Burettes	69.80 €
Bulb pipettes	63.80 €
Graduated pipettes	69.80 €
Volumetric flasks	63.80 €
Graduated cylinders	63.80 €
Oxygen flask, Winkler pattern (150 + 300 cm³)	63.80 €
Density bottle of glass, without thermometer (1 cm³ - 100 cm³)	65.00 €

Standard calibration is done in compliance with DIN EN ISO 8655.

Our calibration laboratory meets the requirements according to ISO 9001 and GLP/GMP guidelines (factory setting).

For DAkkS calibrations the requirements of DIN EN ISO/IEC 17025 are also met (D-K-18572-01-00, accredited by the DAkkS – Deutsche Akkreditierungsstelle GmbH).

Our General Terms and Conditions apply. They are printed in the current price list or are available for download from www.brand.de.

Delivery Address:**Billing address, if different:**

Name

Clinic / Company

Institute / Dept.

Street / P.O.

Postal Code / City

Telephone / Fax

E-Mail / URL

Name

Clinic / Company

Institute / Dept.

Street / P.O.

Postal Code / City

Telephone / Fax

E-Mail / URL

Request:

According to the price list shown on the reverse, please provide:

Calibration: with Calibration Certificate (factory setting) with DAkkS Calibration Certificate

Prepare an additional certificate documenting the metrological status before maintenance and repair.
(According to DIN EN ISO/IEC 17 025 (with DAkkS Calibration Certificate) the metrological status of the instrument has to be documented before a volume relevant intervention.)

Please perform a maintenance and cleaning (BRAND instruments only)

Adjust the instrument, if required (In this instance, you will receive information on measuring results before adjustment.)

for calibrating service only: minor repairs of a value up to € 25,- + VAT will be carried out and invoiced without further queries.

Please repair the instruments, if necessary (BRAND instruments only).

Please repair the instrument. There are defects (please enter in the table below):

Up to EUR 50,- EUR 100,- repair can be done without additional contact.

If costs will exceed this value, please provide an estimate first before proceeding.

Reminder service:

Please send us a reminder of when the next instrument calibration is due based on the information given in the table below.

Calibration is preferably carried out using BRAND disposables (tips, disposables). For instruments from other manufacturers, the corresponding original disposables will be used if available ; for third party instruments, please send consumables e.g. tips with your order if possible.

Single-channel piston-operated pipette (fixed volume)



at least 20 pipette tips for each instrument

Single-channel piston-operated pipette (variable)



at least 60 pipette tips for each instrument

Multichannel piston-operated pipette



at least 400 pipette tips for each instrument

Positive-displacement piston-operated pipette



two corresponding tips / capillaries for each instrument

Repetitive dispensers (steppers)



two corresponding dispenser tips for each instrument

Instrument Description	Serial Number	Calibration Cycle	For Repair: identified defect

Date

Signature

Declaration on the Absence of Health Hazards

We are legally bound to give our staff a maximum of protection from health hazards caused by contaminated instruments. We therefore ask for your understanding that we cannot carry out any calibration/ repair unless this declaration is submitted completed and signed.

Media used:

- | | |
|-----------------------------------|---------------------------------------|
| <input type="checkbox"/> Acids | <input type="checkbox"/> Bases |
| <input type="checkbox"/> Solvents | <input type="checkbox"/> Serum, blood |
| <input type="checkbox"/> | |

Decontamination:

If there is a risk of contamination: What decontamination method was used?

- | | |
|---------------------------------------|---|
| <input type="checkbox"/> Autoclave | <input type="checkbox"/> Ethylene oxide |
| <input type="checkbox"/> Formaldehyde | <input type="checkbox"/> |

Re: Instrument Consignment dated _____ **/ Re: Delivery Note No.** _____

The Undersigned hereby declares:

- that the instruments have been carefully cleaned and decontaminated before shipment.
- that the instruments pose no danger through bacteriological, chemical, radiological and / or viral contamination.
- to be authorised to make declarations on behalf of the Institution represented.

Company / Laboratory (Stamp)

Tel. / Fax / E-Mail

Name

Position

Date, Signature