



Free*
Boro 3.3
bottle

Titrette® promo package

Don't miss our promotional package for Titrette® with Boro 3.3 bottle for free* only from January 1 until March 31, 2022!

BRAND. For lab. For life.®

Bottle-top burettes Titrette®

Titrette® is the first bottle-top burette to meet class A error limits for glass burettes. In addition, it stands out for easy titration, a compact design, excellent precision, an exchangeable dosing unit, and an optional PC interface. You can use the Titrette® to titrate quickly and reliably, even with restricted space and independent of available power connections – in the laboratory, in production, or in the field.

Your benefits:

- + Non-slip hand wheels and a smooth-running precision gear set for fast or drip titration make handling reliable and especially simple.
- + Separate buttons for on/off, pause, and CLEAR to delete the display and select function.
- + Buy a Titrette® and get selected bottle for free

Titrette® promo package

Promotion runs from January 1 until March 31, 2022



Free reservoir bottle!*

Just order on shop.brand.de/en/promotions or fill out the form and send it back to BRAND.



o Yes, I would like to take part in the promotion and order at the special promotional price of as per BRAND GMBH + CO KG General Terms and Conditions for orders from outside shop.brand.de/en:

Titrette® bottle-top burette

DE-M marking, performance certificate, telescoping filling tube, recirculation tube, 2 batteries (AAA/UM4/LR03), 3 PP bottle adapters (GL 32, S 40, NS 29/32), 2 amber colored light shield inspection windows.

Titrette® Nominal volume	Description	Pack of	Reservoir bottle Boro 3.3 NS 29/32	Pack of	Quantity	Special price €*	Total price €*
10 ml	standard	1	clear glass, Ø 160 mm, height 200 mm	1		734,20	
			amber glass, Ø 160 mm, height 200 mm	1		734,20	
25 ml	standard	1	clear glass, Ø 160 mm, height 200 mm	1		734,20	
			amber glass, Ø 160 mm, height 200 mm	1		734,20	
50 ml	standard	1	clear glass, Ø 160 mm, height 200 mm	1		734,20	
			amber glass, Ø 160 mm, height 200 mm	1		734,20	
10 ml	with RS232 interface	1	clear glass, Ø 160 mm, height 200 mm	1		943,20	
			amber glass, Ø 160 mm, height 200 mm	1		943,20	
25 ml	with RS232 interface	1	clear glass, Ø 160 mm, height 200 mm	1		943,20	
			amber glass, Ø 160 mm, height 200 mm	1		943,20	
50 ml	with RS232 interface	1	clear glass, Ø 160 mm, height 200 mm	1		943,20	
			amber glass, Ø 160 mm, height 200 mm	1		943,20	
25 ml	acc. to Soxhlet-Henkel	1	clear glass, Ø 160 mm, height 200 mm	1		734,20	
			amber glass, Ø 160 mm, height 200 mm	1		734,20	

Total:

Date _____ Signature _____

Address:

.....
First name

.....
Last name

.....
Institute / Company

.....
Street

.....
Zip code and City

.....
E-Mail

.....
Tel./Fax

Approval

Yes, I agree that BRAND GMBH + CO KG (hereinafter "BRAND") may contact me by e-mail, telephone and/or mail for the purpose of providing information on products, news, events and free tickets, as well as for conducting satisfaction surveys and for sending samples, and that BRAND may process the data required for this purpose. My data will not be passed on to third parties. This consent is valid until I revoke it. I can revoke my consent at any time with effect for the future and without giving reasons. I can view the contact options for

my revocation in BRAND's data privacy statement at <https://www.brand.de/en/privacy-statement> or I can send my revocation to: privacy@brand.de. I am also entitled to the other rights described in the data privacy statement.

Date

Signature

* All prices are recommended sales prices without VAT. You will receive an order confirmation with two separate positions for the instrument and the accessory. This offer is valid until 31.03.22. No cash alternative available. Warranty claims remain unaffected.

BRAND®, BRAND. For lab. For life.®, Titrette® as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners. Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | P.O. Box 1155 | 97877 Wertheim | Germany
T +49 9342 808-0 | F +49 9342 808-98000 | info@brand.de | www.brand.de

