



Save now
up to
25%*

HandyStep® S

Don't miss our special price for Repetitive pipette HandyStep® S only from June 1 until December 31, 2020!

BRAND. For lab. For life.®

HandyStep® S – the repetitive pipette for so many applications!

Ideal for routine or research applications in diagnostics, molecular biology, environmental analysis, and many more. HandyStep® S and PD-tips work on the direct displacement principle. This offers the highest-precision dispensing of liquid media with high viscosity, high density, or high vapor pressure. Direct displacement permits contamination-free operation, since no aerosols are formed.

Your benefits:

- + Serial pipetting made easy
- + Works with BRAND PD-Tips and third-party tips
- + Robust, stable and durable

Promo HandyStep® S

Promotion runs from June 1 until December 31, 2020

Save up to 25%*!

Just order on shop.brand.de/en/promotions or fill out the form and send it back to BRAND.



Yes, I would like to take part in the promotion and order at the special promotional price of as per BRAND GMBH + CO KG General Terms and Conditions for orders from outside shop.brand.de/en:

HandyStep® S

Cat. No.	Description	Quantity	Special price €*	Total price €*
705110	Repetitive pipette HandyStep® S, shelf/rack mount, DE-M marking, 3 PD-Tips: 0.1 ml, 1 ml and 10 ml.		155,25	

Total:

Date

Signature

Address:

.....
First name

.....
Last name

.....
Institute / Company

.....
Street

.....
Zip code and City

.....
E-Mail

.....
Tel./Fax

Approval

Yes, I agree that BRAND GMBH + CO KG (further referred to as BRAND) may contact me by

e-mail phone mail

for the purpose of product information, news, events, free tickets, satisfaction surveys and to provide product samples (if not applicable, pls. cross out). BRAND is allowed to process my data for the intended purpose. My personal data will not be transferred to third parties.

I can withdraw my consent by short notice at any time without giving any reason and separately for each of the above communication channels. Contact information for my withdrawal is available at the BRAND Privacy Statement on <https://www.brand.de/en/privacy-statement/> or I may withdraw my consent by sending an e-mail to privacy@brand.de.

.....
Date

.....
Signature

* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2020.

BRAND®, BRAND. For lab. For life.®, HandyStep® S as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners. Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | P.O. Box 1155 | 97877 Wertheim | Germany
T +49 9342 808-0 | F +49 9342 808-98000 | info@brand.de | www.brand.de

