



Now save  
up to  
**35%\***

## accu-jet® *pro*

Don't miss our special price for the current model pipetting aid accu-jet® *pro* only from April 1 until June 30, 2021!

**BRAND. For lab. For life.®**

### Pipetting aid accu-jet® *pro*

All with one hand: select the delivery mode (gravity-delivery / blow-out) and adjust the motor speed range with your thumb; use variable button pressure for finest control of filling and delivery speed. The ergonomic handgrip, perfectly balanced design, weight of only 190 g – for fatigue-free pipetting even in prolonged operations.

The accu-jet® *pro* is supplied with two battery compartment covers. The second cover with ridges enables you to place the instrument inverted on the bench, to avoid contamination of the pipette adapter.

### Your benefits:

- + Powerful: Fills a 50 ml pipette in less than 10 seconds.
- + Sensitive: Precise selection of pipetting speed.
- + Long-lasting: up to 8 hours of continuous pipetting without recharging

### From the end of 2021: the new accu-jet® *pro*

Featuring a new design and all the advantages of the current accu-jet® *pro* model.

# Promo accu-jet® pro

Promotion runs from April 1 until June 30, 2021



## Save up to 35%\*!

Simply enter the desired quantity in the table and send the completed form to your distributor.

Yes, I would like to take part in the promotion and order at the special promotional price\* of as per BRAND GMBH + CO KG General Terms and Conditions for orders from outside shop.brand.de/en:

### accu-jet® pro Pipetting aid

with nickel-metal hydride battery, 2 battery compartment covers, wall support, AC adapter (100-240 V; 50/60 Hz) and 2 spare membrane filters 0.2 µm, sterile

Cat. No.	Description	Quantity	Special price €*	Total price €*
26300	dark blue, AC adapter Europe 100-240 V; 50/60 Hz		249,00	
26301	magenta, AC adapter Europe 100-240 V;50/60 Hz		249,00	
26302	green, AC adapter Europe 100-240 V; 50/60 Hz		249,00	
26303	royal blue, AC adapter Europe 100-240 V;50/60Hz		249,00	

**Total:**

Date

Signature

#### Address:

.....  
First name

.....  
Last name

.....  
Institute / Company

.....  
Street

.....  
Zip code and City

.....  
E-Mail

.....  
Tel./Fax

#### Approval

Yes, I agree that BRAND GMBH + CO KG (further referred to as BRAND) may contact me by

e-mail  phone  mail

for the purpose of product information, news, events, free tickets, satisfaction surveys and to provide product samples (if not applicable, pls. cross out). BRAND is allowed to process my data for the intended purpose. My personal data will not be transferred to third parties.

I can withdraw my consent by short notice at any time without giving any reason and separately for each of the above communication channels. Contact information for my withdrawal is available at the BRAND Privacy Statement on <https://www.brand.de/en/privacy-statement/> or I may withdraw my consent by sending an e-mail to [privacy@brand.de](mailto:privacy@brand.de).

.....  
Date

.....  
Signature

\* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2021.

accu-jet®, BRAND®, BRAND. For lab. For life.® as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners. Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | P.O. Box 1155 | 97877 Wertheim | Germany  
T +49 9342 808 0 | F +49 9342 808 98000 | [info@brand.de](mailto:info@brand.de) | [www.brand.de](http://www.brand.de)

