



## Transferpette® S-8-12

Do not miss our special prices for microliter pipette Transferpette® S-8-12 – only from April 1 until June 30, 2022

**BRAND. For lab. For life.®**

### Get to know the Transferpette® S-8-12.

The Transferpette® S air displacement pipette allows you to work efficiently and ergonomically with both small and large volumes. Low operating forces and a short stroke ensure relaxed pipetting, even during long pipetting series. Due to its perfect ergonomics, the Transferpette® S provides a comfortable grip in any hand position, whether you are right- or left-handed, or have big or small hands. Changing your grip is not necessary, since the microliter pipette provides true one-hand operation: set the volume, pipette, eject the tip.

### Your advantages:

- + Efficient: one-hand operation
- + Flexible: ergonomic design for every hand
- + Reliable: volume-change protection prevents accidental adjustments

# Promo Transferpette® S -8-12

Promotion period from April 1 until June 30, 2022

## Save up to 30%\*!

Just order on [shop.brand.de/en/promotions](http://shop.brand.de/en/promotions) or fill out the form and send it back to BRAND.



Yes, I would like to take part in the promotion and order at the special promotional price of as per BRAND GMBH + CO KG General Terms and Conditions for orders from outside [shop.brand.de/en](http://shop.brand.de/en):

### Transferpette® S -8-12

Items supplied: Transferpette® S -8-12, DE-M marking, performance certificate, 1 shelf/rack mount, 1 reagent reservoir, 1 mounting tool, silicone grease and 1 set of sealing rings made of silicone/ FKM.

Cat. No.	Capacity	Quantity	Promo price €*	Total price €*
<b>Transferpette® S-8</b>				
705900	0.5 - 10 µl		420.88	
705906	5 - 50 µl		420.88	
705908	10 - 100 µl		420.88	
705910	20 - 200 µl		420.88	
705912	30 - 300 µl		420.88	

Cat. No.	Capacity	Quantity	Promo price €*	Total price €*
<b>Transferpette® S-12</b>				
705920	0.5 - 10 µl		482.34	
705926	5 - 50 µl		482.34	
705928	10 - 100 µl		482.34	
705930	20 - 200 µl		482.34	
705932	30 - 300 µl		482.34	

**Total:**

Date

Signature

#### Address:

.....  
First name

.....  
Last name

.....  
Institute / Company

.....  
Street

.....  
Zip code and City

.....  
E-Mail

.....  
Tel./Fax

#### Approval

Yes, I agree that BRAND GMBH + CO KG (hereinafter "BRAND") may contact me by e-mail, telephone and/or mail for the purpose of providing information on products, news, events and free tickets, as well as for conducting satisfaction surveys and for sending samples, and that BRAND may process the data required for this purpose. My data will not be passed on to third parties.

This consent is valid until I revoke it. I can revoke my consent at any time with effect for the future and without giving reasons. I can view the contact options for my revocation in BRAND's data privacy statement at <https://www.brand.de/en/privacy-statement> or I can send my revocation to: [privacy@brand.de](mailto:privacy@brand.de). I am also entitled to the other rights described in the data privacy statement.

.....  
Date

.....  
Signature

\* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2022.

BRAND®, BRAND. For lab. For life.®, Transferpette®, as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | P.O. Box 1155 | 97861 Wertheim | Germany  
T +49 9342 808-0 | F +49 9342 808-98000 | [info@brand.de](mailto:info@brand.de) | [www.brand.de](http://www.brand.de)

